

## Value sheet of Mindray BS Measurement System

Русский : Таблица результатов для системы BS компании Mindray

Português : Planilha de valores do Sistema de Medição BS da Mindray

Español : Hoja de valores del sistema de medición Mindray BS

Italiano : Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe : Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский : Данные совпадают во всех группах.

Português : A dados de cada grupo é a mesma.

Español : la datos de cada grupo es la misma.

Italiano : la dati di ogni gruppo è la stessa.

Türkçe : her grubun veri aynıdır.

1.**BS-200E**: BS-200E, BS-220E;

7.**BS-410**: BS-410, BS-470;

2.**BS-300**: BS-300, BS-320;

8.**BS-480**: BS-480, BS-490;

3.**BS-330**: BS-330, BS-350;

9.**BS-600**: BS-600, BS-620;

4.**BS-330E**:**BS-330E**(Serial Number starts with "XQ-"),

10.**BS-600M**: BS-600M, BS-620M;

**BS-350E**(Serial Number starts with "XS-");

11.**BS-800**: BS-800, BS-820, BS-800M, BS-820M,  
BS-1800, BS-1800plus;

5.**BS-380**: BS-380, BS-390;

12.**BS-870**: BS-870;

6.**BS-400**: BS-400, BS-420;

13.**BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M.

14.**BS-2800M**: BS-2600M, BS-2800M

**LOT** : 061622003

**EXP** : 2023-07-11

Abbreviated name	AFU Control(L)					AFU Control(H)						
	Model	Unit	Assay Value	Range (Assay Value $\pm$ 3SD)	1 SD	Model	Unit	Assay Value	Range (Assay Value $\pm$ 3SD)	1 SD		
AFU	<b>BS-200E</b>	<sup>1</sup>	U/L	22.5	18.0 — 27.0	1.5	<b>BS-200E</b>	<sup>1</sup>	U/L	62.0	49.5 — 74.5	4.2
	<b>BS-300</b>	<sup>2</sup>	U/L	25.9	20.7 — 31.1	1.7	<b>BS-300</b>	<sup>2</sup>	U/L	69.2	55.3 — 83.1	4.6
	<b>BS-330</b>	<sup>3</sup>	U/L	26.2	20.9 — 31.5	1.8	<b>BS-330</b>	<sup>3</sup>	U/L	71.0	56.7 — 85.3	4.8
	<b>BS-330E</b>	<sup>4</sup>	U/L	22.5	18.0 — 27.0	1.5	<b>BS-330E</b>	<sup>4</sup>	U/L	62.0	49.5 — 74.5	4.2
	<b>BS-380</b>	<sup>5</sup>	U/L	24.6	19.7 — 29.5	1.6	<b>BS-380</b>	<sup>5</sup>	U/L	67.6	54.0 — 81.2	4.5
	<b>BS-400</b>	<sup>6</sup>	U/L	23.9	19.1 — 28.7	1.6	<b>BS-400</b>	<sup>6</sup>	U/L	64.8	51.8 — 77.8	4.3
	<b>BS-410</b>	<sup>7</sup>	U/L	25.9	20.7 — 31.1	1.7	<b>BS-410</b>	<sup>7</sup>	U/L	70.0	55.9 — 84.1	4.7
	<b>BS-480</b>	<sup>8</sup>	U/L	24.3	19.4 — 29.2	1.6	<b>BS-480</b>	<sup>8</sup>	U/L	65.8	52.6 — 79.0	4.4
	<b>BS-600</b>	<sup>9</sup>	U/L	23.5	18.8 — 28.2	1.6	<b>BS-600</b>	<sup>9</sup>	U/L	63.8	51.0 — 76.6	4.3
	<b>BS-600M</b>	<sup>10</sup>	U/L	28.5	22.8 — 34.2	1.9	<b>BS-600M</b>	<sup>10</sup>	U/L	77.9	62.2 — 93.6	5.2
	<b>BS-800</b>	<sup>11</sup>	U/L	22.7	18.1 — 27.3	1.5	<b>BS-800</b>	<sup>11</sup>	U/L	62.0	49.5 — 74.5	4.2
	<b>BS-870</b>	<sup>12</sup>	U/L	23.5	18.8 — 28.2	1.6	<b>BS-870</b>	<sup>12</sup>	U/L	64.4	51.5 — 77.3	4.3
	<b>BS-2000</b>	<sup>13</sup>	U/L	26.7	21.3 — 32.1	1.8	<b>BS-2000</b>	<sup>13</sup>	U/L	73.2	58.5 — 87.9	4.9
	<b>BS-2800M</b>	<sup>14</sup>	U/L	27.5	22.0 — 33.0	1.8	<b>BS-2800M</b>	<sup>14</sup>	U/L	75.0	59.9 — 90.1	5.0

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value $\pm$ 3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа $\pm$ 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise $\pm$ 3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo $\pm$ 3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione $\pm$ 3SDS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri $\pm$ 3SD)

AFU	
English	$\alpha$ -L-Fucosidase
Русский	$\alpha$ -L-фукозидаза
Português	$\alpha$ -L-Fucosidase
Español	$\alpha$ -L-fucosidasa
Italiano	$\alpha$ -L-Fucosidasi
Türkçe	$\alpha$ -L-Fukosidazdir