

**Value sheet of Mindray BS Measurement System**

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-200: BS-200, BS-220, BS-330, BS-350;

8.BS-480: BS-480, BS-490;

2.BS-200E: BS-200E, BS-220E,;

9.BS-600: BS-600, BS-620;

3.BS-300: BS-300, BS-320;

10.BS-600M: BS-600M, BS-620M;

4.BS-330E:BS-330E(Serial Number starts with "XQ-"), BS-350E(Serial Number starts with "XS-")

11.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

5.BS-380: BS-380, BS-390;

12.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

6.BS-400: BS-400, BS-420;

13.BS-2800M:BS-2600M, BS-2800M;

7.BS-430: BS-430, BS-450,BS-460;

14.For applicable analyzers of the analyte, please subject to the parameter sheet and instrument.

**LOT : 061624003**

**🕒: 2025-03-27**

| English   | Abbreviated name         | Model   | Unit    | Assay Value        | Range (Assay Value±3SD)                    |
|-----------|--------------------------|---------|---------|--------------------|--|
| Русский   | сокращенное наименование | модель  | Прибор  | Результат анализа  | Диапазон(результат анализа ± 3CO)          |
| Português | Nome abreviado           | Modelo  | Unidade | Valores da análise | Faixa(Valores da análise ±3SD)             |
| Español   | nombre abreviado         | modelo  | Unidad  | Valor de ensayo    | Rango(Valor de ensayo ±3SD)                |
| Italiano  | abbreviazione            | modelli | Unità   | Valori di dosaggio | Intervallo (valore di concentrazione±3 DS) |
| Türkçe    | kısaltılmış ad           | model   | Ünite   | Tayin Değeri       | Aralık (Tayin Değeri±3SD)                  |

**AFU**

|           |                  |
|-----------|------------------|
| English   | α-L-Fucosidase   |
| Русский   | α-L-фукозидаза   |
| Português | α-L-Fucosidase   |
| Español   | α-L-fucosidasa   |
| Italiano  | α-L-Fucosidasi   |
| Türkçe    | α-L-Fukosidazdır |

| Abbreviated name | Unit | Level: L/N             |             |                         |   |      | Level: H/P |                        |                         |      |      |      |     |
|------------------|------|------------------------|-------------|-------------------------|---|------|------------|------------------------|-------------------------|------|------|------|-----|
|                  |      | Model                  | Assay Value | Range (Assay Value±3SD) |   | 1 SD | Model      | Assay Value            | Range (Assay Value±3SD) |      | 1 SD |      |     |
| AFU              | U/L  | BS-200 <sup>1</sup>    | 25.8        | 20.7                    | — | 30.9 | 1.7        | BS-200 <sup>1</sup>    | 63.9                    | 51.0 | —    | 76.8 | 4.3 |
|                  |      | BS-200E <sup>2</sup>   | 20.8        | 16.6                    | — | 25.0 | 1.4        | BS-200E <sup>2</sup>   | 52.2                    | 41.7 | —    | 62.7 | 3.5 |
|                  |      | BS-300 <sup>3</sup>    | 25.2        | 20.1                    | — | 30.3 | 1.7        | BS-300 <sup>3</sup>    | 61.7                    | 49.4 | —    | 74.0 | 4.1 |
|                  |      | BS-330E <sup>4</sup>   | 20.8        | 16.6                    | — | 25.0 | 1.4        | BS-330E <sup>4</sup>   | 52.2                    | 41.7 | —    | 62.7 | 3.5 |
|                  |      | BS-380 <sup>5</sup>    | 22.2        | 17.7                    | — | 26.7 | 1.5        | BS-380 <sup>5</sup>    | 55.4                    | 44.3 | —    | 66.5 | 3.7 |
|                  |      | BS-400 <sup>6</sup>    | 23.0        | 18.5                    | — | 27.5 | 1.5        | BS-400 <sup>6</sup>    | 57.4                    | 46.0 | —    | 68.8 | 3.8 |
|                  |      | BS-430 <sup>7</sup>    | 25.1        | 20.0                    | — | 30.2 | 1.7        | BS-430 <sup>7</sup>    | 61.8                    | 49.5 | —    | 74.1 | 4.1 |
|                  |      | BS-480 <sup>8</sup>    | 23.2        | 18.4                    | — | 28.0 | 1.6        | BS-480 <sup>8</sup>    | 58.3                    | 46.6 | —    | 70.0 | 3.9 |
|                  |      | BS-600 <sup>9</sup>    | 23.0        | 18.5                    | — | 27.5 | 1.5        | BS-600 <sup>9</sup>    | 56.7                    | 45.3 | —    | 68.1 | 3.8 |
|                  |      | BS-600M <sup>10</sup>  | 28.1        | 22.4                    | — | 33.8 | 1.9        | BS-600M <sup>10</sup>  | 69.9                    | 55.8 | —    | 84.0 | 4.7 |
|                  |      | BS-800 <sup>11</sup>   | 22.6        | 18.1                    | — | 27.1 | 1.5        | BS-800 <sup>11</sup>   | 55.6                    | 44.5 | —    | 66.7 | 3.7 |
|                  |      | BS-2000 <sup>12</sup>  | 26.2        | 20.8                    | — | 31.6 | 1.8        | BS-2000 <sup>12</sup>  | 64.6                    | 51.7 | —    | 77.5 | 4.3 |
|                  |      | BS-2800M <sup>13</sup> | 26.3        | 20.9                    | — | 31.7 | 1.8        | BS-2800M <sup>13</sup> | 65.5                    | 52.3 | —    | 78.7 | 4.4 |

**mindray**  
AFU Control

For use on: BS-2000, BS-2800M,  
BS-600M, BS-800

LOT 061624003

2025-03-27

