

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

1. **BS-200**: BS-200, BS-220;
2. **BS-200E**: BS-200E, BS-220E;
3. **BS-230**: BS-230, BS-240;
4. **BS-240E**: BS240E, BS240Pro;
5. **BS-300**: BS-300, BS-320;
6. **BS-330E**: BS-330E (Serial Number starts with "XQ-"), BS-350E (Serial Number starts with "XS-");
7. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E (V35.00) (Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E (V35.00) (Serial Number starts with "W9-" and software version starts with "35.00");

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

8. **BS-380**: BS-380, BS-390;
9. **BS-400**: BS-400, BS-420;
10. **BS-430**: BS-430, BS-450, BS-460;
11. **BS-480**: BS-480, BS-490;
12. **BS-600**: BS-600, BS-620;
13. **BS-600M**: BS-600M, BS-620M;
14. **BS-800**: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;
15. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;
16. **BS-2800M**: BS-2600M, BS-2800M.

17. **S0: 0.9% NaCl, Conc. Of S0=0;**

Português: S0: 0,9% NaCl, Conc. de S0=0;

Italiano: S0: 0,9% NaCl, conc. di S0=0;

Русский: S0: 0, 9% NaCl, конц. S0=0;

Español: S0: 0,9% NaCl, Conc. de S0=0;

Türkçe: S0: %0,9 NaCl, S0 Kons.=0.

LOT : 145223001 : 2025-02-01

Abbreviated name			CHE		
Model	Unit	Calibration Value ¹⁷	Model	Unit	Calibration Value ¹⁷
BS-200 ¹	U/L	5114	BS-400 ⁹	U/L	5114
	µkat/L	85.4		µkat/L	85.4
BS-200E ²	U/L	5114	BS-430 ¹⁰	U/L	5114
	µkat/L	85.4		µkat/L	85.4
BS-230 ³	U/L	5114	BS-480 ¹¹	U/L	5114
	µkat/L	85.4		µkat/L	85.4
BS-240E ⁴	U/L	5114	BS-600 ¹²	U/L	5114
	µkat/L	85.4		µkat/L	85.4
BS-300 ⁵	U/L	5114	BS-600M ¹³	U/L	5147
	µkat/L	85.4		µkat/L	86.0
BS-330E ⁶	U/L	5114	BS-800 ¹⁴	U/L	5114
	µkat/L	85.4		µkat/L	85.4
BS-360E ⁷	U/L	5114	BS-2000 ¹⁵	U/L	5112
	µkat/L	85.4		µkat/L	85.4
BS-380 ⁸	U/L	5114	BS-2800M ¹⁶	U/L	5147
	µkat/L	85.4		µkat/L	86.0

English	Abbreviated name	Model	Unit	Calibration Value
Русский	сокращенное наименование	модель	Прибор	Принцип калибровки
Português	Nome abreviado	Modelo	Unidade	Regra de calibração
Español	nombre abreviado	modelo	Unidad	Regla de calibración
Italiano	abbreviazione	modelli	Unità	Regola di calibrazione
Türkçe	kısaltılmış ad	model	Ünite	Kalibrasyon Kuralı

CHE	
English	Cholinesterase
Русский	Холинэстераза
Português	Colinesterase
Español	Colinesterasa
Italiano	colinesterasi
Türkçe	Kolinesteraz

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

1. **BS-200**: BS-200, BS-220;
2. **BS-200E**: BS-200E, BS-220E;
3. **BS-230**: BS-230, BS-240;
4. **BS-240E**: BS240E, BS240Pro;
5. **BS-300**: BS-300, BS-320;
6. **BS-330E**: BS-330E (Serial Number starts with "XQ-"), BS-350E (Serial Number starts with "XS-");
7. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E (V35.00) (Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E (V35.00) (Serial Number starts with "W9-" and software version starts with "35.00");

8. **BS-380**: BS-380, BS-390;
9. **BS-400**: BS-400, BS-420;
10. **BS-430**: BS-430, BS-450, BS-460;
11. **BS-480**: BS-480, BS-490;
12. **BS-600**: BS-600, BS-620;
13. **BS-600M**: BS-600M, BS-620M;
14. **BS-800**: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;
15. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;
16. **BS-2800M**: BS-2600M, BS-2800M.

LOT : 145223001 : 2025-02-01

CHE													
Model	Unit	Assay Value	Range (Assay Value ± 3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value ± 3SD)		1 SD		
BS-200	1 U/L	5757	4600	—	6914	386	BS-400	9 U/L	5804	4637	—	6971	389
	µkat/L	96.1	76.8	—	115.5	6.4		µkat/L	96.9	77.4	—	116.4	6.5
BS-200E	2 U/L	5663	4525	—	6801	379	BS-430	10 U/L	5842	4668	—	7016	391
	µkat/L	94.6	75.6	—	113.6	6.3		µkat/L	97.6	78.0	—	117.2	6.5
BS-230	3 U/L	5797	4632	—	6962	388	BS-480	11 U/L	5774	4613	—	6935	387
	µkat/L	96.8	77.4	—	116.3	6.5		µkat/L	96.4	77.0	—	115.8	6.5
BS-240E	4 U/L	5733	4581	—	6885	384	BS-600	12 U/L	5842	4668	—	7016	391
	µkat/L	95.7	76.5	—	115.0	6.4		µkat/L	97.6	78.0	—	117.2	6.5
BS-300	5 U/L	5677	4536	—	6818	380	BS-600M	13 U/L	5796	4631	—	6961	388
	µkat/L	94.8	75.8	—	113.9	6.3		µkat/L	96.8	77.3	—	116.2	6.5
BS-330E	6 U/L	5663	4525	—	6801	379	BS-800	14 U/L	5842	4668	—	7016	391
	µkat/L	94.6	75.6	—	113.6	6.3		µkat/L	97.6	78.0	—	117.2	6.5
BS-360E	7 U/L	5721	4571	—	6871	383	BS-2000	15 U/L	5830	4658	—	7002	391
	µkat/L	95.5	76.3	—	114.7	6.4		µkat/L	97.4	77.8	—	116.9	6.5
BS-380	8 U/L	5804	4637	—	6971	389	BS-2800M	16 U/L	5796	4631	—	6961	388
	µkat/L	96.9	77.4	—	116.4	6.5		µkat/L	96.8	77.3	—	116.2	6.5

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value ± 3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон (результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa (Valores da análise ± 3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango (Valor de ensayo ± 3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore diconcentrazione ± 3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri ± 3SD)

CHE	
English	Cholinesterase
Русский	Холинэстераза
Português	Colinesterase
Español	Colinesterasa
Italiano	colinesterasi
Türkçe	Kolinesteraz