

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

1. **BS-200**: BS-200, BS-220;
2. **BS-200E**: BS-200E, BS-220E;
3. **BS-230**: BS-230, BS-240;
4. **BS-240E**: BS240E, BS240Pro;
5. **BS-300**: BS-300, BS-320;
6. **BS-330E**: BS-330E (Serial Number starts with "XQ-"), BS-350E (Serial Number starts with "XS-");
7. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E (V35.00) (Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E (V35.00) (Serial Number starts with "W9-" and software version starts with "35.00");

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

8. **BS-380**: BS-380, BS-390;
9. **BS-400**: BS-400, BS-420;
10. **BS-430**: BS-430, BS-450, BS-460;
11. **BS-480**: BS-480, BS-490;
12. **BS-600**: BS-600, BS-620;
13. **BS-600M**: BS-600M, BS-620M;
14. **BS-800**: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;
15. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;
16. **BS-2800M**: BS-2600M, BS-2800M.

17. **S0: 0.9% NaCl, Conc. Of S0=0;**

Português: S0: 0,9% NaCl, Conc. de S0=0;

Italiano: S0: 0,9% NaCl, conc. di S0=0;

Русский: S0: 0, 9% NaCl, конц. S0=0;

Español: S0: 0,9% NaCl, Conc. de S0=0;

Türkçe: S0: %0,9 NaCl, S0 Kons.=0.

LOT : 145223007 : 2025-08-17

Abbreviated name			CHE		
Model	Unit	Calibration Value ¹⁷	Model	Unit	Calibration Value ¹⁷
BS-200 ¹	U/L	5156	BS-400 ⁹	U/L	5156
	µkat/L	86.1		µkat/L	86.1
BS-200E ²	U/L	5156	BS-430 ¹⁰	U/L	5156
	µkat/L	86.1		µkat/L	86.1
BS-230 ³	U/L	5156	BS-480 ¹¹	U/L	5156
	µkat/L	86.1		µkat/L	86.1
BS-240E ⁴	U/L	5156	BS-600 ¹²	U/L	5156
	µkat/L	86.1		µkat/L	86.1
BS-300 ⁵	U/L	5156	BS-600M ¹³	U/L	5188
	µkat/L	86.1		µkat/L	86.6
BS-330E ⁶	U/L	5156	BS-800 ¹⁴	U/L	5156
	µkat/L	86.1		µkat/L	86.1
BS-360E ⁷	U/L	5156	BS-2000 ¹⁵	U/L	5205
	µkat/L	86.1		µkat/L	86.9
BS-380 ⁸	U/L	5156	BS-2800M ¹⁶	U/L	5188
	µkat/L	86.1		µkat/L	86.6

English	Abbreviated name	Model	Unit	Calibration Value
Русский	сокращенное наименование	модель	Прибор	Принцип калибровки
Português	Nome abreviado	Modelo	Unidade	Regra de calibração
Español	nombre abreviado	modelo	Unidad	Regla de calibración
Italiano	abbreviazione	modelli	Unità	Regola di calibrazione
Türkçe	kısaltılmış ad	model	Ünite	Kalibrasyon Kuralı

CHE	
English	Cholinesterase
Русский	Холинэстераза
Português	Colinesterase
Español	Colinesterasa
Italiano	colinesterasi
Türkçe	Kolinesteraz

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

1. **BS-200**: BS-200, BS-220;

2. **BS-200E**: BS-200E, BS-220E;

3. **BS-230**: BS-230, BS-240;

4. **BS-240E**: BS240E, BS240Pro;

5. **BS-300**: BS-300, BS-320;

6. **BS-330E**: BS-330E (Serial Number starts with "XQ-"),
BS-350E (Serial Number starts with "XS-");

7. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,
BS-330E (V35.00) (Serial Number starts with "W8-"
and software version starts with "35.00"),
BS-350E (V35.00) (Serial Number starts with "W9-"
and software version starts with "35.00");

8. **BS-380**: BS-380, BS-390;

9. **BS-400**: BS-400, BS-420;

10. **BS-430**: BS-430, BS-450, BS-460;

11. **BS-480**: BS-480, BS-490;

12. **BS-600**: BS-600, BS-620;

13. **BS-600M**: BS-600M, BS-620M;

14. **BS-800**: BS-800, BS-820, BS-800M, BS-820M,
BS-1800, BS-1800plus;

15. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

16. **BS-2800M**: BS-2600M, BS-2800M.

LOT : 145223007 : 2025-08-17

CHE													
Model	Unit	Assay Value	Range (Assay Value ± 3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value ± 3SD)		1 SD		
BS-200	1 U/L	5880	4698	—	7062	394	BS-400	9 U/L	5905	4717	—	7093	396
	µkat/L	98.2	78.5	—	117.9	6.6		µkat/L	98.6	78.8	—	118.5	6.6
BS-200E	2 U/L	5709	4560	—	6858	383	BS-430	10 U/L	5934	4740	—	7128	398
	µkat/L	95.3	76.2	—	114.5	6.4		µkat/L	99.1	79.2	—	119.0	6.6
BS-230	3 U/L	5901	4716	—	7086	395	BS-480	11 U/L	5932	4741	—	7123	397
	µkat/L	98.5	78.8	—	118.3	6.6		µkat/L	99.1	79.2	—	119.0	6.6
BS-240E	4 U/L	5775	4614	—	6936	387	BS-600	12 U/L	5934	4740	—	7128	398
	µkat/L	96.4	77.1	—	115.8	6.5		µkat/L	99.1	79.2	—	119.0	6.6
BS-300	5 U/L	5912	4724	—	7100	396	BS-600M	13 U/L	5915	4727	—	7103	396
	µkat/L	98.7	78.9	—	118.6	6.6		µkat/L	98.8	78.9	—	118.6	6.6
BS-330E	6 U/L	5709	4560	—	6858	383	BS-800	14 U/L	5934	4740	—	7128	398
	µkat/L	95.3	76.2	—	114.5	6.4		µkat/L	99.1	79.2	—	119.0	6.6
BS-360E	7 U/L	5780	4619	—	6941	387	BS-2000	15 U/L	5956	4759	—	7153	399
	µkat/L	96.5	77.1	—	115.9	6.5		µkat/L	99.5	79.5	—	119.5	6.7
BS-380	8 U/L	5905	4717	—	7093	396	BS-2800M	16 U/L	5915	4727	—	7103	396
	µkat/L	98.6	78.8	—	118.5	6.6		µkat/L	98.8	78.9	—	118.6	6.6

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value ± 3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон (результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa (Valores da análise ± 3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango (Valor de ensayo ± 3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore diconcentrazione ± 3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri ± 3SD)

CHE	
English	Cholinesterase
Русский	Холинэстераза
Português	Colinesterase
Español	Colinesterasa
Italiano	colinesterasi
Türkçe	Kolinesteraz