

# Coagulation Control

## Reference Value Table

Product Name	Batch code And Use-by date	Applicable Instruments	Abbreviation/Analyte	Unit	Reference Value	Scope
Coagulation control (level 1)	LOT	CX-9000, CX-9010, CX-9200 and CX-9210 Auto Coagulation Analyzers of Shenzhen Mindray Bio-medical Electronics Co., Ltd	PT	s	11.84	9.47 ~ 14.21
	8623060118		APTT	s	30.71	24.57 ~ 36.85
			TT	s	15.18	12.14 ~ 18.22
	2024-06-09		FIB	g/L	3.17	2.54 ~ 3.80
				mg/dL	317	254 ~ 380

Product Name	Batch code And Use-by date	Applicable Instruments	Abbreviation/Analyte	Unit	Reference Value	Scope
Coagulation control (level 1)	LOT	C3510 and ExC810 Auto Coagulation Analyzers of Beijing Mindray Medical Instrument Co., Ltd	PT	s	13.59	10.87 ~ 16.31
	8623060118		APTT	s	32.19	25.75 ~ 38.63
			TT	s	14.50	11.60 ~ 17.40
	2024-06-09		FIB	g/L	3.02	2.42 ~ 3.62
				mg/dL	302	242 ~ 362