

# Coagulation Control

## Reference Value Table

Product Name	Batch code And Use-by date	Applicable Instruments	Abbreviation/Analyte	Unit	Reference Value	Scope
Coagulation control (level 1)	LOT	CX-9000, CX-9010, CX-9200 and CX-9210 Auto Coagulation Analyzers of Shenzhen Mindray Bio-medical Electronics Co., Ltd	PT	s	11.42	9.14 ~ 13.70
	8623080118		APTT	s	28.27	22.62 ~ 33.92
			TT	s	15.49	12.39 ~ 18.59
	2024-08-11		FIB	g/L	2.87	2.30 ~ 3.44
				mg/dL	287	230 ~ 344

Product Name	Batch code And Use-by date	Applicable Instruments	Abbreviation/Analyte	Unit	Reference Value	Scope
Coagulation control (level 1)	LOT	C3510 and ExC810 Auto Coagulation Analyzers of Beijing Mindray Medical Instrument Co., Ltd	PT	s	13.58	10.86 ~ 16.30
	8623080118		APTT	s	29.75	23.80 ~ 35.70
			TT	s	14.82	11.86 ~ 17.78
	2024-08-11		FIB	g/L	2.72	2.18 ~ 3.26
				mg/dL	272	218 ~ 326