

Coagulation Control

Reference Value Table

Product Name	Batch code And Use-by date	Applicable Instruments	Abbreviation/Analyte	Unit	Reference Value	Scope
Coagulation control (level 2)	LOT	CX-9000, CX-9010, CX-9200 and CX-9210 Auto Coagulation Analyzers of Shenzhen Mindray Bio-medical Electronics Co., Ltd	PT	s	36.48	29.18 ~ 43.78
	8723040118		APTT	s	62.53	50.02 ~ 75.04
			TT	s	/	/
	2024-04-10		FIB	g/L	/	/
		mg/dL		/	/	

Product Name	Batch code And Use-by date	Applicable Instruments	Abbreviation/Analyte	Unit	Reference Value	Scope
Coagulation control (level 2)	LOT	C3510 and ExC810 Auto Coagulation Analyzers of Beijing Mindray Medical Instrument Co., Ltd	PT	s	32.21	25.77 ~ 38.65
	8723040118		APTT	s	67.98	54.38 ~ 81.58
			TT	s	/	/
	2024-04-10		FIB	g/L	/	/
		mg/dL		/	/	