

## Reference Value Table

Product Name	Batch code And Use-by date	Applicable Instruments	Abbreviation /Analyte	Unit	Reference Value	Scope
Coagulation control (level 2)	LOT	CX series Auto Coagulation Analyzers of Shenzhen Mindray Bio-medical Electronics Co., Ltd	PT	s	33.77	27.02 ~ 40.52
	8724020218		APTT	s	56.47	45.18 ~ 67.76
			TT	s	/	/
	2025-02-27		FIB	g/L	/	/
				mg/dL	/	/

Product Name	Batch code And Use-by date	Applicable Instruments	Abbreviation /Analyte	Unit	Reference Value	Scope
Coagulation control (level 2)	LOT	C3510 and ExC810 Auto Coagulation Analyzers of Beijing Mindray Medical Instrument Co., Ltd	PT	s	30.56	24.45 ~ 36.67
	8724020218		APTT	s	58.37	46.70 ~ 70.04
			TT	s	/	/
	2025-02-27		FIB	g/L	/	/
				mg/dL	/	/