

# C-reaction Protein Control



## Value sheet of Mindray Measurement System



<b>LOT</b> : 061423003 <b>EXP</b> : 2024-03-18					
Parameter	Unit	Analyzer	Concentration Level	Assay Value	Range
C-reaction Protein	mg/L	BC-5390 CRP,BC-5180 CRP, BC-5310 CRP, BC-5100 CRP	I	2.86	2.00 - 3.72
			II	43.59	30.51 - 56.67
		BC-7500[B] CRP,BC-7500[N] CRP, BC-7500[R] CRP,BC-7500[NR] CRP, BC-7500[B] CS,BC-7500[N] CS, BC-7500[R] CS,BC-7500[NR] CS, BC-7300[B] CRP,BC-7300[N] CRP, BC-7300[R] CRP,BC-7300[NR] CRP	I	3.84	2.69 - 4.99
			II	49.02	34.31 - 63.73
		CRP-M100	I	2.86	2.00 - 3.72
			II	43.59	30.51 - 56.67
		BC-700[B] CS BC-760[B] CS BC-700[R] CS BC-760[R] CS	I	3.74	2.62 - 4.86
			II	49.84	34.89 - 64.79
		BP 200n BP 260n	I	/	/ - /
			II	/	/ - /