

C-reactive Protein (CRP) Calibrator

Value sheet of Mindray Measurement System



LOT: 158323004		🕒: 2024-07-02					
Parameter	Level	Unit	Assay Value				
			BC-5390 CRP/ BC-5180 CRP/ BC-5310 CRP/ BC-5100 CRP	CRP-M100	BC-7500[B] CRP/BC-7500[N] CRP/BC-7500[R] CRP/ BC-7500[NR] CRP/ BC-7300[B] CRP/BC-7300[N] CRP/ BC-7300[R] CRP/BC-7300[NR] CR /BC-7500[B] CS/ BC-7500[N] CS/BC-7500[R] CS/BC-7500[NR] CS	BC-700[B] CS/ BC-700[R] CS/ BC-760[B] CS/ BC-760[R] CS	BP 200n/BP 260n
C-reactive Protein (CRP)	Cal a	mg/L	5.87	5.87	7.48	7.60	7.48
	Cal b		20.61	20.61	23.46	25.39	23.46
	Cal c		39.49	39.49	43.46	45.99	43.46
	Cal d		83.31	83.31	86.79	93.61	86.79
	Cal e		151.75	151.75	153.16	161.95	153.16