

C-reactive Protein (CRP) Calibrator

Value sheet of Mindray Measurement System



<div><div>LOT: 158323005</div><div> 2024-08-25</div></div>							
Parameter	Level	Unit	Assay Value				
			BC-5390 CRP/ BC-5180 CRP/ BC-5310 CRP/ BC-5100 CRP	CRP-M100	BC-7500[B] CRP/ BC-7500[N] CRP/ BC-7500[R] CRP/ BC-7500[NR] CRP/ BC-7300[B] CRP/ BC-7300[N] CRP/ BC-7300[R] CRP/ BC-7300[NR] CR / BC-7500[B] CS/ BC-7500[N] CS/ BC-7500[R] CS/ BC-7500[NR] CS	BC-700[B] CS/ BC-700[R] CS/ BC-760[B] CS/ BC-760[R] CS	BP 200n/BP 260n
C-reactive Protein (CRP)	Cal a	mg/L	5.77	5.77	7.39	7.54	7.39
	Cal b		20.64	20.64	23.18	24.78	23.18
	Cal c		39.81	39.81	42.57	44.74	42.57
	Cal d		84.39	84.39	85.24	90.06	85.24
	Cal e		155.75	155.75	148.94	157.76	148.94