

C-reactive Protein (CRP) Calibrator

Value sheet of Mindray Measurement System



LOT: 158323008							
2024-11-22							
Parameter	Level	Unit	Assay Value				
			BC-5390 CRP/ BC-5180 CRP/ BC-5310 CRP/ BC-5100 CRP	CRP-M100	BC-7500[B]CRP/ BC-7500[N] CRP/ BC-7500[R] CRP/ BC-7500[NR] CRP/ BC-7300[B] CRP/ BC-7300[N] CRP/ BC-7300[R] CRP/ BC-7300[NR] CR / BC-7500[B] CS/ BC-7500[N]CS/ BC-7500[R]CS/ BC-7500[NR] CS	BC-700[B] CS/ BC-700[R] CS/ BC-760[B] CS/ BC-760[R] CS	BP 200n/BP 260n
C-reactive Protein (CRP)	Cal a	mg/L	5.72	5.72	7.30	7.71	7.30
	Cal b		20.56	20.56	23.70	26.62	23.70
	Cal c		38.93	38.93	41.54	45.79	41.54
	Cal d		83.90	83.90	83.43	92.95	83.43
	Cal e		151.20	151.20	147.77	159.45	147.77