

C-reactive Protein (CRP) Calibrator

Value sheet of Mindray Measurement System



<div><div>LOT: 158324001</div><div> 2025-01-22</div></div>							
Parameter	Level	Unit	Assay Value				
			BC-5390 CRP/ BC-5180 CRP/ BC-5310 CRP/ BC-5100 CRP	CRP-M100	BC-7500[B] CRP/ BC-7500[N] CRP/ BC-7500[R] CRP/ BC-7500[NR] CRP/ BC-7300[B] CRP/ BC-7300[N] CRP/ BC-7300[R] CRP/ BC-7300[NR] CR / BC-7500[B] CS/ BC-7500[N] CS/ BC-7500[R] CS/ BC-7500[NR] CS	BC-700[B] CS/ BC-700[R] CS/ BC-760[B] CS/ BC-760[R] CS	BP 200n/BP 260n
C-reactive Protein (CRP)	Cal a	mg/L	5.81	5.81	7.42	7.67	7.42
	Cal b		20.78	20.78	23.90	25.95	23.90
	Cal c		39.40	39.40	42.17	45.21	42.17
	Cal d		85.29	85.29	84.60	93.17	84.60
	Cal e		153.08	153.08	146.53	162.32	146.53