

# C-reaction Protein Control



## Value sheet of Mindray Measurement System



| <b>LOT</b> : 061423002 <b>EXP</b> : 2024-02-17 |      |  |                     |             |               |
|--|------|--|---------------------|-------------|---------------|
| Parameter                                      | Unit | Analyzer   | Concentration Level | Assay Value | Range         |
| C-reaction Protein                             | mg/L | BC-5390 CRP,BC-5180 CRP,<br>BC-5310 CRP, BC-5100 CRP   | I                   | 2.70        | 1.89 - 3.51   |
|  |      |  | II                  | 43.72       | 30.60 - 56.84 |
|  |      | BC-7500[B] CRP,BC-7500[N] CRP,<br>BC-7500[R] CRP,BC-7500[NR] CRP,<br>BC-7500[B] CS,BC-7500[N] CS,<br>BC-7500[R] CS,BC-7500[NR] CS,<br>BC-7300[B] CRP,BC-7300[N] CRP,<br>BC-7300[R] CRP,BC-7300[NR] CRP | I                   | 3.54        | 2.48 - 4.60   |
|  |      |  | II                  | 50.45       | 35.31 - 65.59 |
|  |      | CRP-M100   | I                   | 2.70        | 1.89 - 3.51   |
|  |      |  | II                  | 43.72       | 30.60 - 56.84 |
|  |      | BC-700[B] CS<br>BC-760[B] CS<br>BC-700[R] CS<br>BC-760[R] CS   | I                   | 3.50        | 2.45 - 4.55   |
|  |      |  | II                  | 50.74       | 35.52 - 65.96 |
|  |      | BP 200n<br>BP 260n   | I                   | /           | / - /         |
|  |      |  | II                  | /           | / - /         |