


mindray

BeneHeart C1A/C2

AED Trainer

Smarter & Faster



C1A/C2 AED Trainer

C1A/C2 AED Trainer is designed specific for training first-aid personnel before they use BeneHeart C1A/C2 AED. The trainer fully simulates the BeneHeart C1A/C2 AED in the training guide and drill program, with a view to help first-aid personnel to get familiar with the AED and grasp required basic operation skills when using it in the first-aid process.



High Level Simulation

- Appearance, interface and operation process identical to the 'real' C-series
- Preset 6 basic scenes and 4 simulation scenarios to better meet training needs
- Meet the latest AHA/ERC guidelines
- Support semi-auto and fully auto modes

Easy to Maintain

- The AED trainer is equipped with high quality training pads. These high quality pads will last for multiple uses
- The trainer uses standard AA batteries, which can be purchased conveniently

Specifications

AED Trainer	
Size (W x D x H)	21.0 cm x 28.6 cm x 7.8 cm
Weight	1.4 kg (C1A AED trainer, without screen, with batteries) 1.7 kg (C2 AED trainer, with screen, with batteries)
Batteries *	DC 12V, 8 AA batteries
Display	7 inch (for C2 AED trainer only)
Buttons	Language selection (up to 26 languages), adult/child mode switch, shock button
CPR metronome	30:2, 15:2
Remote Control	
Batteries *	DC 3V, 2 AAA batteries
Standard scenario	6
Pads connected	Yes, no
Shock advised	Yes, no
Discharge mode	Semi-auto, Fully auto

* Note: The batteries are not included in the trainer & acc. kits.

www.mindray.com

P/N:ENG-C1A/C2 AED Trainer-210285X2P-20200221

©2020 Shenzhen Mindray Bio-Medical Electronics Co.,Ltd. All rights reserved.

mindray
healthcare within reach