

The materials and information provided via this brochure are distributed internationally, but not all the products or services found on this brochure are necessarily available in your country or area. Please contact local office or agents to make sure of the availability of the relevant products or services.

**Mindray Medical Netherlands B.V.**  
Drs.W.van Royenstraat 8, 3871 AN Hoevelaken, The Netherlands  
Voice: (31-33) 254-4911  
Fax: (31-33) 253-4280

**Mindray Medical Italy S.r.l.**  
Via Leonardo DA Vinci, 158, 20090 Trezzano Sul Naviglio (MI), Italy  
Voice: (39-02) 5737-401  
Fax: (39-02) 5523-2018

**Mindray Medical España S.L.**  
Avenida Manoteras 38, Edificio B, Bajo 8008. 28050-Madrid Spain  
Voice: (34-91) 392-3754  
Fax: (34-91) 392-3756

**Mindray Medical Germany GmbH**  
Goebelstrasse 21, 3. OG, 64293 Darmstadt  
Voice: (49 6151) 3910 0  
Fax: (49 6151) 3910 300

**Mindray (UK) Limited**  
Mindray House, Kingfisher Way, Hinchingsbrooke Business Park, Huntingdon, Cambs, PE29 6FN  
Voice: (44) 1480-416-840  
Fax: (44) 1480-436-588

**Mindray Medical France SARL**  
Europarc Creteil 1 allée des cerisiers 94000 CRETEIL, France  
Voice: (33-1) 4513-9150  
Fax: (33-1) 4513-9151

**Mindray Medical Rus Co.Ltd**  
2 Zvenigorodskaya Street, 13, bld 41, 123022 Moscow, Russia  
Voice: (7-499) 553-6036  
Fax: (7-499) 553-6039

**Mindray Medical Colombia S.A.S**  
Room No.702, Torre UNIKA, Carrera 9 No.77-67, Bogota, Colombia  
Voice: (57-1) 3130892 / 3210916

**Mindray do Brasil – Comércio e Distribuição de Equipamentos Médicos Ltda.**  
Rua Tavares Bastos, 329, Perdizes, São Paulo - SP, CEP 05012-020  
Voice: (55-11) 3124-8026  
Fax: (55-11) 3078-8035

**Mindray Medical Mexico S. de R.L. de C.V.**  
Félix Parra # 175 Colonia. San José Insurgentes Sur Delegación Benito Juárez 03900 México, D. F.  
Voice: (52-55) 5661-9450 / (52-55) 5662-6620  
Fax: (52-55) 5662-6597

**PT.Mindray Medical Indonesia**  
The East 11th Floor, Unit 08 Jl. Dr. Ide Anak Agung Gde Agung Kav. E 3.2 No.1 Jakarta 12950  
Voice: (6221) 2902-7280  
Fax: (6221) 2902-7283

**Mindray Medical Egypt**  
Ground Floor -11 El Imam Metwally El Sharawi Street, Sheraton Heliopolis - El Nozha, Cairo, Egypt  
Voice: (2) 02 22691047  
Fax: (2) 02 22691046

**Mindray Medical Australia**  
ADD: 19 Harker St, Burwood VIC 3125, Australia  
Voice: (61-3) 9830 6610  
Fax: (61-2) 9550 2841

**Mindray Medical Nigeria**  
ADD: 1004 Estate,Victoria Island, Lagos, Nigeria  
Voice: (+234) 7030682053

**Mindray Medical India Pvt. Ltd. (Delhi Office)**  
16th floor, Building No. 9-B, DLF Cyber City, DLF Phase- 3, Gurgaon- 122002, Haryana  
Voice: (91-124) 4632-488  
Fax: (91-124) 4632499

**Mindray Medical India Pvt. Ltd. (Mumbai office)**  
B-404, City Point, Andheri Kurla Road, Andheri East, Mumbai, 400 059, India  
Voice: (91-22) 4020-0000  
Fax: (91-22) 4020-0011


**Mindray Medikal Teknoloji Istanbul Ltd. Sti.**  
KÖYALTI MEVKİLİ 29 EKİM CAD.N.1 KUYUMCUKENT 2.NO.PLZ.3.KT. D.1- 9 34520 BAĞÇELİEVLER/YENİBOSNA, İSTANBUL, TÜRKİYE  
Voice: (90-212) 482 08 77  
Fax: (90-212) 482 08 78

**Mindray Medikal Teknoloji Istanbul Ltd. Sti. Branch Office Ankara**  
TURAN GÜNEŞ BULVARI ALEKSANDER DÜBÇEK CAD. NO:18/A ÇANKAYA/ANKARA, TÜRKİYE  
Voice: (90-312) 441 96 92  
Fax: (90-312) 441 96 93

**Mindray Medical Vietnam Co.,Ltd.**  
Unit 2, 4th Floor, "B" Building, Broadway Office Park, 102 Nguyen Luong Bang Blvd., Tan Phu Ward, Dist. 7, HCMC, Vietnam  
Voice: (84) 8 5413-6970  
Fax: (84) 8 5413-6971

**Mindray Medical Thailand Limited**  
1768 Thai summit Tower Room 407, 4th floor, New petchaburee Road, Kwaeng Bangkokpi, Khet Huai Khwang, Bangkok 10320, Thailand  
Voice: (66) 2652-7288  
Fax: (66) 2652-7266

**Mindray Medical Korea**  
KNI Tower 9F 823-19 Yeoksam dong, Seoul-ganum, Korea  
Tel: 008225688040  
Fax: 008225688042



# BeneHeartD3

## Defibrillator /Monitor

Technical Specifications	
<b>Physical</b>	
Dimensions:	288mm(w)x203mm(d)x275mm(h)
<b>Weight :</b>	
Main unit:	4.7 kg
Battery package (each):	0.54kg
External paddle set:	0.86 kg
<b>Environmental and Physical Requirements</b>	
Water Resistance:	IPX4 (without external power)
Solids Resistance:	IP4X
Temperature:	Operating: 0 to 45 °C Storage: -30 to 70 °C
Humidity:	Operating / Storage: 10 to 95%.non-condensation
Altitude:	Operating / Storage: -381m to +4575 m
Shock and Vibration:	Meets the requirements of 21.102, ISO9919 (Shock and vibration for transport)
Bump:	Meets the requirements of 6.3.4.2, EN1789 (Medical devices for use in road ambulances)
Free fall:	Meets the requirements of 6.3.4.3, EN1789 (Height of fall: 0.75 m)
EMC:	Meets IEC60601-1-2
Safety:	Meets EN/IEC 60601-1
<b>Display</b>	
Type:	TFT Color LCD
Dimensions:	7 inch
Resolution:	800x480 pixels
Display Waveforms:	Max. 3 channels
Wave Viewing Time:	Max. 16 s (ECG)
<b>Power</b>	
<b>AC Power</b>	
Line voltage:	100 to 240 V~ (±10%)
Current:	1.8 to 0.8 A
Frequency:	50/60 Hz (±3 Hz)
<b>DC Power (through DC-AC Inverter)</b>	
Input voltage:	12 VDC
Power Consumption:	190 W
<b>Battery</b>	
Type:	15.1V, 5600mAh, rechargeable lithium ion battery pack
Number:	1
Charge time:	Less than 3 hours to 90% and less than 4 hours to 100% with equipment power of
Capacity Indicator:	5-segment led indicator for fast battery capacity evaluation
Capacity (new, fully charged battery):	Monitoring Mode: 6 hours, Monitoring with a 5-lead ECG,Resp, SpO <sub>2</sub> , CO <sub>2</sub> and NIBP measurements set at an interval of 15 minutes.
	Wi-Fi is disabled

Mindray Building, Keji 12th Road South,  
High-tech Industrial Park, Nanshan, Shenzhen 518057, P.R. China  
Tel: +86 755 8188 8998 Fax: +86 755 26582680  
E-mail: intl-market@mindray.com www.mindray.com

**mindray**

**mindray** (healthcare with us) are registered trademarks or trademarks owned by Shenzhen Mindray Bio-medical Electronics Co., LTD.  
 ©2015 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved. Specifications subject to changes without prior notice.  
 P/N:ENG-BeneHeart D3 datasheet-210285x4P-20170808

	Defib Mode: 200 times, 360J discharge at intervals of 1minute without recording
	Pacing Mode: 4.5 hours, 50 Ohm load impedance, Pacing rate: 80bpm,Pacing output: 60mA
<b>Recorder</b>	
Method:	High-resolution thermal dot array
Waveforms:	Max. 3 channels
Speed:	6.25mm/s, 12.5mm/s, 25mm/s, 50mm/s
Paper width:	50 mm
Reports:	Real time waveforms,Event Summary, Tabular Trends, Frozen Waveforms, Review , User test, and Configuration
Auto Recording:	Recorder can be configured to record marked events, charge, shock , Alarm, auto test
<b>Data Storage</b>	
Patient profiles:	Max. 100 patients
Events:	Up to 1000 events for one patient
Waveform Storage:	Up to 24 hours of consecutive ECG waveform
Tabular Trends:	72 hours, resolution: 1 min.
Voice recording:	Max. 180 min in total; max. 60 min for each patient
Data Export:	Data can be exported to PC through USB flash memory
<b>Defibrillator</b>	
Waveform:	Biphasic truncated exponential waveform , with impedance compensation
Energy accuracy:	±2 J or 15% of setting, whichever is greater, into 50 Ohm
Power on time:	Less than 2 seconds with a new, fully charged battery
Charge Time:	Less than 3 seconds to 200 Joules with a new, fully charged battery
	Less than 7 seconds to 360 Joules with a new, fully charged battery
ECG recovery time:	Less than 2.5 seconds
Shock Delivery:	Via multifunction defib electrode pads, or paddles
Patient Impedance Range:	25 to 300 Ω (External defibrillation)
<b>Manual Mode</b>	
Output Energy:	1,2,3,4,5,6,7,8,9,10,15,20,30,50,70,100,150,170,200,300,360 J
Synchronous Cardioversion:	Energy transfer begins within 60ms of the QRS peak
	Energy transfer begins within 25ms of the External Sync Pulse
<b>AED Mode</b>	
Output Energy:	User configurable
AED Shock Series:	Energy level: 100 to 360 J, configurable
	Shocks series: 1, 2, 3, configurable
	Default configuration meets 2015 AHA Guidelines
	CPRmode with 1-channel ECG monitoring
Sensitivity and Specificity:	Meets AAMI DF-80
<b>Noninvasive Pacing</b>	
Waveform:	Monophasic square wave pulse
Pulse Width:	20 ms or 40 ms, ±5%
Refractory period:	200 to 300 ms, ±3% (function of rate)
Pacing Mode:	Demand or fixed
Pacing rate:	30 ppm to 210 ppm, ±1.5%
Pacing output:	0 mA to 200 mA, ±5% or 5mA, whichever is greater
4:1 pacing:	Pacing pulse frequency reduced by factor of 4 when activated
<b>ECG Monitoring</b>	
Lead type:	3 leads ECG, 5 leads ECG
Lead Selection:	3 leads ECG: I, II, III; 5 leads ECG: I, II, III, aVR, aVL, aVF, V
Heart Rate Display:	Adult: 15 to 300 bpm
	Pediatric: 15 to 350 bpm
	Neonate: 15 to 350 bpm
Resolution:	1 bpm

Arrhythmia :	Yes
Alarms:	Yes
ECG size:	2.5 mm/mV(×0.25), 5 mm/mV(×0.5), 10 mm/mV(×1), 20 mm/mV(×2), 40 mm/mV(×4), Auto
Sweep speed:	6.25 mm/s, 12.5 mm/s, 25 mm/s, 50 mm/s
Patient Isolation (defibrillation proof):	
	Type CF: ECG, RESP, SpO <sub>2</sub> , NIBP
	Type BF: CO <sub>2</sub>
<b>Respiration</b>	
Method:	Trans-thoracic impedance
Range:	Adult: 0 to 200 rpm
	Pediatric, neonate: 0 to 200 rpm
Resolution:	1 rpm
<b>SpO<sub>2</sub>Pulse Oximetry</b>	
<b>Mindray SpO<sub>2</sub></b>	
Range:	0 to 100%
Resolution:	1%
PR Range:	20 to 300 bpm
<b>Nellcor SpO<sub>2</sub></b>	
Range:	1 to 100%
Resolution:	1%
PR Range:	20 to 300 bpm.
<b>NIBP</b>	
Operating mode:	Manual, Auto, STAT
Static pressure range:	0 to 300 mmHg
Displayed Pressures:	Systolic, Diastolic, Mean
Cuff inflation pressure(Default):Adult:	160±5 mmHg
	Pediatric: 140±5 mmHg
	Neonate: 90 ± 5 mmHg
<b>CO<sub>2</sub></b>	
Measurement range:	0 to 150 mmHg
Resolution:	1 mmHg
awRR measurement range:	0 to 150 rpm
awRR accuracy:	0<60 rpm: ±1 rpm
	60 to 150 rpm: ±2 rpm