

## Product Comparison


		A7 Advantage	A5 Advantage	A4 Advantage
<b>Hardware</b>	Auxiliary power outlets	4	4	3
	Central locking brake	✓	✓	
	Integrated cable sweep system	✓	✓	
	Backup touchscreen control	Touchpad	Touchpad	Optional Mouse
	Backup Cylinder Yokes	O <sub>2</sub> , Air, N <sub>2</sub> O	O <sub>2</sub> , Air, N <sub>2</sub> O	O <sub>2</sub> , Air
	Integrated gas analysis module	Optional	Optional	Optional
<b>Ventilation Modes</b>	Manual/Spontaneous	✓	✓	✓
	VCV	✓	✓	✓
	CPAP/PS	✓	✓	✓
	PCV	✓	✓	Optional
	PCV – Volume Guarantee (PCV-VG)	✓	✓	Optional
	SIMV-VC	✓	✓	Optional
	SIMV-PC	✓	✓	Optional
	SIMV-VG	Optional	Optional	Optional
	APRV	Optional	Optional	
Lung Recruitment	Optional	Optional		
<b>Advanced Features</b>	Optimizer low flow support*	Optional	Optional	Optional
	Spirometry loops display	Optional	Optional	
	Cardiac bypass mode	Optional	Optional	
	Gas analysis with agent ID, age MAC	Optional	Optional	Optional
	High pressure O <sub>2</sub> port	✓	✓	
	Electronic fresh gas flow control	✓		
	Agent usage calculation	Optional	Optional	Optional
	Quick release APL valve	✓	Optional	Optional
	Integrated suction	✓		
	Electronically controlled Aux Common Gas Outlet	✓		
	Pull out work table	✓		
	<b>Safety Features</b>	Auxiliary O <sub>2</sub> /Air supply	✓	✓
Lithium-ion battery backup		90 min. (2)	150 min. (2)	150 min. (2)
Fresh gas with agent during power outage		O <sub>2</sub> /Air	O <sub>2</sub> /N <sub>2</sub> O/Air	O <sub>2</sub> /N <sub>2</sub> O/Air
<b>Optional Accessories</b>	Flexible bag arm	✓	✓	✓
	Three (in-line) vaporizer mount	✓	✓	✓
	Utility tray (side mount)	✓	✓	✓
	Various mounting arms	✓	✓	✓
<b>Cost Saving Features</b>	Non-proprietary vaporizer	✓	✓	✓
	Non-proprietary sodalime cannister	✓	✓	✓
	Only one PM required annually	✓	✓	✓

\*Optimizer functionality requires Gas Module (AGM).

# A Series ADVANTAGE

ADVANCED ANESTHESIA SYSTEMS

## Standard Components



A7 Advantage model shown. Not all features listed.

Mindray DS USA, Inc., 800 MacArthur Blvd., Mahwah, NJ 07430  
Tel: 800.288.2121 Tel: 201.995.8000 Fax: 800.926.4275 [www.mindray.com](http://www.mindray.com)

©2018 Mindray DS USA, Inc. Mindray® is a registered trademark of Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved.  
10/18 P/N: 0002-08-1597 Rev A

**mindray**  
NORTH AMERICA