

TE X-Series Ultrasound Systems

Flex Editions

Smart Fluid Management

Smart Nerve

Auto Gastric Antrum

Flexibility Redefined

The TE X Flex Edition Ultrasound Systems are scaled to address the diverse and evolving needs of today's healthcare landscape. Offering a blend of flexibility and cutting-edge technology, the TE X Flex Acute Care and TE X Flex Procedural Editions offer clinicians a versatile tool tailored to a wide range of clinical scenarios, from routine examinations to more complex diagnostic challenges.

AI Innovations

The Mindray TE X Flex Editions harness the transformative potential of artificial intelligence (AI) to redefine Point of Care imaging. Discover their cutting-edge capabilities, including:

- Smart Fluid Management
- Smart Nerve
- Auto Gastric Antrum




TE X-Series Comparison


Introducing the newest additions to the TE X-Series: the Acute Care and Procedural Editions. Tailored for emergency, critical care, anesthesia, and interventional radiology, these new editions provide specialized and adaptable solutions for key medical fields. Building upon the TE X foundation, these systems feature streamlined designs and user-centric interfaces to elevate patient care within specialized clinical settings.

	TE X Flex Acute Care	TE X Flex Procedural	TE X
System			
23.8" Rotatable HD Touchscreen	✓	✓	✓
Lithium Ion Batteries	4	4	6
4 Active Transducer Ports	✓	✓	✓
B-Mode, Color, PW, M-Mode, PD	✓	✓	✓
TDI	✓	N/A	✓
CW	✓	N/A	✓
Features			
Fluid Management Package	✓	N/A	✓
X-Pilot Application-based Workflow	✓	N/A	✓
Smart EchoVue	N/A	N/A	✓
Advance Cardiac Suite	N/A	N/A	✓
Smart Nerve	N/A	✓	✓
iNeedle ⁺	✓	✓	✓
Auto Gastric Antrum (Auto GA)	N/A	✓	✓
Smart Fetal Heart Rate (Smart FHR)	✓	N/A	✓
iWorks™ Protocols	✓	Optional	✓


Clinical Gallery



2D PSAX View




Pleural Effusion




4C Tricuspid Regurgitation



Ascites



Quadratus Lumborum



Brachial Plexus with Smart Nerve

Mindray North America

800 MacArthur Blvd., Mahwah, NJ 07430

Tel: 800.288.2121 Support: 877.913.9663 Fax: 800.926.4275 www.mindray.com

Mindray® is a registered trademark of Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All brands and product names are trademarks of their respective owners. ©2024 Mindray DS USA, Inc. Subject to change. 03/24 P/N: 0002-08-40761 Rev 1.0

mindray

healthcare within reach