

M9 Premium

Ultrasound System




Premium Capability, Maximum Mobility


Mindray's M9 Premium Ultrasound System provides a powerful and innovative platform that exceeds the demands of today's challenging healthcare environment. Powered by Mindray's revolutionary mQuadro+ system architecture and cutting-edge 3T™ transducer technology, the M9 Premium Ultrasound System provides market-leading performance and cutting-edge software features to ensure clinicians have access to the tools they need to provide excellent patient care and elevate outcomes.

Enhanced Image Clarity


With new single crystal 3T transducers and the powerful mQuadro beamformer platform, premium image quality has never been more accessible on a wide variety of patient exams.




Morrison's Pouch




Gallbladder with Stones




Carotid Artery with Pulse Wave Doppler



Kidney Perfusion




Parasternal Long




4 Chamber with Color

Advanced Technologies


The M9 Premium offers a full suite of market-leading cardiac analysis and Artificial Intelligence (AI)-driven Smart Tools to help expedite clinical decision making and achieve reproducibility between exams.




Smart B-Line
Automated B-Line Analysis




Smart IVC
Automated IVC Measurement




Smart VTI
Automated VTI and Cardiac Output




TTQA
Myocardial Strain Analysis




RV Strain
Right Heart Strain Analysis




Auto DFR
Automated Diastolic Function Ratio



Auto EF
Automated Ejection Fraction



TEE
Transesophageal Echo



Auto GA
Automated Gastric Antrum Area

Innovative Design

The newly designed cart on the M9 Premium System is inspired by user feedback and designed to protect the system, enabling convenient transport and storage.

15.6" HD display with 1920 x 1080 resolution LED display

Two storage baskets provide ample storage

Small footprint for easy maneuvering in tight spaces



Three active transducer ports



Adjustable height ranging from 32" - 42"



Additional battery pack on cart for up to 3 hours of continuous scanning

Impressive Value

The Industry's Best Investment Protection

Living Technology

Living Technology™ is a continually evolving software-based approach to providing our customers with easily upgradeable enhancements made possible by our core imaging technologies. These upgrades secure product investment protection by ensuring that Mindray Ultrasound Systems remain at the cutting-edge of imaging performance throughout their entire life cycle.

Standard warranty* includes:

Five-years of coverage for:

- Main system and parts, and technical phone support Monday – Friday, 8 am – 8 pm ET (excluding U.S. holidays)
- Standard transducers (normal wear and failure)
- Software performance updates with Living Technology

One-year of coverage for:

- Specialty transducers including, but not limited to, TEE and 4D transducers

* Standard five (5) year warranty is available for systems purchased directly from Mindray North America. Warranties for systems purchased from Mindray authorized agents can vary.

Mindray North America

800 MacArthur Blvd., Mahwah, NJ 07430

Tel: 800.288.2121 Support: 877.913.9663 Fax: 800.926.4275 www.mindray.com

Mindray® is a registered trademark of Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All brands and product names are trademarks of their respective owners. ©2022 Mindray DS USA, Inc. Subject to change. 11/22 P/N: 0002-08-40617 Rev B

mindray

healthcare within reach