

Metabolic Multi Control (H)



Target Value Sheet

Product	Lot Number and Expiration Date	Analytes	Instrument Model	Unit	Target Values	Target Values and Ranges (Target Values±3SD)			ISD	International Unit	Target Values	Target Values and Ranges (Target Values±3SD)			ISD
Metabolic Multi Control (H)	LOT 2025060101	VD-T	CL-900i Series	ng/mL	69.86	41.93	~	97.79	9.31	nmol/L	174.65	104.83	~	244.48	23.28
			CL-8000i Series		69.86	41.93	~	97.79	9.31		174.65	104.83	~	244.48	23.28
			CL-2000i/CL-1000i/CL-6000i Series		69.86	41.93	~	97.79	9.31		174.65	104.83	~	244.48	23.28
		CT	CL-900i Series	pg/mL	158.40	95.04	~	221.76	21.12	pmol/L	46.35	27.81	~	64.89	6.18
			CL-8000i Series		158.40	95.04	~	221.76	21.12		46.35	27.81	~	64.89	6.18
			CL-2000i/CL-1000i/CL-6000i Series		158.40	95.04	~	221.76	21.12		46.35	27.81	~	64.89	6.18
		Folate	CL-900i Series	ng/mL	13.36	8.02	~	18.70	1.78	nmol/L	30.27	18.17	~	42.37	4.03
			CL-8000i Series		13.36	8.02	~	18.70	1.78		30.27	18.17	~	42.37	4.03
			CL-2000i/CL-1000i/CL-6000i Series		13.36	8.02	~	18.70	1.78		30.27	18.17	~	42.37	4.03
		VB12	CL-900i Series	pg/mL	931.48	558.88	~	1304.08	124.20	pmol/L	687.43	412.45	~	962.41	91.66
			CL-8000i Series		931.48	558.88	~	1304.08	124.20		687.43	412.45	~	962.41	91.66
			CL-2000i/CL-1000i/CL-6000i Series		931.48	558.88	~	1304.08	124.20		687.43	412.45	~	962.41	91.66
		PTH	CL-900i Series	pg/mL	808.03	484.81	~	1131.25	107.74	pmol/L	85.65	51.39	~	119.91	11.42
			CL-8000i Series		808.03	484.81	~	1131.25	107.74		85.65	51.39	~	119.91	11.42
			CL-2000i/CL-1000i/CL-6000i Series		808.03	484.81	~	1131.25	107.74		85.65	51.39	~	119.91	11.42
		FERR	CL-900i Series	ng/mL	306.34	183.79	~	428.89	40.85	µg/L	306.34	183.79	~	428.89	40.85
			CL-8000i Series		306.34	183.79	~	428.89	40.85		306.34	183.79	~	428.89	40.85
			CL-2000i/CL-1000i/CL-6000i Series		306.34	183.79	~	428.89	40.85		306.34	183.79	~	428.89	40.85
		VD-T II	CL-900i Series	ng/mL	57.82	35.26	~	80.38	7.52	nmol/L	144.55	88.15	~	200.95	18.80
			CL-8000i Series		57.82	35.26	~	80.38	7.52		144.55	88.15	~	200.95	18.80
			CL-2000i/CL-1000i/CL-6000i Series		57.82	35.26	~	80.38	7.52		144.55	88.15	~	200.95	18.80
	2026/12/3	/	/	/	/	/	~	/	/	/	/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
		/	/	/	/	/	~	/	/	/	/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
		/	/	/	/	/	~	/	/	/	/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
		/	/	/	/	/	~	/	/	/	/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
		/	/	/	/	/	~	/	/	/	/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
		/	/	/	/	/	~	/	/	/	/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
		/	/	/	/	/	~	/	/	/	/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
		/	/	/	/	/	~	/	/	/	/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/

Remark: