

AMH Control



Target Value Sheet

Product	Lot Number and Expiration Date	Analytes	Instrument Model	Unit	Target Values	Target Values and Ranges (Target Values±3SD)	1SD	International Unit	Target Values	Target Values and Ranges (Target Values±3SD)	1SD
AMH Control (L)	<div><div>LOT</div><div>2022110111</div><div></div><div>2024-06-25</div></div>	AMH	Mindray CL-series Chemiluminescence Immunoassay Analyzer	ng/mL	1.00	0.61 ~ 1.39	0.13	pmoL/L	7.14	4.36 ~ 9.92	0.93