

# AMH Control



## Target Value Sheet

| Product         | Lot Number and Expiration Date   | Analytes | Instrument Model   | Unit  | Target Values | Target Values and Ranges (Target Values±3SD) | 1SD  | International Unit | Target Values | Target Values and Ranges (Target Values±3SD) | 1SD  |
|-----------------|--|----------|--|-------|---------------|--|------|--------------------|---------------|--|------|
| AMH Control (H) | <div>LOT</div> <div>2022110111</div> <div></div> <div>2024-06-25</div> | AMH      | Mindray CL-series Chemiluminescence Immunoassay Analyzer | ng/mL | 4.69          | 2.80 ~ 6.58                                  | 0.63 | pmoL/L             | 33.49         | 19.99 ~ 46.98                                | 4.50 |