

Thyroid Function Multi Control (H)

Target Value Sheet 1

Product	Lot Number and Expiration Date	Analytes	Instrument Model	Unit	Target Values	Target Values and Ranges (Target Values±3SD)			ISD	International Unit	Target Values	Target Values and Ranges (Target Values±3SD)			ISD
Thyroid Function Multi Control (H)	LOT 2025040211	FT3	CL-900i Series	pg/mL	8.46	5.07	~	11.85	1.13	pmol/L	12.99	7.79	~	18.20	1.74
			CL-8000i Series		8.46	5.07	~	11.85	1.13		12.99	7.79	~	18.20	1.74
			CL-2000i/CL-1000i/CL-6000i Series		8.46	5.07	~	11.85	1.13		12.99	7.79	~	18.20	1.74
		FT4	CL-900i Series	ng/dL	3.68	2.21	~	5.15	0.49	pmol/L	47.36	28.44	~	66.28	6.31
			CL-8000i Series		3.68	2.21	~	5.15	0.49		47.36	28.44	~	66.28	6.31
			CL-2000i/CL-1000i/CL-6000i Series		3.68	2.21	~	5.15	0.49		47.36	28.44	~	66.28	6.31
		T3	CL-900i Series	ng/mL	2.23	1.33	~	3.13	0.30	nmol/L	3.43	2.04	~	4.81	0.46
			CL-8000i Series		2.23	1.33	~	3.13	0.30		3.43	2.04	~	4.81	0.46
			CL-2000i/CL-1000i/CL-6000i Series		2.23	1.33	~	3.13	0.30		3.43	2.04	~	4.81	0.46
		T4	CL-900i Series	µg/dL	14.97	8.97	~	20.97	2.00	nmol/L	192.66	115.44	~	269.88	25.74
			CL-8000i Series		14.97	8.97	~	20.97	2.00		192.66	115.44	~	269.88	25.74
			CL-2000i/CL-1000i/CL-6000i Series		14.97	8.97	~	20.97	2.00		192.66	115.44	~	269.88	25.74
	2026/4/28	TSH	CL-900i Series	µIU/mL	24.893	14.936	~	34.850	3.319	mIU/L	24.893	14.936	~	34.850	3.319
			CL-8000i Series		24.893	14.936	~	34.850	3.319		24.893	14.936	~	34.850	3.319
			CL-2000i/CL-1000i/CL-6000i Series		24.893	14.936	~	34.850	3.319		24.893	14.936	~	34.850	3.319
		Tg	CL-900i Series	ng/mL	271.11	162.66	~	379.56	36.15	µg/L	271.11	162.66	~	379.56	36.15
			CL-8000i Series		271.11	162.66	~	379.56	36.15		271.11	162.66	~	379.56	36.15
			CL-2000i/CL-1000i/CL-6000i Series		271.11	162.66	~	379.56	36.15		271.11	162.66	~	379.56	36.15
		/	/	/	/	/	~	/	/	/	/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
		/	/	/	/	/	~	/	/	/	/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
		/	/	/	/	/	~	/	/	/	/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
		/	/	/	/	/	~	/	/	/	/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
		/	/	/	/	/	~	/	/	/	/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
		/	/	/	/	/	~	/	/	/	/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
		/	/	/	/	/	~	/	/	/	/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/

Remark:

The target value and range of
The target value and range of
The target value and range of
The target value and range of

FT4
T4
TSH
Tg

can only be used for reagent after lot#
can only be used for reagent after lot#
can only be used for reagent after lot#
can only be used for reagent after lot#

20230109
20230109
20230809
20240309

(excluding lot#20230109)
(excluding lot#20230109)
(excluding lot#20230809)
(excluding lot#20240309)

Thyroid Function Multi Control (H)

Target Value Sheet 2

Product	Lot Number and Expiration Date	Analytes	Instrument Model	Unit	Target Values	Target Values and Ranges (Target Values±3SD)			ISD	International Unit	Target Values	Target Values and Ranges (Target Values±3SD)			ISD
Thyroid Function Multi Control (H)	LOT 2025040211	FT3	CL-900i Series	pg/mL	8.46	5.07	~	11.85	1.13	pmol/L	12.99	7.79	~	18.20	1.74
			CL-8000i Series		8.46	5.07	~	11.85	1.13		12.99	7.79	~	18.20	1.74
			CL-2000i/CL-1000i/CL-6000i Series		8.46	5.07	~	11.85	1.13		12.99	7.79	~	18.20	1.74
		FT4	CL-900i Series	ng/dL	2.79	1.68	~	3.90	0.37	pmol/L	35.91	21.62	~	50.19	4.76
			CL-8000i Series		2.79	1.68	~	3.90	0.37		35.91	21.62	~	50.19	4.76
			CL-2000i/CL-1000i/CL-6000i Series		2.79	1.68	~	3.90	0.37		35.91	21.62	~	50.19	4.76
		T3	CL-900i Series	ng/mL	2.23	1.33	~	3.13	0.30	nmol/L	3.43	2.04	~	4.81	0.46
			CL-8000i Series		2.23	1.33	~	3.13	0.30		3.43	2.04	~	4.81	0.46
			CL-2000i/CL-1000i/CL-6000i Series		2.23	1.33	~	3.13	0.30		3.43	2.04	~	4.81	0.46
		T4	CL-900i Series	µg/dL	15.32	9.20	~	21.44	2.04	nmol/L	197.17	118.40	~	275.93	26.25
			CL-8000i Series		15.32	9.20	~	21.44	2.04		197.17	118.40	~	275.93	26.25
			CL-2000i/CL-1000i/CL-6000i Series		15.32	9.20	~	21.44	2.04		197.17	118.40	~	275.93	26.25
		TSH	CL-900i Series	µIU/mL	34.681	20.809	~	48.553	4.624	mIU/L	34.681	20.809	~	48.553	4.624
			CL-8000i Series		32.118	19.272	~	44.964	4.282		32.118	19.272	~	44.964	4.282
			CL-2000i/CL-1000i/CL-6000i Series		34.681	20.809	~	48.553	4.624		34.681	20.809	~	48.553	4.624
		Tg	CL-900i Series	ng/mL	27.12	16.26	~	37.98	3.62	µg/L	27.12	16.26	~	37.98	3.62
			CL-8000i Series		27.12	16.26	~	37.98	3.62		27.12	16.26	~	37.98	3.62
			CL-2000i/CL-1000i/CL-6000i Series		27.12	16.26	~	37.98	3.62		27.12	16.26	~	37.98	3.62
		/	/	/	/	/	~	/	/	/	/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
		/	/	/	/	/	~	/	/	/	/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
		/	/	/	/	/	~	/	/	/	/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
		/	/	/	/	/	~	/	/	/	/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
		/	/	/	/	/	~	/	/	/	/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
		/	/	/	/	/	~	/	/	/	/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
		/	/	/	/	/	~	/	/	/	/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
		/	/	/	/	/	~	/	/	/	/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/
			/		/	/	~	/	/		/	/	~	/	/

Remark:

The target value and range of
The target value and range of
The target value and range of
The target value and range of

FT4
T4
TSH
Tg

can only be used for reagent lot#
can only be used for reagent lot#
can only be used for reagent lot#
can only be used for reagent lot#

20230109
20230109
20230809
20240309

or earlier.
or earlier.
or earlier.
or earlier.