

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов измерительной системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадает во всех группах.

Português : A dados de cada grupo é a mesma.

Español : la datos de cada grupo es la misma.

Italiano : la dati di ogni gruppo è la stessa.

Türkçe : her grubun veri aynıdır.

1.BS-200E: BS-200E, BS-220E,;

9.BS-600: BS-600, BS-620;

2.BS-240: BS-230, BS-240;

10.BS-600M: BS-600M, BS-620M;

3.BS-240E: BS240E, BS240Pro;

11.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

4.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00");

12.BS-1000M: BS-1000M, BS-1100M;

5.BS-380: BS-380, BS-390;

13.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

6.BS-400: BS-400, BS-420;

14.BS-2800M:BS-2600M, BS-2800M;

7.BS-430: BS-430, BS-450,BS-460;

15.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction.

8.BS-480: BS-480, BS-490;

LOT : 046125007

: 2026-08-08

English	Abbreviated name	Calibration Value	Level	Model
Русский	сокращенное наименование	Эталонное значение	Уровень	модель
Português	Nome abreviado	Valor de calibração	Nível	Modelo
Español	nombre abreviado	Valor de calibración	Nivel	modelo
Italiano	abbreviazione	Valore di calibrazione	Livello	modelli
Türkçe	kisaltılmış ad	Kalibrasyon Değeri	Düzey	model

English	Expanded Uncertainty
Русский	Расширенная неопределенность
Português	Incerteza Expandida
Español	Incertidumbre expandida
Italiano	Incertezza estesa
Türkçe	Daha Uzun Süreli Belirsizlik

	FER
English	Ferritin
Русский	ферритин
Português	Ferritina
Español	Ferritina
Italiano	Ferritina
Türkçe	Ferritin

Abbreviated name		FER		Calibration Rule		spline			
Model		BS-200E ¹		BS-240 ²		BS-240E ³		BS-360E ⁴	
Level	Unit	Calibrator Value	Expanded Uncertainty	Calibrator Value	Expanded Uncertainty	Calibrator Value	Expanded Uncertainty	Calibrator Value	Expanded Uncertainty
0.9%NaCl		0	/	0	/	0	/	0	/
a	ng/mL	121	9	121	9	118	8	118	8
b		235	5	235	5	231	5	229	5
c		573	25	573	25	564	24	549	24

d		1066	37	1066	37	1085	38	1055	37
0.9%NaCl		0	/	0	/	0	/	0	/
a		272	19	272	19	265	19	265	19
b	pmol/L	528	12	528	12	519	11	515	11
c		1288	56	1288	56	1267	55	1234	53
d		2395	84	2395	84	2438	86	2371	83
Model		BS-380 ⁵		BS-400 ⁶		BS-430 ⁷		BS-480 ⁸	
Level	Unit	Calibrator Value	Expanded Uncertainty	Calibrator Value	Expanded Uncertainty	Calibrator Value	Expanded Uncertainty	Calibrator Value	Expanded Uncertainty
0.9%NaCl		0	/	0	/	0	/	0	/
a		114	8	114	8	119	9	114	8
b	ng/mL	223	5	223	5	235	5	223	5
c		538	23	538	23	568	25	538	23
d		1051	37	1051	37	1095	38	1051	37
0.9%NaCl		0	/	0	/	0	/	0	/
a		256	18	256	18	267	19	256	18
b	pmol/L	501	11	501	11	528	12	501	11
c		1209	52	1209	52	1276	55	1209	52
d		2362	83	2362	83	2460	86	2362	83
Model		BS-600 ⁹		BS-600M ¹⁰		BS-800 ¹¹		BS-1000M ¹²	
Level	Unit	Calibrator Value	Expanded Uncertainty	Calibrator Value	Expanded Uncertainty	Calibrator Value	Expanded Uncertainty	Calibrator Value	Expanded Uncertainty
0.9%NaCl		0	/	0	/	0	/	0	/
a		114	8	114	8	118	8	114	8
b	ng/mL	223	5	223	5	229	5	223	5
c		538	23	538	23	549	24	538	23
d		1051	37	1051	37	1055	37	1051	36
0.9%NaCl		0	/	0	/	0	/	0	/
a		256	18	256	18	265	19	256	18
b	pmol/L	501	11	501	11	515	11	501	11
c		1209	52	1209	52	1234	53	1209	52
d		2362	83	2362	83	2371	83	2362	81
Model		BS-2000 ¹³		BS-2800M ¹⁴					
Level	Unit	Calibrator Value	Expanded Uncertainty	Calibrator Value	Expanded Uncertainty				
0.9%NaCl		0	/	0	/				
a		116	8	114	8				
b	ng/mL	223	5	223	5				
c		534	23	538	23				
d		1036	36	1051	36				
0.9%NaCl		0	/	0	/				
a		261	19	256	18				
b	pmol/L	501	11	501	11				
c		1200	52	1209	52				
d		2328	80	2362	81				

mindray

FER Calibrator

For use on: BS-1000M, BS-2000,
BS-2800M, BS-600M,
BS-800

LOT 046125007

 2026-08-08

