

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-200: BS-200, BS-220, BS-330, BS-350;

11.BS-480: BS-480, BS-490;

2.BS-200E: BS-200E, BS-220E,;

12.BS-600: BS-600, BS-620;

3.BS-240: BS-230, BS-240;

13.BS-600M: BS-600M, BS-620M;

4.BS-240E: BS240E, BS240Pro;

14.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

5.BS-300: BS-300, BS-320;

15.BS-1000M: BS-1000M, BS-1100M;

6.BS-330E:BS-330E(Serial Number starts with "XQ-"), BS-350E(Serial Number starts with "XS-")

16.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

7.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with "W9-" and software version starts with "35.00");

17.BS-2800M:BS-2600M, BS-2800M.

8.BS-380: BS-380, BS-390;

18.0.9% NaCl is needed, Conc. Of 0.9% NaCl =0;

9.BS-400: BS-400, BS-420;

19.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction.

10.BS-430: BS-430, BS-450, BS-460;



: 145225005



: 2026-11-30

English	Abbreviated name	Unit	Calibration Rule	Model
Русский	сокращенное наименование	Прибор	Принцип калибровки	модель
Português	Nome abreviado	Unidade	Regra de calibração	Modelo
Español	nombre abreviado	Unidad	Regla de calibración	modelo
Italiano	abbreviazione	Unità	Regola di calibrazione	modelli
Türkçe	kısaltılmış ad	Ünite	Kalibrasyon Kuralı	model

English	Expanded Uncertainty
Русский	Расширенная неопределенность
Português	Incerteza Expandida
Español	Incertidumbre expandida
Italiano	Incertezza estesa
Türkçe	Daha Uzun Süreli Belirsizlik

CHE	
English	Cholinesterase
Русский	Холинэстераза
Português	Colinesterase
Español	Colinesterasa
Italiano	Colinesterasi
Türkçe	Kolinesteraz

Abbreviated name	Unit	Calibration Rule			Two-point linear		
		Model	Calibrator Value	Expanded Uncertainty	Model	Calibrator Value	Expanded Uncertainty
CHE	U/L	BS-200 ¹	5016	201	BS-430 ¹⁰	5016	195
		BS-200E ²	5016	195	BS-480 ¹¹	5031	201
		BS-240 ³	5016	201	BS-600 ¹²	5031	201
		BS-240E ⁴	5016	195	BS-600M ¹³	5031	195
		BS-300 ⁵	5031	201	BS-800 ¹⁴	5016	195
		BS-330E ⁶	5016	195	BS-1000M ¹⁵	5031	194
		BS-360E ⁷	5016	195	BS-2000 ¹⁶	5060	196
		BS-380 ⁸	5031	201	BS-2800M ¹⁷	5031	194
		BS-400 ⁹	5031	201			
	μkat/L	BS-200 ¹	83.8	3	BS-430 ¹⁰	83.8	3
		BS-200E ²	83.8	3	BS-480 ¹¹	84.0	3
		BS-240 ³	83.8	3	BS-600 ¹²	84.0	3
		BS-240E ⁴	83.8	3	BS-600M ¹³	84.0	3
		BS-300 ⁵	84.0	3	BS-800 ¹⁴	83.8	3
		BS-330E ⁶	83.8	3	BS-1000M ¹⁵	84.0	3
		BS-360E ⁷	83.8	3	BS-2000 ¹⁶	84.5	3
		BS-380 ⁸	84.0	3	BS-2800M ¹⁷	84.0	3
		BS-400 ⁹	84.0	3			