

## Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-240: BS-230, BS-240;

9.BS-600M: BS-600M, BS-620M;

2.BS-240E: BS240E, BS240Pro;

10.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

3.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with "W9-" and software version starts with "35.00");

11.BS-1000M: BS-1000M, BS-1100M;

12.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

4.BS-380: BS-380, BS-390;

13.BS-2800M:BS-2600M, BS-2800M.

5.BS-400: BS-400, BS-420;

14.0.9% NaCl is needed, Conc. Of 0.9% NaCl =0;

6.BS-430: BS-430, BS-450, BS-460;

15.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction.

7.BS-480: BS-480, BS-490;

8.BS-600: BS-600, BS-620;



: 146723011



: 2025-01-24

| English   | Abbreviated name                           | Unit    | Calibration Rule       | Model   |
|-----------|--------------------------------------------|---------|------------------------|---------|
| Русский   | сокращенное наименование                   | Прибор  | Принцип калибровки     | модель  |
| Português | Nome abreviado                             | Unidade | Regra de calibração    | Modelo  |
| Español   | nombre abreviado                           | Unidad  | Regla de calibración   | modelo  |
| Italiano  | abbreviazione                              | Unità   | Regola di calibrazione | modelli |
| Türkçe    | kısaltılmış ad                             | Ünite   | Kalibrasyon Kuralı     | model   |
| English   | Expanded Uncertainty                       |         |                        |         |
| Русский   | Расширенная неопределенность               |         |                        |         |
| Português | Incerteza Expandida                        |         |                        |         |
| Español   | Incertidumbre expandida                    |         |                        |         |
| Italiano  | Incertezza estesa                          |         |                        |         |
| Türkçe    | Daha Uzun Süreli Belirsizlik               |         |                        |         |
| UIBC      |                                            |         |                        |         |
| English   | Unsaturated Iron Binding Capacity          |         |                        |         |
| Русский   | ненасыщенная железосвязывающая способность |         |                        |         |
| Português | Capacidade de ligação de ferro insaturado  |         |                        |         |
| Español   | Capacidad de unión de hierro no saturado   |         |                        |         |
| Italiano  | Capacità di legame del ferro insaturo      |         |                        |         |
| Türkçe    | Unsature Demir Bağlama Kapasitesi          |         |                        |         |

| Abbreviated name | Unit   | Calibration Rule     |                  |                      | Two-point linear       |                      |
|------------------|--------|----------------------|------------------|----------------------|------------------------|----------------------|
|                  |        | Model                | Calibrator Value | Expanded Uncertainty | Model                  | Expanded Uncertainty |
| UIBC             | μmol/L | BS-240 <sup>1</sup>  | /                | /                    | BS-600 <sup>8</sup>    | /                    |
|                  |        | BS-240E <sup>2</sup> | /                | /                    | BS-600M <sup>9</sup>   | /                    |
|                  |        | BS-360E <sup>3</sup> | /                | /                    | BS-800 <sup>10</sup>   | /                    |
|                  |        | BS-380 <sup>4</sup>  | /                | /                    | BS-1000M <sup>11</sup> | 114                  |
|                  |        | BS-400 <sup>5</sup>  | /                | /                    | BS-2000 <sup>12</sup>  | /                    |
|                  |        | BS-430 <sup>6</sup>  | /                | /                    | BS-2800M <sup>13</sup> | /                    |
|                  |        | BS-480 <sup>7</sup>  | /                | /                    |                        |                      |
|                  | μg/dL  | BS-240 <sup>1</sup>  | /                | /                    | BS-600 <sup>8</sup>    | /                    |
|                  |        | BS-240E <sup>2</sup> | /                | /                    | BS-600M <sup>9</sup>   | /                    |
|                  |        | BS-360E <sup>3</sup> | /                | /                    | BS-800 <sup>10</sup>   | /                    |
|                  |        | BS-380 <sup>4</sup>  | /                | /                    | BS-1000M <sup>11</sup> | 637                  |
|                  |        | BS-400 <sup>5</sup>  | /                | /                    | BS-2000 <sup>12</sup>  | /                    |
|                  |        | BS-430 <sup>6</sup>  | /                | /                    | BS-2800M <sup>13</sup> | /                    |
|                  |        | BS-480 <sup>7</sup>  | /                | /                    |                        |                      |

# **mindray**

## **UIBC Calibrator**

For use on: BS-1000M, BS-2000,  
BS-2800M, BS-600M,  
BS-800

**LOT** 146723011

 **2025-01-24**

