

Value sheet of Mindray BS Measurement System

Русский : Таблица результатов для системы BS компании Mindray

Português : Planilha de valores do Sistema de Medição BS da Mindray

Español : Hoja de valores del sistema de medición Mindray BS

Italiano : Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe : Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский : Данные совпадают во всех группах.

Português : A dados de cada grupo é a mesma.

Español : la datos de cada grupo es la misma.

Italiano : la dati di ogni gruppo è la stessa.

Türkçe : her grubun veri aynıdır.

1.BS-120: BS-120, BS-130, BS-180, BS-190;

11.BS-480: BS-480, BS-490;

2.BS-200: BS-200, BS-220, BS-330, BS-350;

12.BS-600: BS-600, BS-620;

3.BS-200E: BS-200E, BS-220E,;

13.BS-600M: BS-600M, BS-620M;

4.BS-240: BS-230, BS-240;

14.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

5.BS-240E: BS240E, BS240Pro;

6.BS-330E:BS-330E(Serial Number starts with "XQ-"),
BS-350E(Serial Number starts with "XS-")

15.BS-1000M: BS-1000M, BS-1100M;

7.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with "W9-" and software version starts with "35.00");

16.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

8.BS-380: BS-380, BS-390;

18.0.9% NaCl is needed, Conc. Of 0.9% NaCl =0;

9.BS-400: BS-400, BS-420;

19.For applicable models of the analyte, please refer to the reagent

10.BS-430: BS-430, BS-450, BS-460;

parameter sheet and instruction.



LOT : 147925005



: 2026-07-09

English	Abbreviated name	Unit	Calibration Rule	Model
Русский	сокращенное наименование	Прибор	Принцип калибровки	модель
Português	Nome abreviado	Unidade	Regra de calibração	Modelo
Español	nombre abreviado	Unidad	Regla de calibración	modelo
Italiano	abbreviazione	Unità	Regola di calibrazione	modelli
Türkçe	kısaltılmış ad	Ünite	Kalibrasyon Kuralı	model

English	Expanded Uncertainty
Русский	Расширенная неопределенность
Português	Incerteza Expandida
Español	Incertidumbre expandida
Italiano	Incertezza estesa
Türkçe	Daha Uzun Süreli Belirsizlik

TPUC	
English	Total Protein in Urine/CSF (TPUC)
Русский	Общий белок в моче/СМЖ (TPUC)
Português	Proteína total em urina e LCR (PTUL)
Español	proteína total en orina/LCR (TPUC)
Italiano	Proteine Totali Urinarie/CSF (TPUC)
Türkçe	BOS/İdrar Protein (TPUC)

Abbreviated name	Unit	Calibration Rule				Two-point linear	
		Model	Calibrator Value	Expanded Uncertainty	Model	Calibrator Value	Expanded Uncertainty
TPUC	mg/L	BS-120 ¹	1260	17	BS-430 ¹⁰	1312	18
		BS-200 ²	1260	17	BS-480 ¹¹	1312	18
		BS-200E ³	1312	18	BS-600 ¹²	1322	18
		BS-240 ⁴	1260	17	BS-600M ¹³	1322	18
		BS-240E ⁵	1312	18	BS-800 ¹⁴	1312	18
		BS-330E ⁶	1312	18	BS-1000M ¹⁵	1322	18
		BS-360E ⁷	1312	18	BS-2000 ¹⁶	1336	19
		BS-380 ⁸	1322	18	BS-2800M ¹⁷	1322	18
		BS-400 ⁹	1322	18			

mindray

TPUC Calibrator

For use on: BS-1000M, BS-2000,
BS-2800M, BS-600M,
BS-800

LOT 147925005

 2026-07-09

