

Value sheet of Mindray BS Measurement System

Русский : Таблица результатов для системы BS компании Mindray  
Português : Planilha de valores do Sistema de Medição BS da Mindray  
Español : Hoja de valores del sistema de medición Mindray BS  
Italiano : Scheda dei valori del sistema di misurazione BS di Mindray  
Türkçe : Mindray BS Ölçüm Sistemi'nin değer sayfası



|                                                                                                                                                                                                                                      |  |  |                                                                                                       |  |  |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|--|-------------------------------------------------------------------------------------------------------|--|--|
| The data of each group is same.                                                                                                                                                                                                      |  |  | Русский: Данные совпадают во всех группах.                                                            |  |  |
| Português : A dados de cada grupo é a mesma.                                                                                                                                                                                         |  |  | Español : la datos de cada grupo es la misma.                                                         |  |  |
| Italiano : la dati di ogni gruppo è la stessa.                                                                                                                                                                                       |  |  | Türkçe : her grubun veri aynıdır.                                                                     |  |  |
| 1.BS-120: BS-120, BS-130, BS-180, BS-190;                                                                                                                                                                                            |  |  | 11.BS-430: BS-430, BS-450,BS-460;                                                                     |  |  |
| 2.BS-200: BS-200, BS-220, BS-330, BS-350;                                                                                                                                                                                            |  |  | 12.BS-480: BS-480, BS-490;                                                                            |  |  |
| 3.BS-200E: BS-200E, BS-220E,;                                                                                                                                                                                                        |  |  | 13.BS-600: BS-600, BS-620;                                                                            |  |  |
| 4.BS-240: BS-230, BS-240;                                                                                                                                                                                                            |  |  | 14.BS-600M: BS-600M, BS-620M;                                                                         |  |  |
| 5.BS-240E: BS240E, BS240Pro;                                                                                                                                                                                                         |  |  | 15.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;                                    |  |  |
| 6.BS-300: BS-300, BS-320;                                                                                                                                                                                                            |  |  | 16.BS-1000M: BS-1000M, BS-1100M;                                                                      |  |  |
| 7.BS-330E:BS-330E(Serial Number starts with "XQ-"), BS-350E(Serial Number starts with "XS-")                                                                                                                                         |  |  | 17.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;                                                     |  |  |
| 8.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00"); |  |  | 18.BS-2800M:BS-2600M, BS-2800M;                                                                       |  |  |
| 9.BS-380: BS-380, BS-390;                                                                                                                                                                                                            |  |  | 19.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction. |  |  |
| 10.BS-400: BS-400, BS-420;                                                                                                                                                                                                           |  |  |                                                                                                       |  |  |

LOT : 044824005      : 2025-07-31

| English   | Abbreviated name         | Model   | Unit    | Assay Value        | Range (Assay Value±3SD)                    |
|-----------|--------------------------|---------|---------|--------------------|--------------------------------------------|
| Русский   | сокращенное наименование | модель  | Прибор  | Результат анализа  | Диапазон(результат анализа ± 3CO)          |
| Português | Nome abreviado           | Modelo  | Unidade | Valores da análise | Faixa(Valores da análise ±3SD)             |
| Español   | nombre abreviado         | modelo  | Unidad  | Valor de ensayo    | Rango(Valor de ensayo ±3SD)                |
| Italiano  | abbreviazione            | modelli | Unità   | Valori di dosaggio | Intervallo (valore di concentrazione±3 DS) |
| Türkçe    | kısaltılmış ad           | model   | Ünite   | Tayin Değeri       | Aralık (Tayin Değeri±3SD)                  |
| LIP       |                          |         |         |                    |                                            |
| English   |                          |         | Lipase  |                    |                                            |
| Русский   |                          |         | Липаза  |                    |                                            |
| Português |                          |         | Lipase  |                    |                                            |
| Español   |                          |         | Lipasa  |                    |                                            |
| Italiano  |                          |         | Lipasi  |                    |                                            |
| Türkçe    |                          |         | Lipaz   |                    |                                            |

| Abbreviated name     | Unit  | Model                | Assay Value         | Range(Assay Value±3SD) |       |       | 1 SD                  | Model                  | Assay Value          | Range(Assay Value±3SD) |       |       | 1 SD  |
|----------------------|-------|----------------------|---------------------|------------------------|-------|-------|-----------------------|------------------------|----------------------|------------------------|-------|-------|-------|
| LIP                  | U/L   | BS-120 <sup>1</sup>  | 47.2                | 37.6                   | —     | 56.8  | 3.2                   | BS-400 <sup>10</sup>   | 47.0                 | 37.7                   | —     | 56.3  | 3.1   |
|                      |       | BS-200 <sup>2</sup>  | 47.2                | 37.6                   | —     | 56.8  | 3.2                   | BS-430 <sup>11</sup>   | 45.0                 | 36.0                   | —     | 54.0  | 3.0   |
|                      |       | BS-200E <sup>3</sup> | 43.6                | 34.9                   | —     | 52.3  | 2.9                   | BS-480 <sup>12</sup>   | 47.0                 | 37.7                   | —     | 56.3  | 3.1   |
|                      |       | BS-240 <sup>4</sup>  | 47.2                | 37.6                   | —     | 56.8  | 3.2                   | BS-600 <sup>13</sup>   | 47.0                 | 37.7                   | —     | 56.3  | 3.1   |
|                      |       | BS-240E <sup>5</sup> | 46.6                | 37.3                   | —     | 55.9  | 3.1                   | BS-600M <sup>14</sup>  | 46.0                 | 36.7                   | —     | 55.3  | 3.1   |
|                      |       | BS-300 <sup>6</sup>  | 47.0                | 37.7                   | —     | 56.3  | 3.1                   | BS-800 <sup>15</sup>   | 46.5                 | 37.2                   | —     | 55.8  | 3.1   |
|                      |       | BS-330E <sup>7</sup> | 43.6                | 34.9                   | —     | 52.3  | 2.9                   | BS-1000M <sup>16</sup> | 47.0                 | 37.7                   | —     | 56.3  | 3.1   |
|                      |       | BS-360E <sup>8</sup> | 43.6                | 34.9                   | —     | 52.3  | 2.9                   | BS-2000 <sup>17</sup>  | 47.1                 | 37.5                   | —     | 56.7  | 3.2   |
|                      |       | BS-380 <sup>9</sup>  | 47.0                | 37.7                   | —     | 56.3  | 3.1                   | BS-2800M <sup>18</sup> | 47.0                 | 37.7                   | —     | 56.3  | 3.1   |
|                      |       | µkat/L               | BS-120 <sup>1</sup> | 0.788                  | 0.628 | —     | 0.949                 | 0.053                  | BS-400 <sup>10</sup> | 0.785                  | 0.630 | —     | 0.940 |
| BS-200 <sup>2</sup>  | 0.788 |                      | 0.628               | —                      | 0.949 | 0.053 | BS-430 <sup>11</sup>  | 0.752                  | 0.601                | —                      | 0.902 | 0.050 |       |
| BS-200E <sup>3</sup> | 0.728 |                      | 0.583               | —                      | 0.873 | 0.048 | BS-480 <sup>12</sup>  | 0.785                  | 0.630                | —                      | 0.940 | 0.052 |       |
| BS-240 <sup>4</sup>  | 0.788 |                      | 0.628               | —                      | 0.949 | 0.053 | BS-600 <sup>13</sup>  | 0.785                  | 0.630                | —                      | 0.940 | 0.052 |       |
| BS-240E <sup>5</sup> | 0.778 |                      | 0.623               | —                      | 0.934 | 0.052 | BS-600M <sup>14</sup> | 0.768                  | 0.613                | —                      | 0.924 | 0.052 |       |

| Abbreviated name | Unit | Model                | Assay Value | Range(Assay Value±3SD) |   |       | 1 SD  | Model                  | Assay Value | Range(Assay Value±3SD) |   |       | 1 SD  |
|------------------|------|----------------------|-------------|------------------------|---|-------|-------|------------------------|-------------|------------------------|---|-------|-------|
|                  |      | BS-300 <sup>6</sup>  | 0.785       | 0.630                  | — | 0.940 | 0.052 | BS-800 <sup>15</sup>   | 0.777       | 0.621                  | — | 0.932 | 0.052 |
|                  |      | BS-330E <sup>7</sup> | 0.728       | 0.583                  | — | 0.873 | 0.048 | BS-1000M <sup>16</sup> | 0.785       | 0.630                  | — | 0.940 | 0.052 |
|                  |      | BS-360E <sup>8</sup> | 0.728       | 0.583                  | — | 0.873 | 0.048 | BS-2000 <sup>17</sup>  | 0.787       | 0.626                  | — | 0.947 | 0.053 |
|                  |      | BS-380 <sup>9</sup>  | 0.785       | 0.630                  | — | 0.940 | 0.052 | BS-2800M <sup>18</sup> | 0.785       | 0.630                  | — | 0.940 | 0.052 |