

Value sheet of Mindray BS Measurement System

Русский : Таблица результатов для системы BS компании Mindray  
Português : Planilha de valores do Sistema de Medição BS da Mindray  
Español : Hoja de valores del sistema de medición Mindray BS  
Italiano : Scheda dei valori del sistema di misurazione BS di Mindray  
Türkçe : Mindray BS Ölçüm Sistemi'nin değer sayfası



|                                                                                                                                                                                                                                      |  |  |                                                                                                       |  |  |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|--|-------------------------------------------------------------------------------------------------------|--|--|
| The data of each group is same.                                                                                                                                                                                                      |  |  | Русский: Данные совпадают во всех группах.                                                            |  |  |
| Português : A dados de cada grupo é a mesma.                                                                                                                                                                                         |  |  | Español : la datos de cada grupo es la misma.                                                         |  |  |
| Italiano : la dati di ogni gruppo è la stessa.                                                                                                                                                                                       |  |  | Türkçe : her grubun veri aynıdır.                                                                     |  |  |
| 1.BS-120: BS-120, BS-130, BS-180, BS-190;                                                                                                                                                                                            |  |  | 11.BS-430: BS-430, BS-450,BS-460;                                                                     |  |  |
| 2.BS-200: BS-200, BS-220, BS-330, BS-350;                                                                                                                                                                                            |  |  | 12.BS-480: BS-480, BS-490;                                                                            |  |  |
| 3.BS-200E: BS-200E, BS-220E,;                                                                                                                                                                                                        |  |  | 13.BS-600: BS-600, BS-620;                                                                            |  |  |
| 4.BS-240: BS-230, BS-240;                                                                                                                                                                                                            |  |  | 14.BS-600M: BS-600M, BS-620M;                                                                         |  |  |
| 5.BS-240E: BS240E, BS240Pro;                                                                                                                                                                                                         |  |  | 15.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;                                    |  |  |
| 6.BS-300: BS-300, BS-320;                                                                                                                                                                                                            |  |  | 16.BS-1000M: BS-1000M, BS-1100M;                                                                      |  |  |
| 7.BS-330E:BS-330E(Serial Number starts with "XQ-"), BS-350E(Serial Number starts with "XS-")                                                                                                                                         |  |  | 17.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;                                                     |  |  |
| 8.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00"); |  |  | 18.BS-2800M:BS-2600M, BS-2800M;                                                                       |  |  |
| 9.BS-380: BS-380, BS-390;                                                                                                                                                                                                            |  |  | 19.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction. |  |  |
| 10.BS-400: BS-400, BS-420;                                                                                                                                                                                                           |  |  |                                                                                                       |  |  |

LOT: 044824013      : 2025-11-30

| English   | Abbreviated name         | Model   | Unit    | Assay Value        | Range (Assay Value±3SD)                    |
|-----------|--------------------------|---------|---------|--------------------|--------------------------------------------|
| Русский   | сокращенное наименование | модель  | Прибор  | Результат анализа  | Диапазон(результат анализа ± 3CO)          |
| Português | Nome abreviado           | Modelo  | Unidade | Valores da análise | Faixa(Valores da análise ±3SD)             |
| Español   | nombre abreviado         | modelo  | Unidad  | Valor de ensayo    | Rango(Valor de ensayo ±3SD)                |
| Italiano  | abbreviazione            | modelli | Unità   | Valori di dosaggio | Intervallo (valore di concentrazione±3 DS) |
| Türkçe    | kısaltılmış ad           | model   | Ünite   | Tayin Değeri       | Aralık (Tayin Değeri±3SD)                  |
| LIP       |                          |         |         |                    |                                            |
| English   |                          |         | Lipase  |                    |                                            |
| Русский   |                          |         | Липаза  |                    |                                            |
| Português |                          |         | Lipase  |                    |                                            |
| Español   |                          |         | Lipasa  |                    |                                            |
| Italiano  |                          |         | Lipasi  |                    |                                            |
| Türkçe    |                          |         | Lipaz   |                    |                                            |

| Abbreviated name | Unit | Model                | Assay Value | Range(Assay Value±3SD) | 1 SD  | Model                  | Assay Value | Range(Assay Value±3SD) | 1 SD  |
|------------------|------|----------------------|-------------|------------------------|-------|------------------------|-------------|------------------------|-------|
| LIP              | U/L  | BS-120 <sup>1</sup>  | 46.5        | 37.2 — 55.8            | 3.1   | BS-400 <sup>10</sup>   | 46.7        | 37.4 — 56.0            | 3.1   |
|                  |      | BS-200 <sup>2</sup>  | 46.5        | 37.2 — 55.8            | 3.1   | BS-430 <sup>11</sup>   | 45.1        | 36.1 — 54.1            | 3.0   |
|                  |      | BS-200E <sup>3</sup> | 44.1        | 35.1 — 53.1            | 3.0   | BS-480 <sup>12</sup>   | 46.7        | 37.4 — 56.0            | 3.1   |
|                  |      | BS-240 <sup>4</sup>  | 46.5        | 37.2 — 55.8            | 3.1   | BS-600 <sup>13</sup>   | 46.7        | 37.4 — 56.0            | 3.1   |
|                  |      | BS-240E <sup>5</sup> | 46.2        | 36.9 — 55.5            | 3.1   | BS-600M <sup>14</sup>  | 45.7        | 36.4 — 55.0            | 3.1   |
|                  |      | BS-300 <sup>6</sup>  | 46.7        | 37.4 — 56.0            | 3.1   | BS-800 <sup>15</sup>   | 46.3        | 37.0 — 55.6            | 3.1   |
|                  |      | BS-330E <sup>7</sup> | 44.1        | 35.1 — 53.1            | 3.0   | BS-1000M <sup>16</sup> | 46.7        | 37.4 — 56.0            | 3.1   |
|                  |      | BS-360E <sup>8</sup> | 44.1        | 35.1 — 53.1            | 3.0   | BS-2000 <sup>17</sup>  | 47.4        | 37.8 — 57.0            | 3.2   |
|                  |      | BS-380 <sup>9</sup>  | 46.7        | 37.4 — 56.0            | 3.1   | BS-2800M <sup>18</sup> | 46.7        | 37.4 — 56.0            | 3.1   |
|                  |      | BS-120 <sup>1</sup>  | 0.777       | 0.621 — 0.932          | 0.052 | BS-400 <sup>10</sup>   | 0.780       | 0.625 — 0.935          | 0.052 |
|                  |      | BS-200 <sup>2</sup>  | 0.777       | 0.621 — 0.932          | 0.052 | BS-430 <sup>11</sup>   | 0.753       | 0.603 — 0.903          | 0.050 |
|                  |      | BS-200E <sup>3</sup> | 0.736       | 0.586 — 0.887          | 0.050 | BS-480 <sup>12</sup>   | 0.780       | 0.625 — 0.935          | 0.052 |
|                  |      | BS-240 <sup>4</sup>  | 0.777       | 0.621 — 0.932          | 0.052 | BS-600 <sup>13</sup>   | 0.780       | 0.625 — 0.935          | 0.052 |

| Abbreviated name | Unit   | Model                | Assay Value | Range(Assay Value±3SD) |   |       | 1 SD  | Model                  | Assay Value | Range(Assay Value±3SD) |   |       | 1 SD  |
|------------------|--------|----------------------|-------------|------------------------|---|-------|-------|------------------------|-------------|------------------------|---|-------|-------|
|                  | μkat/L | BS-240E <sup>5</sup> | 0.772       | 0.616                  | — | 0.927 | 0.052 | BS-600M <sup>14</sup>  | 0.763       | 0.608                  | — | 0.919 | 0.052 |
|                  |        | BS-300 <sup>6</sup>  | 0.780       | 0.625                  | — | 0.935 | 0.052 | BS-800 <sup>15</sup>   | 0.773       | 0.618                  | — | 0.929 | 0.052 |
|                  |        | BS-330E <sup>7</sup> | 0.736       | 0.586                  | — | 0.887 | 0.050 | BS-1000M <sup>16</sup> | 0.780       | 0.625                  | — | 0.935 | 0.052 |
|                  |        | BS-360E <sup>8</sup> | 0.736       | 0.586                  | — | 0.887 | 0.050 | BS-2000 <sup>17</sup>  | 0.792       | 0.631                  | — | 0.952 | 0.053 |
|                  |        | BS-380 <sup>9</sup>  | 0.780       | 0.625                  | — | 0.935 | 0.052 | BS-2800M <sup>18</sup> | 0.780       | 0.625                  | — | 0.935 | 0.052 |