

Value sheet of Mindray BS Measurement System

Русский : Таблица результатов для системы BS компании Mindray
Português : Planilha de valores do Sistema de Medição BS da Mindray
Español : Hoja de valores del sistema de medición Mindray BS
Italiano : Scheda dei valori del sistema di misurazione BS di Mindray
Türkçe : Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Português : A dados de cada grupo é a mesma.

Italiano : la dati di ogni gruppo è la stessa.

1.BS-120: BS-120, BS-130, BS-180, BS-190;
2.BS-200: BS-200, BS-220, BS-330, BS-350;
3.BS-200E: BS-200E, BS-220E,;
4.BS-240: BS-230, BS-240;
5.BS-240E: BS240E, BS240Pro;
6.BS-330E:BS-330E(Serial Number starts with "XQ-"), BS-350E(Serial Number starts with "XS-")
7.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00");
8.BS-380: BS-380, BS-390;
9.BS-400: BS-400, BS-420;

Русский: Данные совпадают во всех группах.

Español : la datos de cada grupo es la misma.

Türkçe : her grubun veri aynıdır.

10.BS-430: BS-430, BS-450,BS-460;
11.BS-480: BS-480, BS-490;
12.BS-600: BS-600, BS-620;
13.BS-600M: BS-600M, BS-620M;
14.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;
15.BS-1000M: BS-1000M, BS-1100M;
16.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;
17.BS-2800M:BS-2600M, BS-2800M;
18.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction.

LOT : 053125010

: 2027-03-22

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

MALB					
English	Microalbumin				
Русский	микроальбумина				
Português	Microalbuminúria				
Español	Microalbúmina				
Italiano	Microalbumina				
Türkçe	Mikroalbumin				

Abbreviated name	Unit	Model	Assay Value	Range(Assay Value±3SD)			1 SD	Model	Assay Value	Range(Assay Value±3SD)			1 SD
MALB	mg/L	BS-120 ¹	196	166	—	226	10	BS-430 ¹⁰	193	163	—	223	10
		BS-200 ²	196	166	—	226	10	BS-480 ¹¹	192	162	—	222	10
		BS-200E ³	191	161	—	221	10	BS-600 ¹²	192	162	—	222	10
		BS-240 ⁴	196	166	—	226	10	BS-600M ¹³	192	162	—	222	10
		BS-240E ⁵	191	161	—	221	10	BS-800 ¹⁴	191	161	—	221	10
		BS-360E ⁷	191	161	—	221	10	BS-1000M ¹⁵	192	162	—	222	10
		BS-380 ⁸	192	162	—	222	10	BS-2000 ¹⁶	192	162	—	222	10
		BS-400 ⁹	192	162	—	222	10	BS-2800M ¹⁷	192	162	—	222	10
	μmol/L	BS-120 ¹	2.98	2.52	—	3.44	0.15	BS-430 ¹⁰	2.93	2.48	—	3.39	0.15
		BS-200 ²	2.98	2.52	—	3.44	0.15	BS-480 ¹¹	2.92	2.46	—	3.37	0.15
		BS-200E ³	2.90	2.45	—	3.36	0.15	BS-600 ¹²	2.92	2.46	—	3.37	0.15
		BS-240 ⁴	2.98	2.52	—	3.44	0.15	BS-600M ¹³	2.92	2.46	—	3.37	0.15
		BS-240E ⁵	2.90	2.45	—	3.36	0.15	BS-800 ¹⁴	2.90	2.45	—	3.36	0.15
		BS-360E ⁷	2.90	2.45	—	3.36	0.15	BS-1000M ¹⁵	2.92	2.46	—	3.37	0.15

Abbreviated name	Unit	Model	Assay Value	Range(Assay Value±3SD)			1 SD	Model	Assay Value	Range(Assay Value±3SD)			1 SD
		BS-380 ⁸	2.92	2.46	—	3.37	0.15	BS-2000 ¹⁶	2.92	2.46	—	3.37	0.15
		BS-400 ⁹	2.92	2.46	—	3.37	0.15	BS-2800M ¹⁷	2.92	2.46	—	3.37	0.15

mindray

MALB Control

For use on: BS-1000M, BS-2000,
BS-2800M, BS-600M,
BS-800

LOT 053125010

 **2027-03-22**

