

Value sheet of Mindray BS Measurement System

Русский : Таблица результатов для системы BS компании Mindray
Português : Planilha de valores do Sistema de Medição BS da Mindray
Español : Hoja de valores del sistema de medición Mindray BS
Italiano : Scheda dei valori del sistema di misurazione BS di Mindray
Türkçe : Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Português : A dados de cada grupo é a mesma.

Italiano : la dati di ogni gruppo è la stessa.

1.BS-240: BS-230, BS-240;

2.BS-240E: BS240E, BS240Pro;

3.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00");

4.BS-380: BS-380, BS-390;

5.BS-400: BS-400, BS-420;

6.BS-430: BS-430, BS-450,BS-460;

7.BS-480: BS-480, BS-490;

Русский: Данные совпадают во всех группах.

Español : la datos de cada grupo es la misma.

Türkçe : her grubun veri aynıdır.

8.BS-600: BS-600, BS-620;

9.BS-600M: BS-600M, BS-620M;

10.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

11.BS-1000M: BS-1000M, BS-1100M;

12.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

13.BS-2800M:BS-2600M, BS-2800M;

14.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction.

LOT : 053524001

🕒 : 2025-04-16

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)
UIBC					
English	Unsaturated Iron Binding Capacity				
Русский	ненасыщенная железосвязывающая способность				
Português	Capacidade de ligação de ferro insaturado				
Español	Capacidad de unión de hierro no saturado				
Italiano	Capacità di legame del ferro insaturo				
Türkçe	Unsature Demir Bağlama Kapasitesi				

Abbreviated name	Unit	Model	Assay Value	Range(Assay Value±3SD)			1 SD	Model	Assay Value	Range(Assay Value±3SD)			1 SD
UIBC	μmol/L	BS-240 ¹	/	/	—	/	/	BS-600 ⁸	/	/	—	/	/
		BS-240E ²	/	/	—	/	/	BS-600M ⁹	/	/	—	/	/
		BS-360E ³	/	/	—	/	/	BS-800 ¹⁰	/	/	—	/	/
		BS-380 ⁴	/	/	—	/	/	BS-1000M ¹¹	32.7	26.1	—	39.3	2.2
		BS-400 ⁵	/	/	—	/	/	BS-2000 ¹²	/	/	—	/	/
		BS-430 ⁶	/	/	—	/	/	BS-2800M ¹³	/	/	—	/	/
		BS-480 ⁷	/	/	—	/	/						
	μg/dL	BS-240 ¹	/	/	—	/	/	BS-600 ⁸	/	/	—	/	/
		BS-240E ²	/	/	—	/	/	BS-600M ⁹	/	/	—	/	/
		BS-360E ³	/	/	—	/	/	BS-800 ¹⁰	/	/	—	/	/
		BS-380 ⁴	/	/	—	/	/	BS-1000M ¹¹	183	146	—	220	12
		BS-400 ⁵	/	/	—	/	/	BS-2000 ¹²	/	/	—	/	/
		BS-430 ⁶	/	/	—	/	/	BS-2800M ¹³	/	/	—	/	/
		BS-480 ⁷	/	/	—	/	/						

mindray
UIBC Control

For use on: BS-1000M, BS-2000,
BS-2800M, BS-600M,
BS-800

LOT 053524001

2025-04-16

