

Value sheet of Mindray BS Measurement System

Русский : Таблица результатов для системы BS компании Mindray
Português : Planilha de valores do Sistema de Medição BS da Mindray
Español : Hoja de valores del sistema de medición Mindray BS
Italiano : Scheda dei valori del sistema di misurazione BS di Mindray
Türkçe : Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.			Русский: Данные совпадают во всех группах.		
Português : A dados de cada grupo é a mesma.			Español : la datos de cada grupo es la misma.		
Italiano : la dati di ogni gruppo è la stessa.			Türkçe : her grubun veri aynıdır.		
1.BS-200E: BS-200E, BS-220E,;			9.BS-600: BS-600, BS-620;		
2.BS-240: BS-230, BS-240;			10.BS-600M: BS-600M, BS-620M;		
3.BS-240E: BS240E, BS240Pro;			11.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;		
4.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00");			12.BS-1000M: BS-1000M, BS-1100M;		
5.BS-380: BS-380, BS-390;			13.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;		
6.BS-400: BS-400, BS-420;			14.BS-2800M:BS-2600M, BS-2800M;		
7.BS-430: BS-430, BS-450,BS-460;			15.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction.		
8.BS-480: BS-480, BS-490;					

LOT : 054625006



: 2026-06-19

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)
RBP					
English	Retinol Binding Protein				
Русский	Ретинол-связывающий белок				
Português	Proteína transportadora do retinol				
Español	Proteína fijadora de retinol				
Italiano	Proteina legante il retinolo				
Türkçe	Retinol Bağlayıcı Protein				

Abbreviated name	Unit	Model	Assay Value	Range(Assay Value±3SD)			1 SD	Model	Assay Value	Range(Assay Value±3SD)			1 SD
RBP	mg/dL	BS-200E ¹	4.71	3.75	—	5.67	0.32	BS-480 ⁸	4.77	3.81	—	5.73	0.32
		BS-240 ²	4.79	3.83	—	5.75	0.32	BS-600 ⁹	4.77	3.81	—	5.73	0.32
		BS-240E ³	4.71	3.75	—	5.67	0.32	BS-600M ¹⁰	4.77	3.81	—	5.73	0.32
		BS-360E ⁴	4.71	3.75	—	5.67	0.32	BS-800 ¹¹	4.71	3.75	—	5.67	0.32
		BS-380 ⁵	4.77	3.81	—	5.73	0.32	BS-1000M ¹²	4.77	3.81	—	5.73	0.32
		BS-400 ⁶	4.77	3.81	—	5.73	0.32	BS-2000 ¹³	4.79	3.83	—	5.75	0.32
		BS-430 ⁷	4.71	3.75	—	5.67	0.32	BS-2800M ¹⁴	4.77	3.81	—	5.73	0.32

mindray
RBP Control

For use on: BS-1000M, BS-2000,
BS-2800M, BS-600M,
BS-800

LOT 054625006

2026-06-19

