

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası

**The data of each group is same.****Português: A dados de cada grupo é a mesma.****Italiano: la dati di ogni gruppo è la stessa.**

1.BS-200E: BS-200E, BS-220E,;

2.BS-240: BS-230, BS-240;

3.BS-240E: BS240E, BS240Pro;

4.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-

330E(V35.00)(Serial Number starts with "W8-" and software version

starts with "35.00"), BS-350E(V35.00) (Serial Number starts with "W9-"

and software version starts with "35.00");

5.BS-380: BS-380, BS-390;

6.BS-400: BS-400, BS-420;

7.BS-430: BS-430, BS-450,BS-460;

8.BS-480: BS-480, BS-490;

Русский: Данные совпадают во всех группах.**Español: la datos de cada grupo es la misma.****Türkçe: her grubun veri aynıdır.**

9.BS-600: BS-600, BS-620;

10.BS-600M: BS-600M, BS-620M;

11.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

12.BS-1000M: BS-1000M, BS-1100M;

13.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

14.BS-2800M:BS-2600M, BS-2800M;

15.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction.

LOT: 054625007**: 2026-08-25**

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value $\pm$ 3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа $\pm$ 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise $\pm$ 3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo $\pm$ 3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione $\pm$ 3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri $\pm$ 3SD)

RBP

English	Retinol Binding Protein
Русский	Ретинол-связывающий белок
Português	Proteína transportadora do retinol
Español	Proteína fijadora de retinol
Italiano	Proteina legante il retinolo
Türkçe	Retinol Bağlayıcı Protein

Abbreviated name	Unit	Model	Assay Value	Range(Assay Value $\pm$ 3SD)	1 SD	Model	Assay Value	Range(Assay Value $\pm$ 3SD)	1 SD
RBP	mg/dL	BS-200E¹	4.70	3.77 — 5.63	0.31	BS-480⁸	4.69	3.76 — 5.62	0.31
		BS-240²	4.70	3.77 — 5.63	0.31	BS-600⁹	4.69	3.76 — 5.62	0.31
		BS-240E³	4.70	3.77 — 5.63	0.31	BS-600M¹⁰	4.69	3.76 — 5.62	0.31
		BS-360E⁴	4.70	3.77 — 5.63	0.31	BS-800¹¹	4.70	3.77 — 5.63	0.31
		BS-380⁵	4.69	3.76 — 5.62	0.31	BS-1000M¹²	4.69	3.76 — 5.62	0.31
		BS-400⁶	4.69	3.76 — 5.62	0.31	BS-2000¹³	4.69	3.76 — 5.62	0.31
		BS-430⁷	4.70	3.77 — 5.63	0.31	BS-2800M¹⁴	4.69	3.76 — 5.62	0.31

mindray
RBP Control

For use on: BS-1000M, BS-2000,
BS-2800M, BS-600M,
BS-800

LOT 054625007

2026-08-25

