

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-120: BS-120, BS-130, BS-180, BS-190;

9.BS-480: BS-480, BS-490;

2.BS-200: BS-200, BS-220, BS-330, BS-350;

10.BS-600: BS-600, BS-620;

3.BS-200E: BS-200E, BS-220E,;

11.BS-600M: BS-600M, BS-620M;

4.BS-240: BS-230, BS-240;

12.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

5.BS-240E: BS240E, BS240Pro;

13.BS-1000M: BS-1000M, BS-1100M;

6.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-

14.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with "W9-" and software version starts with "35.00");

15.BS-2800M: BS-2600M, BS-2800M;

7.BS-380: BS-380, BS-390;

16.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction.

8.BS-400: BS-400, BS-420;

LOT: 057523006

EXP: 2025-02-10

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value $\pm$ 3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа $\pm$ 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise $\pm$ 3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo $\pm$ 3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione $\pm$ 3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri $\pm$ 3SD)

G6PD

English	Glucose-6-Phosphate Dehydrogenase
Русский	глюкозо-6-фосфатдегидрогеназа
Português	Glicose-6-fosfato desidrogenase
Español	Glucosa-6-fosfato deshidrogenasa
Italiano	Glucosio-6-fosfato deidrogenasi
Türkçe	Glikoz 6 Fosfat Dehidrogenaz

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
G6PD	U/L	BS-120 ¹	/	/	—	/	/	BS-120 ¹	/	/	—	/	/
		BS-200 ²	/	/	—	/	/	BS-200 ²	/	/	—	/	/
		BS-200E ³	/	/	—	/	/	BS-200E ³	/	/	—	/	/
		BS-240 ⁴	/	/	—	/	/	BS-240 ⁴	/	/	—	/	/
		BS-240E ⁵	/	/	—	/	/	BS-240E ⁵	/	/	—	/	/
		BS-360E ⁶	/	/	—	/	/	BS-360E ⁶	/	/	—	/	/
		BS-380 ⁷	/	/	—	/	/	BS-380 ⁷	/	/	—	/	/
		BS-400 ⁸	/	/	—	/	/	BS-400 ⁸	/	/	—	/	/
		BS-480 ⁹	/	/	—	/	/	BS-480 ⁹	/	/	—	/	/
		BS-600 ¹⁰	/	/	—	/	/	BS-600 ¹⁰	/	/	—	/	/
		BS-600M ¹¹	/	/	—	/	/	BS-600M ¹¹	/	/	—	/	/
		BS-800 ¹²	/	/	—	/	/	BS-800 ¹²	/	/	—	/	/
		BS-1000M ¹³	1005	804	—	1206	67	BS-1000M ¹³	2067	1653	—	2481	138
		BS-2000 ¹⁴	/	/	—	/	/	BS-2000 ¹⁴	/	/	—	/	/
		BS-2800M ¹⁵	/	/	—	/	/	BS-2800M ¹⁵	/	/	—	/	/
		BS-120 ¹	/	/	—	/	/	BS-120 ¹	/	/	—	/	/
		BS-200 ²	/	/	—	/	/	BS-200 ²	/	/	—	/	/

μkat/L	BS-200E ³	/	/	—	/	/	BS-200E ³	/	/	—	/	/
	BS-240 ⁴	/	/	—	/	/	BS-240 ⁴	/	/	—	/	/
	BS-240E ⁵	/	/	—	/	/	BS-240E ⁵	/	/	—	/	/
	BS-360E ⁶	/	/	—	/	/	BS-360E ⁶	/	/	—	/	/
	BS-380 ⁷	/	/	—	/	/	BS-380 ⁷	/	/	—	/	/
	BS-400 ⁸	/	/	—	/	/	BS-400 ⁸	/	/	—	/	/
	BS-480 ⁹	/	/	—	/	/	BS-480 ⁹	/	/	—	/	/
	BS-600 ¹⁰	/	/	—	/	/	BS-600 ¹⁰	/	/	—	/	/
	BS-600M ¹¹	/	/	—	/	/	BS-600M ¹¹	/	/	—	/	/
	BS-800 ¹²	/	/	—	/	/	BS-800 ¹²	/	/	—	/	/
	BS-1000M ¹³	16.8	13.4	—	20.1	1.1	BS-1000M ¹³	34.5	27.6	—	41.4	2.3
	BS-2000 ¹⁴	/	/	—	/	/	BS-2000 ¹⁴	/	/	—	/	/
	BS-2800M ¹⁵	/	/	—	/	/	BS-2800M ¹⁵	/	/	—	/	/

mindray

G6PD Control

For use on: BS-1000M, BS-2000,
BS-2800M, BS-600M,
BS-800

LOT 057523006

 2025-02-10

