

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-120: BS-120, BS-130, BS-180, BS-190;

10.BS-480: BS-480, BS-490;

2.BS-200: BS-200, BS-220, BS-330, BS-350;

11.BS-600: BS-600, BS-620;

3.BS-200E: BS-200E, BS-220E,;

12.BS-600M: BS-600M, BS-620M;

4.BS-240: BS-230, BS-240;

13.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

5.BS-240E: BS240E, BS240Pro;

14.BS-1000M: BS-1000M, BS-1100M;

6.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-

15.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with "W9-" and software version starts with "35.00");

16.BS-2800M:BS-2600M, BS-2800M;

7.BS-380: BS-380, BS-390;

17.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction.

8.BS-400: BS-400, BS-420;

9.BS-430: BS-430, BS-450,BS-460;

LOT : 057525003

: 2027-01-06

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value $\pm$ 3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа $\pm$ 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise $\pm$ 3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo $\pm$ 3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione $\pm$ 3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri $\pm$ 3SD)

G6PD

English	Glucose-6-Phosphate Dehydrogenase
Русский	глюкозо-6-фосфатдегидрогеназа
Português	Glicose-6-fosfato desidrogenase
Español	Glucosa-6-fosfato deshidrogenasa
Italiano	Glucosio-6-fosfato deidrogenasi
Türkçe	Glikoz 6 Fosfat Dehidrogenaz

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
G6PD	U/L	BS-120 ¹	1112	887	—	1337	75	BS-120 ¹	2233	1783	—	2683	150
		BS-200 ²	1147	916	—	1378	77	BS-200 ²	2276	1820	—	2732	152
		BS-200E ³	996	795	—	1197	67	BS-200E ³	2081	1664	—	2498	139
		BS-240 ⁴	1036	829	—	1243	69	BS-240 ⁴	2081	1664	—	2498	139
		BS-240E ⁵	1098	876	—	1320	74	BS-240E ⁵	2186	1748	—	2624	146
		BS-360E ⁶	1033	826	—	1240	69	BS-360E ⁶	2081	1664	—	2498	139
		BS-380 ⁷	1079	863	—	1295	72	BS-380 ⁷	2280	1821	—	2739	153
		BS-400 ⁸	1053	840	—	1266	71	BS-400 ⁸	2202	1758	—	2646	148
		BS-430 ⁹	1102	880	—	1324	74	BS-430 ⁹	2199	1758	—	2640	147
		BS-480 ¹⁰	1163	929	—	1397	78	BS-480 ¹⁰	2280	1821	—	2739	153
		BS-600 ¹¹	1163	929	—	1397	78	BS-600 ¹¹	2280	1821	—	2739	153
		BS-600M ¹²	1163	929	—	1397	78	BS-600M ¹²	2280	1821	—	2739	153
		BS-800 ¹³	1114	889	—	1339	75	BS-800 ¹³	2202	1758	—	2646	148
		BS-1000M ¹⁴	1050	840	—	1260	70	BS-1000M ¹⁴	2102	1679	—	2525	141
		BS-2000 ¹⁵	1164	930	—	1398	78	BS-2000 ¹⁵	2272	1816	—	2728	152
		BS-2800M ¹⁶	1163	929	—	1397	78	BS-2800M ¹⁶	2280	1821	—	2739	153

μkat/L	BS-120¹	18.6	14.8	—	22.3	1.3	BS-120¹	37.3	29.8	—	44.8	2.5
	BS-200²	19.2	15.3	—	23.0	1.3	BS-200²	38.0	30.4	—	45.6	2.5
	BS-200E³	16.6	13.3	—	20.0	1.1	BS-200E³	34.8	27.8	—	41.7	2.3
	BS-240⁴	17.3	13.8	—	20.8	1.2	BS-240⁴	34.8	27.8	—	41.7	2.3
	BS-240E⁵	18.3	14.6	—	22.0	1.2	BS-240E⁵	36.5	29.2	—	43.8	2.4
	BS-360E⁶	17.3	13.8	—	20.7	1.2	BS-360E⁶	34.8	27.8	—	41.7	2.3
	BS-380⁷	18.0	14.4	—	21.6	1.2	BS-380⁷	38.1	30.4	—	45.7	2.6
	BS-400⁸	17.6	14.0	—	21.1	1.2	BS-400⁸	36.8	29.4	—	44.2	2.5
	BS-430⁹	18.4	14.7	—	22.1	1.2	BS-430⁹	36.7	29.4	—	44.1	2.5
	BS-480¹⁰	19.4	15.5	—	23.3	1.3	BS-480¹⁰	38.1	30.4	—	45.7	2.6
	BS-600¹¹	19.4	15.5	—	23.3	1.3	BS-600¹¹	38.1	30.4	—	45.7	2.6
	BS-600M¹²	19.4	15.5	—	23.3	1.3	BS-600M¹²	38.1	30.4	—	45.7	2.6
	BS-800¹³	18.6	14.8	—	22.4	1.3	BS-800¹³	36.8	29.4	—	44.2	2.5
	BS-1000M¹⁴	17.5	14.0	—	21.0	1.2	BS-1000M¹⁴	35.1	28.0	—	42.2	2.4
	BS-2000¹⁵	19.4	15.5	—	23.3	1.3	BS-2000¹⁵	37.9	30.3	—	45.6	2.5
	BS-2800M¹⁶	19.4	15.5	—	23.3	1.3	BS-2800M¹⁶	38.1	30.4	—	45.7	2.6

mindray

G6PD Control

For use on: BS-1000M, BS-2000,
BS-2800M, BS-600M,
BS-800

LOT 057525003

 2027-01-06

