

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray
Português: Planilha de valores do Sistema de Medição BS da Mindray
Español: Hoja de valores del sistema de medición Mindray BS
Italiano: Scheda dei valori del sistema di misurazione BS di Mindray
Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.
Português: A dados de cada grupo é a mesma.
Italiano: I dati di ogni gruppo è la stessa.
1.BS-120: BS-120, BS-130, BS-180, BS-190;
2.BS-200: BS-200, BS-220, BS-330, BS-350;
3.BS-200E: BS-200E, BS-220E,;
4.BS-240: BS-230, BS-240;
5.BS-240E: BS240E, BS240Pro;
6.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00");
7.BS-380: BS-380, BS-390;
8.BS-400: BS-400, BS-420;
9.BS-430: BS-430, BS-450,BS-460;

Русский: Данные совпадают во всех группах.
Español: los datos de cada grupo es la misma.
Türkçe: her grubun veri aynıdır.
10.BS-480: BS-480, BS-490;
11.BS-600: BS-600, BS-620;
12.BS-600M: BS-600M, BS-620M;
13.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;
14.BS-1000M: BS-1000M, BS-1100M;
15.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;
16.BS-2800M:BS-2600M, BS-2800M;
17.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction.

LOT : 057625004

2026-11-20

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

β-HB	
English	β-Hydroxybutyrate
Русский	β-гидроксибутирата
Português	β-Hidroxibutirato
Español	β-hidroxibutirato
Italiano	β-idrossibutirato
Türkçe	β-Hidroksibótrat

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
β-HB	mmol/L	BS-120 ¹	0.285	0.228	—	0.342	0.019	BS-120 ¹	1.13	0.89	—	1.37	0.08
		BS-200 ²	0.285	0.228	—	0.342	0.019	BS-200 ²	1.13	0.89	—	1.37	0.08
		BS-200E ³	0.289	0.232	—	0.346	0.019	BS-200E ³	1.14	0.90	—	1.38	0.08
		BS-240 ⁴	0.285	0.228	—	0.342	0.019	BS-240 ⁴	1.13	0.89	—	1.37	0.08
		BS-240E ⁵	0.289	0.232	—	0.346	0.019	BS-240E ⁵	1.14	0.90	—	1.38	0.08
		BS-360E ⁶	0.289	0.232	—	0.346	0.019	BS-360E ⁶	1.14	0.90	—	1.38	0.08
		BS-380 ⁷	0.291	0.234	—	0.348	0.019	BS-380 ⁷	1.15	0.91	—	1.39	0.08
		BS-400 ⁸	0.291	0.234	—	0.348	0.019	BS-400 ⁸	1.15	0.91	—	1.39	0.08
		BS-430 ⁹	0.289	0.232	—	0.346	0.019	BS-430 ⁹	1.14	0.90	—	1.38	0.08
		BS-480 ¹⁰	0.291	0.234	—	0.348	0.019	BS-480 ¹⁰	1.15	0.91	—	1.39	0.08
		BS-600 ¹¹	0.291	0.234	—	0.348	0.019	BS-600 ¹¹	1.15	0.91	—	1.39	0.08
		BS-600M ¹²	0.291	0.234	—	0.348	0.019	BS-600M ¹²	1.15	0.91	—	1.39	0.08
		BS-800 ¹³	0.289	0.232	—	0.346	0.019	BS-800 ¹³	1.14	0.90	—	1.38	0.08
		BS-1000M ¹⁴	0.291	0.234	—	0.348	0.019	BS-1000M ¹⁴	1.15	0.91	—	1.39	0.08
		BS-2000 ¹⁵	0.291	0.234	—	0.348	0.019	BS-2000 ¹⁵	1.14	0.90	—	1.38	0.08
		BS-2800M ¹⁶	0.291	0.234	—	0.348	0.019	BS-2800M ¹⁶	1.15	0.91	—	1.39	0.08

mindray
β-HB Control

For use on: BS-1000M, BS-2000,
BS-2800M, BS-600M,
BS-800

LOT 057625004

2026-11-20

