

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.	Русский: Данные совпадают во всех группах.
Português: A dados de cada grupo é a mesma.	Español: la datos de cada grupo es la misma.
Italiano: la dati di ogni gruppo è la stessa.	Türkçe: her grubun veri aynıdır.
1.BS-120: BS-120, BS-130, BS-180, BS-190;	10.BS-480: BS-480, BS-490;
2.BS-200: BS-200, BS-220, BS-330, BS-350;	11.BS-600: BS-600, BS-620;
3.BS-200E: BS-200E, BS-220E,;	12.BS-600M: BS-600M, BS-620M;
4.BS-240: BS-230, BS-240;	13.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;
5.BS-240E: BS240E, BS240Pro;	14.BS-1000M: BS-1000M, BS-1100M;
6.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00");	15.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;
7.BS-380: BS-380, BS-390;	16.BS-2800M:BS-2600M, BS-2800M;
8.BS-400: BS-400, BS-420;	17.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction.
9.BS-430: BS-430, BS-450,BS-460;	

LOT : 057723007

: 2024-11-13

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

ACE	
English	Angiotensin Converting Enzyme
Русский	Ангиотензин-превращающего фермента
Português	Enzima de conversão da angiotensina
Español	Enzima convertidora de angiotensina
Italiano	Enzima convertitore dell'angiotensina
Türkçe	Anjiotensin Dönüştürücü Enzim

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
ACE	U/L	BS-120 ¹	/	/	—	/	/	BS-120 ¹	/	/	—	/	/
		BS-200 ²	/	/	—	/	/	BS-200 ²	/	/	—	/	/
		BS-200E ³	/	/	—	/	/	BS-200E ³	/	/	—	/	/
		BS-240 ⁴	/	/	—	/	/	BS-240 ⁴	/	/	—	/	/
		BS-240E ⁵	/	/	—	/	/	BS-240E ⁵	/	/	—	/	/
		BS-360E ⁶	/	/	—	/	/	BS-360E ⁶	/	/	—	/	/
		BS-380 ⁷	/	/	—	/	/	BS-380 ⁷	/	/	—	/	/
		BS-400 ⁸	/	/	—	/	/	BS-400 ⁸	/	/	—	/	/
		BS-430 ⁹	/	/	—	/	/	BS-430 ⁹	/	/	—	/	/
		BS-480 ¹⁰	/	/	—	/	/	BS-480 ¹⁰	/	/	—	/	/
		BS-600 ¹¹	/	/	—	/	/	BS-600 ¹¹	/	/	—	/	/
		BS-600M ¹²	/	/	—	/	/	BS-600M ¹²	/	/	—	/	/
		BS-800 ¹³	/	/	—	/	/	BS-800 ¹³	/	/	—	/	/
		BS-1000M ¹⁴	48.2	36.2	—	60.2	4.0	BS-1000M ¹⁴	142	106	—	178	12
		BS-2000 ¹⁵	/	/	—	/	/	BS-2000 ¹⁵	/	/	—	/	/

BS-2800M¹⁶ / / — / / BS-2800M¹⁶ / / — / /

μkat/L	BS-120 ¹	/	/	—	/	/	BS-120 ¹	/	/	—	/	/
	BS-200 ²	/	/	—	/	/	BS-200 ²	/	/	—	/	/
	BS-200E ³	/	/	—	/	/	BS-200E ³	/	/	—	/	/
	BS-240 ⁴	/	/	—	/	/	BS-240 ⁴	/	/	—	/	/
	BS-240E ⁵	/	/	—	/	/	BS-240E ⁵	/	/	—	/	/
	BS-360E ⁶	/	/	—	/	/	BS-360E ⁶	/	/	—	/	/
	BS-380 ⁷	/	/	—	/	/	BS-380 ⁷	/	/	—	/	/
	BS-400 ⁸	/	/	—	/	/	BS-400 ⁸	/	/	—	/	/
	BS-430 ⁹	/	/	—	/	/	BS-430 ⁹	/	/	—	/	/
	BS-480 ¹⁰	/	/	—	/	/	BS-480 ¹⁰	/	/	—	/	/
	BS-600 ¹¹	/	/	—	/	/	BS-600 ¹¹	/	/	—	/	/
	BS-600M ¹²	/	/	—	/	/	BS-600M ¹²	/	/	—	/	/
	BS-800 ¹³	/	/	—	/	/	BS-800 ¹³	/	/	—	/	/
	BS-1000M ¹⁴	0.805	0.605	—	1.005	0.067	BS-1000M ¹⁴	2.37	1.77	—	2.97	0.20
	BS-2000 ¹⁵	/	/	—	/	/	BS-2000 ¹⁵	/	/	—	/	/
	BS-2800M ¹⁶	/	/	—	/	/	BS-2800M ¹⁶	/	/	—	/	/

mindray

ACE Control

For use on: BS-1000M, BS-2000,
BS-2800M, BS-600M,
BS-800

LOT 057723007

 2024-11-13

