

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-120: BS-120, BS-130, BS-180, BS-190;

10.BS-480: BS-480, BS-490;

2.BS-200: BS-200, BS-220, BS-330, BS-350;

11.BS-600: BS-600, BS-620;

3.BS-200E: BS-200E, BS-220E,;

12.BS-600M: BS-600M, BS-620M;

4.BS-240: BS-230, BS-240;

13.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

5.BS-240E: BS240E, BS240Pro;

14.BS-1000M: BS-1000M, BS-1100M;

6.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-

15.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with "W9-" and software version starts with "35.00");

16.BS-2800M:BS-2600M, BS-2800M;

7.BS-380: BS-380, BS-390;

17.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction.

8.BS-400: BS-400, BS-420;

9.BS-430: BS-430, BS-450,BS-460;

LOT : 057724005

: 2025-04-07

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value $\pm$ 3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа $\pm$ 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise $\pm$ 3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo $\pm$ 3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione $\pm$ 3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri $\pm$ 3SD)

ACE

English	Angiotensin Converting Enzyme
Русский	Ангиотензин-превращающего фермента
Português	Enzima de conversão da angiotensina
Español	Enzima convertidora de angiotensina
Italiano	Enzima convertitore dell'angiotensina
Türkçe	Anjiotensin Dönüştürücü Enzim

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
ACE	U/L	BS-120 ¹	48.9	36.6	—	61.2	4.1	BS-120 ¹	126	96	—	156	10
		BS-200 ²	49.8	37.5	—	62.1	4.1	BS-200 ²	130	97	—	163	11
		BS-200E ³	43.9	33.1	—	54.7	3.6	BS-200E ³	115	85	—	145	10
		BS-240 ⁴	48.9	36.6	—	61.2	4.1	BS-240 ⁴	126	96	—	156	10
		BS-240E ⁵	40.7	30.5	—	50.9	3.4	BS-240E ⁵	107	80	—	134	9
		BS-360E ⁶	40.7	30.5	—	50.9	3.4	BS-360E ⁶	107	80	—	134	9
		BS-380 ⁷	43.6	32.8	—	54.4	3.6	BS-380 ⁷	109	82	—	136	9
		BS-400 ⁸	43.6	32.8	—	54.4	3.6	BS-400 ⁸	115	85	—	145	10
		BS-430 ⁹	40.7	30.5	—	50.9	3.4	BS-430 ⁹	107	80	—	134	9
		BS-480 ¹⁰	40.7	30.5	—	50.9	3.4	BS-480 ¹⁰	107	80	—	134	9
		BS-600 ¹¹	40.7	30.5	—	50.9	3.4	BS-600 ¹¹	107	80	—	134	9
		BS-600M ¹²	44.8	33.7	—	55.9	3.7	BS-600M ¹²	119	89	—	149	10
		BS-800 ¹³	40.7	30.5	—	50.9	3.4	BS-800 ¹³	107	80	—	134	9
		BS-1000M ¹⁴	43.6	32.8	—	54.4	3.6	BS-1000M ¹⁴	115	85	—	145	10
		BS-2000 ¹⁵	42.6	32.1	—	53.1	3.5	BS-2000 ¹⁵	113	86	—	140	9
		BS-2800M ¹⁶	43.6	32.8	—	54.4	3.6	BS-2800M ¹⁶	115	85	—	145	10

μkat/L	BS-120¹	0.817	0.611	—	1.022	0.068	BS-120¹	2.10	1.60	—	2.61	0.17
	BS-200²	0.832	0.626	—	1.037	0.068	BS-200²	2.17	1.62	—	2.72	0.18
	BS-200E³	0.733	0.553	—	0.913	0.060	BS-200E³	1.92	1.42	—	2.42	0.17
	BS-240⁴	0.817	0.611	—	1.022	0.068	BS-240⁴	2.10	1.60	—	2.61	0.17
	BS-240E⁵	0.680	0.509	—	0.850	0.057	BS-240E⁵	1.79	1.34	—	2.24	0.15
	BS-360E⁶	0.680	0.509	—	0.850	0.057	BS-360E⁶	1.79	1.34	—	2.24	0.15
	BS-380⁷	0.728	0.548	—	0.908	0.060	BS-380⁷	1.82	1.37	—	2.27	0.15
	BS-400⁸	0.728	0.548	—	0.908	0.060	BS-400⁸	1.92	1.42	—	2.42	0.17
	BS-430⁹	0.680	0.509	—	0.850	0.057	BS-430⁹	1.79	1.34	—	2.24	0.15
	BS-480¹⁰	0.680	0.509	—	0.850	0.057	BS-480¹⁰	1.79	1.34	—	2.24	0.15
	BS-600¹¹	0.680	0.509	—	0.850	0.057	BS-600¹¹	1.79	1.34	—	2.24	0.15
	BS-600M¹²	0.748	0.563	—	0.934	0.062	BS-600M¹²	1.99	1.49	—	2.49	0.17
	BS-800¹³	0.680	0.509	—	0.850	0.057	BS-800¹³	1.79	1.34	—	2.24	0.15
	BS-1000M¹⁴	0.728	0.548	—	0.908	0.060	BS-1000M¹⁴	1.92	1.42	—	2.42	0.17
	BS-2000¹⁵	0.711	0.536	—	0.887	0.058	BS-2000¹⁵	1.89	1.44	—	2.34	0.15
	BS-2800M¹⁶	0.728	0.548	—	0.908	0.060	BS-2800M¹⁶	1.92	1.42	—	2.42	0.17

mindray

ACE Control

For use on: BS-1000M, BS-2000,
BS-2800M, BS-600M,
BS-800

LOT 057724005

 **2025-04-07**

