

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-120: BS-120, BS-130, BS-180, BS-190;

10.BS-480: BS-480, BS-490;

2.BS-200: BS-200, BS-220, BS-330, BS-350;

11.BS-600: BS-600, BS-620;

3.BS-200E: BS-200E, BS-220E,;

12.BS-600M: BS-600M, BS-620M;

4.BS-240: BS-230, BS-240;

13.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

5.BS-240E: BS240E, BS240Pro;

14.BS-1000M: BS-1000M, BS-1100M;

6.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-

15.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with "W9-" and software version starts with "35.00");

7.BS-380: BS-380, BS-390;

16.BS-2800M:BS-2600M, BS-2800M;

8.BS-400: BS-400, BS-420;

17.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction.

9.BS-430: BS-430, BS-450,BS-460;

LOT : 057724006

: 2025-09-11

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value $\pm$ 3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа $\pm$ 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise $\pm$ 3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo $\pm$ 3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione $\pm$ 3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri $\pm$ 3SD)

ACE

English	Angiotensin Converting Enzyme
Русский	Ангиотензин-превращающего фермента
Português	Enzima de conversão da angiotensina
Español	Enzima convertidora de angiotensina
Italiano	Enzima convertitore dell'angiotensina
Türkçe	Anjiotensin Dönüştürücü Enzim

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
ACE	U/L	BS-120 ¹	55.5	41.7	—	69.3	4.6	BS-120 ¹	132	99	—	165	11
		BS-200 ²	56.5	42.4	—	70.6	4.7	BS-200 ²	137	104	—	170	11
		BS-200E ³	49.7	37.4	—	62.0	4.1	BS-200E ³	113	86	—	140	9
		BS-240 ⁴	55.5	41.7	—	69.3	4.6	BS-240 ⁴	132	99	—	165	11
		BS-240E ⁵	46.1	34.7	—	57.5	3.8	BS-240E ⁵	105	78	—	132	9
		BS-360E ⁶	46.1	34.7	—	57.5	3.8	BS-360E ⁶	105	78	—	132	9
		BS-380 ⁷	47.2	35.5	—	58.9	3.9	BS-380 ⁷	109	82	—	136	9
		BS-400 ⁸	47.2	35.5	—	58.9	3.9	BS-400 ⁸	115	85	—	145	10
		BS-430 ⁹	46.1	34.7	—	57.5	3.8	BS-430 ⁹	105	78	—	132	9
		BS-480 ¹⁰	46.1	34.7	—	57.5	3.8	BS-480 ¹⁰	105	78	—	132	9
		BS-600 ¹¹	46.1	34.7	—	57.5	3.8	BS-600 ¹¹	105	78	—	132	9
		BS-600M ¹²	49.5	37.2	—	61.8	4.1	BS-600M ¹²	118	88	—	148	10
		BS-800 ¹³	46.1	34.7	—	57.5	3.8	BS-800 ¹³	105	78	—	132	9
		BS-1000M ¹⁴	47.2	35.5	—	58.9	3.9	BS-1000M ¹⁴	115	85	—	145	10
		BS-2000 ¹⁵	45.9	34.5	—	57.3	3.8	BS-2000 ¹⁵	111	84	—	138	9
		BS-2800M ¹⁶	47.2	35.5	—	58.9	3.9	BS-2800M ¹⁶	115	85	—	145	10

μkat/L	BS-120¹	0.927	0.696	—	1.157	0.077	BS-120¹	2.20	1.65	—	2.76	0.18
	BS-200²	0.944	0.708	—	1.179	0.078	BS-200²	2.29	1.74	—	2.84	0.18
	BS-200E³	0.830	0.625	—	1.035	0.068	BS-200E³	1.89	1.44	—	2.34	0.15
	BS-240⁴	0.927	0.696	—	1.157	0.077	BS-240⁴	2.20	1.65	—	2.76	0.18
	BS-240E⁵	0.770	0.579	—	0.960	0.063	BS-240E⁵	1.75	1.30	—	2.20	0.15
	BS-360E⁶	0.770	0.579	—	0.960	0.063	BS-360E⁶	1.75	1.30	—	2.20	0.15
	BS-380⁷	0.788	0.593	—	0.984	0.065	BS-380⁷	1.82	1.37	—	2.27	0.15
	BS-400⁸	0.788	0.593	—	0.984	0.065	BS-400⁸	1.92	1.42	—	2.42	0.17
	BS-430⁹	0.770	0.579	—	0.960	0.063	BS-430⁹	1.75	1.30	—	2.20	0.15
	BS-480¹⁰	0.770	0.579	—	0.960	0.063	BS-480¹⁰	1.75	1.30	—	2.20	0.15
	BS-600¹¹	0.770	0.579	—	0.960	0.063	BS-600¹¹	1.75	1.30	—	2.20	0.15
	BS-600M¹²	0.827	0.621	—	1.032	0.068	BS-600M¹²	1.97	1.47	—	2.47	0.17
	BS-800¹³	0.770	0.579	—	0.960	0.063	BS-800¹³	1.75	1.30	—	2.20	0.15
	BS-1000M¹⁴	0.788	0.593	—	0.984	0.065	BS-1000M¹⁴	1.92	1.42	—	2.42	0.17
	BS-2000¹⁵	0.767	0.576	—	0.957	0.063	BS-2000¹⁵	1.85	1.40	—	2.30	0.15
	BS-2800M¹⁶	0.788	0.593	—	0.984	0.065	BS-2800M¹⁶	1.92	1.42	—	2.42	0.17

mindray

ACE Control

For use on: BS-1000M, BS-2000,
BS-2800M, BS-600M,
BS-800

LOT 057724006

 2025-09-11

