

## Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-120: BS-120, BS-130, BS-180, BS-190;

10.BS-480: BS-480, BS-490;

2.BS-200: BS-200, BS-220, BS-330, BS-350;

11.BS-600: BS-600, BS-620;

3.BS-200E: BS-200E, BS-220E,;

12.BS-600M: BS-600M, BS-620M;

4.BS-240: BS-230, BS-240;

13.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

5.BS-240E: BS240E, BS240Pro;

14.BS-1000M: BS-1000M, BS-1100M;

6.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-

15.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with "W9-" and software version starts with "35.00");

16.BS-2800M:BS-2600M, BS-2800M;

7.BS-380: BS-380, BS-390;

17.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction.

8.BS-400: BS-400, BS-420;

9.BS-430: BS-430, BS-450,BS-460;

LOT : 057724007

: 2025-09-11

| English   | Abbreviated name         | Model   | Unit    | Assay Value        | Range (Assay Value $\pm$ 3SD)                       |
|-----------|--------------------------|---------|---------|--------------------|-----------------------------------------------------|
| Русский   | сокращенное наименование | модель  | Прибор  | Результат анализа  | Диапазон(результат анализа $\pm$ 3CO)               |
| Português | Nome abreviado           | Modelo  | Unidade | Valores da análise | Faixa(Valores da análise $\pm$ 3SD)                 |
| Español   | nombre abreviado         | modelo  | Unidad  | Valor de ensayo    | Rango(Valor de ensayo $\pm$ 3SD)                    |
| Italiano  | abbreviazione            | modelli | Unità   | Valori di dosaggio | Intervallo<br>(valore di concentrazione $\pm$ 3 DS) |
| Türkçe    | kısaltılmış ad           | model   | Ünite   | Tayin Değeri       | Aralık (Tayin Değeri $\pm$ 3SD)                     |

## ACE

|           |                                       |
|-----------|---------------------------------------|
| English   | Angiotensin Converting Enzyme         |
| Русский   | Ангиотензин-превращающего фермента    |
| Português | Enzima de conversão da angiotensina   |
| Español   | Enzima convertidora de angiotensina   |
| Italiano  | Enzima convertitore dell'angiotensina |
| Türkçe    | Anjiotensin Dönüştürücü Enzim         |

| Abbreviated name | Unit | Level: L/N             |             |                         |   |      | Level: H/P |                        |                         |     |      |     |    |
|------------------|------|------------------------|-------------|-------------------------|---|------|------------|------------------------|-------------------------|-----|------|-----|----|
|                  |      | Model                  | Assay Value | Range (Assay Value±3SD) |   | 1 SD | Model      | Assay Value            | Range (Assay Value±3SD) |     | 1 SD |     |    |
| ACE              | U/L  | BS-120 <sup>1</sup>    | 55.5        | 41.7                    | — | 69.3 | 4.6        | BS-120 <sup>1</sup>    | 132                     | 99  | —    | 165 | 11 |
|                  |      | BS-200 <sup>2</sup>    | 56.5        | 42.4                    | — | 70.6 | 4.7        | BS-200 <sup>2</sup>    | 137                     | 104 | —    | 170 | 11 |
|                  |      | BS-200E <sup>3</sup>   | 49.7        | 37.4                    | — | 62.0 | 4.1        | BS-200E <sup>3</sup>   | 113                     | 86  | —    | 140 | 9  |
|                  |      | BS-240 <sup>4</sup>    | 55.5        | 41.7                    | — | 69.3 | 4.6        | BS-240 <sup>4</sup>    | 132                     | 99  | —    | 165 | 11 |
|                  |      | BS-240E <sup>5</sup>   | 46.1        | 34.7                    | — | 57.5 | 3.8        | BS-240E <sup>5</sup>   | 105                     | 78  | —    | 132 | 9  |
|                  |      | BS-360E <sup>6</sup>   | 46.1        | 34.7                    | — | 57.5 | 3.8        | BS-360E <sup>6</sup>   | 105                     | 78  | —    | 132 | 9  |
|                  |      | BS-380 <sup>7</sup>    | 47.2        | 35.5                    | — | 58.9 | 3.9        | BS-380 <sup>7</sup>    | 109                     | 82  | —    | 136 | 9  |
|                  |      | BS-400 <sup>8</sup>    | 47.2        | 35.5                    | — | 58.9 | 3.9        | BS-400 <sup>8</sup>    | 115                     | 85  | —    | 145 | 10 |
|                  |      | BS-430 <sup>9</sup>    | 46.1        | 34.7                    | — | 57.5 | 3.8        | BS-430 <sup>9</sup>    | 105                     | 78  | —    | 132 | 9  |
|                  |      | BS-480 <sup>10</sup>   | 46.1        | 34.7                    | — | 57.5 | 3.8        | BS-480 <sup>10</sup>   | 105                     | 78  | —    | 132 | 9  |
|                  |      | BS-600 <sup>11</sup>   | 46.1        | 34.7                    | — | 57.5 | 3.8        | BS-600 <sup>11</sup>   | 105                     | 78  | —    | 132 | 9  |
|                  |      | BS-600M <sup>12</sup>  | 49.5        | 37.2                    | — | 61.8 | 4.1        | BS-600M <sup>12</sup>  | 118                     | 88  | —    | 148 | 10 |
|                  |      | BS-800 <sup>13</sup>   | 46.1        | 34.7                    | — | 57.5 | 3.8        | BS-800 <sup>13</sup>   | 105                     | 78  | —    | 132 | 9  |
|                  |      | BS-1000M <sup>14</sup> | 47.2        | 35.5                    | — | 58.9 | 3.9        | BS-1000M <sup>14</sup> | 115                     | 85  | —    | 145 | 10 |
|                  |      | BS-2000 <sup>15</sup>  | 45.9        | 34.5                    | — | 57.3 | 3.8        | BS-2000 <sup>15</sup>  | 111                     | 84  | —    | 138 | 9  |
|                  |      | BS-2800M <sup>16</sup> | 47.2        | 35.5                    | — | 58.9 | 3.9        | BS-2800M <sup>16</sup> | 115                     | 85  | —    | 145 | 10 |

|        |                              |       |       |   |       |       |                              |      |      |   |      |      |
|--------|------------------------------|-------|-------|---|-------|-------|------------------------------|------|------|---|------|------|
| μkat/L | <b>BS-120<sup>1</sup></b>    | 0.927 | 0.696 | — | 1.157 | 0.077 | <b>BS-120<sup>1</sup></b>    | 2.20 | 1.65 | — | 2.76 | 0.18 |
|        | <b>BS-200<sup>2</sup></b>    | 0.944 | 0.708 | — | 1.179 | 0.078 | <b>BS-200<sup>2</sup></b>    | 2.29 | 1.74 | — | 2.84 | 0.18 |
|        | <b>BS-200E<sup>3</sup></b>   | 0.830 | 0.625 | — | 1.035 | 0.068 | <b>BS-200E<sup>3</sup></b>   | 1.89 | 1.44 | — | 2.34 | 0.15 |
|        | <b>BS-240<sup>4</sup></b>    | 0.927 | 0.696 | — | 1.157 | 0.077 | <b>BS-240<sup>4</sup></b>    | 2.20 | 1.65 | — | 2.76 | 0.18 |
|        | <b>BS-240E<sup>5</sup></b>   | 0.770 | 0.579 | — | 0.960 | 0.063 | <b>BS-240E<sup>5</sup></b>   | 1.75 | 1.30 | — | 2.20 | 0.15 |
|        | <b>BS-360E<sup>6</sup></b>   | 0.770 | 0.579 | — | 0.960 | 0.063 | <b>BS-360E<sup>6</sup></b>   | 1.75 | 1.30 | — | 2.20 | 0.15 |
|        | <b>BS-380<sup>7</sup></b>    | 0.788 | 0.593 | — | 0.984 | 0.065 | <b>BS-380<sup>7</sup></b>    | 1.82 | 1.37 | — | 2.27 | 0.15 |
|        | <b>BS-400<sup>8</sup></b>    | 0.788 | 0.593 | — | 0.984 | 0.065 | <b>BS-400<sup>8</sup></b>    | 1.92 | 1.42 | — | 2.42 | 0.17 |
|        | <b>BS-430<sup>9</sup></b>    | 0.770 | 0.579 | — | 0.960 | 0.063 | <b>BS-430<sup>9</sup></b>    | 1.75 | 1.30 | — | 2.20 | 0.15 |
|        | <b>BS-480<sup>10</sup></b>   | 0.770 | 0.579 | — | 0.960 | 0.063 | <b>BS-480<sup>10</sup></b>   | 1.75 | 1.30 | — | 2.20 | 0.15 |
|        | <b>BS-600<sup>11</sup></b>   | 0.770 | 0.579 | — | 0.960 | 0.063 | <b>BS-600<sup>11</sup></b>   | 1.75 | 1.30 | — | 2.20 | 0.15 |
|        | <b>BS-600M<sup>12</sup></b>  | 0.827 | 0.621 | — | 1.032 | 0.068 | <b>BS-600M<sup>12</sup></b>  | 1.97 | 1.47 | — | 2.47 | 0.17 |
|        | <b>BS-800<sup>13</sup></b>   | 0.770 | 0.579 | — | 0.960 | 0.063 | <b>BS-800<sup>13</sup></b>   | 1.75 | 1.30 | — | 2.20 | 0.15 |
|        | <b>BS-1000M<sup>14</sup></b> | 0.788 | 0.593 | — | 0.984 | 0.065 | <b>BS-1000M<sup>14</sup></b> | 1.92 | 1.42 | — | 2.42 | 0.17 |
|        | <b>BS-2000<sup>15</sup></b>  | 0.767 | 0.576 | — | 0.957 | 0.063 | <b>BS-2000<sup>15</sup></b>  | 1.85 | 1.40 | — | 2.30 | 0.15 |
|        | <b>BS-2800M<sup>16</sup></b> | 0.788 | 0.593 | — | 0.984 | 0.065 | <b>BS-2800M<sup>16</sup></b> | 1.92 | 1.42 | — | 2.42 | 0.17 |

# **mindray**

## **ACE Control**

For use on: BS-1000M, BS-2000,  
BS-2800M, BS-600M,  
BS-800

**LOT** 057724007

 2025-09-11

