

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-120: BS-120, BS-130, BS-180, BS-190;

10.BS-480: BS-480, BS-490;

2.BS-200: BS-200, BS-220, BS-330, BS-350;

11.BS-600: BS-600, BS-620;

3.BS-200E: BS-200E, BS-220E,;

12.BS-600M: BS-600M, BS-620M;

4.BS-240: BS-230, BS-240;

13.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

5.BS-240E: BS240E, BS240Pro;

14.BS-1000M: BS-1000M, BS-1100M;

6.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-

15.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with "W9-" and software version starts with "35.00");

16.BS-2800M:BS-2600M, BS-2800M;

7.BS-380: BS-380, BS-390;

17.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction.

8.BS-400: BS-400, BS-420;

9.BS-430: BS-430, BS-450,BS-460;

LOT : 057724009

: 2025-12-24

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value $\pm$ 3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа $\pm$ 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise $\pm$ 3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo $\pm$ 3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione $\pm$ 3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri $\pm$ 3SD)

ACE

English	Angiotensin Converting Enzyme
Русский	Ангиотензин-превращающего фермента
Português	Enzima de conversão da angiotensina
Español	Enzima convertidora de angiotensina
Italiano	Enzima convertitore dell'angiotensina
Türkçe	Anjiotensin Dönüştürücü Enzim

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
ACE	U/L	BS-120 ¹	51.7	38.8	—	64.6	4.3	BS-120 ¹	128	95	—	161	11
		BS-200 ²	52.6	39.4	—	65.8	4.4	BS-200 ²	133	100	—	166	11
		BS-200E ³	48.5	36.5	—	60.5	4.0	BS-200E ³	116	86	—	146	10
		BS-240 ⁴	51.7	38.8	—	64.6	4.3	BS-240 ⁴	128	95	—	161	11
		BS-240E ⁵	45.0	33.9	—	56.1	3.7	BS-240E ⁵	108	81	—	135	9
		BS-360E ⁶	45.0	33.9	—	56.1	3.7	BS-360E ⁶	108	81	—	135	9
		BS-380 ⁷	46.2	34.8	—	57.6	3.8	BS-380 ⁷	112	85	—	139	9
		BS-400 ⁸	46.2	34.8	—	57.6	3.8	BS-400 ⁸	119	89	—	149	10
		BS-430 ⁹	45.0	33.9	—	56.1	3.7	BS-430 ⁹	108	81	—	135	9
		BS-480 ¹⁰	45.0	33.9	—	56.1	3.7	BS-480 ¹⁰	108	81	—	135	9
		BS-600 ¹¹	45.0	33.9	—	56.1	3.7	BS-600 ¹¹	108	81	—	135	9
		BS-600M ¹²	47.1	35.4	—	58.8	3.9	BS-600M ¹²	117	87	—	147	10
		BS-800 ¹³	45.0	33.9	—	56.1	3.7	BS-800 ¹³	108	81	—	135	9
		BS-1000M ¹⁴	46.2	34.8	—	57.6	3.8	BS-1000M ¹⁴	119	89	—	149	10
		BS-2000 ¹⁵	45.7	34.3	—	57.1	3.8	BS-2000 ¹⁵	115	85	—	145	10
		BS-2800M ¹⁶	46.2	34.8	—	57.6	3.8	BS-2800M ¹⁶	119	89	—	149	10

μkat/L	BS-120¹	0.863	0.648	—	1.079	0.072	BS-120¹	2.14	1.59	—	2.69	0.18
	BS-200²	0.878	0.658	—	1.099	0.073	BS-200²	2.22	1.67	—	2.77	0.18
	BS-200E³	0.810	0.610	—	1.010	0.067	BS-200E³	1.94	1.44	—	2.44	0.17
	BS-240⁴	0.863	0.648	—	1.079	0.072	BS-240⁴	2.14	1.59	—	2.69	0.18
	BS-240E⁵	0.752	0.566	—	0.937	0.062	BS-240E⁵	1.80	1.35	—	2.25	0.15
	BS-360E⁶	0.752	0.566	—	0.937	0.062	BS-360E⁶	1.80	1.35	—	2.25	0.15
	BS-380⁷	0.772	0.581	—	0.962	0.063	BS-380⁷	1.87	1.42	—	2.32	0.15
	BS-400⁸	0.772	0.581	—	0.962	0.063	BS-400⁸	1.99	1.49	—	2.49	0.17
	BS-430⁹	0.752	0.566	—	0.937	0.062	BS-430⁹	1.80	1.35	—	2.25	0.15
	BS-480¹⁰	0.752	0.566	—	0.937	0.062	BS-480¹⁰	1.80	1.35	—	2.25	0.15
	BS-600¹¹	0.752	0.566	—	0.937	0.062	BS-600¹¹	1.80	1.35	—	2.25	0.15
	BS-600M¹²	0.787	0.591	—	0.982	0.065	BS-600M¹²	1.95	1.45	—	2.45	0.17
	BS-800¹³	0.752	0.566	—	0.937	0.062	BS-800¹³	1.80	1.35	—	2.25	0.15
	BS-1000M¹⁴	0.772	0.581	—	0.962	0.063	BS-1000M¹⁴	1.99	1.49	—	2.49	0.17
	BS-2000¹⁵	0.763	0.573	—	0.954	0.063	BS-2000¹⁵	1.92	1.42	—	2.42	0.17
	BS-2800M¹⁶	0.772	0.581	—	0.962	0.063	BS-2800M¹⁶	1.99	1.49	—	2.49	0.17

mindray

ACE Control

For use on: BS-1000M, BS-2000,
BS-2800M, BS-600M,
BS-800

LOT 057724009

 2025-12-24

