

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-120: BS-120, BS-130, BS-180, BS-190;

10.BS-480: BS-480, BS-490;

2.BS-200: BS-200, BS-220, BS-330, BS-350;

11.BS-600: BS-600, BS-620;

3.BS-200E: BS-200E, BS-220E,;

12.BS-600M: BS-600M, BS-620M;

4.BS-240: BS-230, BS-240;

13.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

5.BS-240E: BS240E, BS240Pro;

14.BS-1000M: BS-1000M, BS-1100M;

6.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-

15.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with "W9-" and software version starts with "35.00");

16.BS-2800M:BS-2600M, BS-2800M;

7.BS-380: BS-380, BS-390;

17.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction.

8.BS-400: BS-400, BS-420;

9.BS-430: BS-430, BS-450,BS-460;

LOT : 057725001

: 2026-03-10

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value $\pm$ 3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа $\pm$ 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise $\pm$ 3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo $\pm$ 3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione $\pm$ 3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri $\pm$ 3SD)

ACE

English	Angiotensin Converting Enzyme
Русский	Ангиотензин-превращающего фермента
Português	Enzima de conversão da angiotensina
Español	Enzima convertidora de angiotensina
Italiano	Enzima convertitore dell'angiotensina
Türkçe	Anjiotensin Dönüştürücü Enzim

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
ACE	U/L	BS-120 ¹	51.9	39.0	—	64.8	4.3	BS-120 ¹	133	100	—	166	11
		BS-200 ²	52.8	39.6	—	66.0	4.4	BS-200 ²	138	105	—	171	11
		BS-200E ³	49.0	36.7	—	61.3	4.1	BS-200E ³	118	88	—	148	10
		BS-240 ⁴	51.9	39.0	—	64.8	4.3	BS-240 ⁴	133	100	—	166	11
		BS-240E ⁵	45.4	34.0	—	56.8	3.8	BS-240E ⁵	110	83	—	137	9
		BS-360E ⁶	45.4	34.0	—	56.8	3.8	BS-360E ⁶	110	83	—	137	9
		BS-380 ⁷	45.9	34.5	—	57.3	3.8	BS-380 ⁷	114	87	—	141	9
		BS-400 ⁸	45.9	34.5	—	57.3	3.8	BS-400 ⁸	121	91	—	151	10
		BS-430 ⁹	45.4	34.0	—	56.8	3.8	BS-430 ⁹	110	83	—	137	9
		BS-480 ¹⁰	45.4	34.0	—	56.8	3.8	BS-480 ¹⁰	110	83	—	137	9
		BS-600 ¹¹	45.4	34.0	—	56.8	3.8	BS-600 ¹¹	110	83	—	137	9
		BS-600M ¹²	45.9	34.5	—	57.3	3.8	BS-600M ¹²	121	91	—	151	10
		BS-800 ¹³	45.4	34.0	—	56.8	3.8	BS-800 ¹³	110	83	—	137	9
		BS-1000M ¹⁴	45.9	34.5	—	57.3	3.8	BS-1000M ¹⁴	121	91	—	151	10
		BS-2000 ¹⁵	46.9	35.2	—	58.6	3.9	BS-2000 ¹⁵	121	91	—	151	10
		BS-2800M ¹⁶	45.9	34.5	—	57.3	3.8	BS-2800M ¹⁶	121	91	—	151	10

μkat/L	BS-120¹	0.867	0.651	—	1.082	0.072	BS-120¹	2.22	1.67	—	2.77	0.18
	BS-200²	0.882	0.661	—	1.102	0.073	BS-200²	2.30	1.75	—	2.86	0.18
	BS-200E³	0.818	0.613	—	1.024	0.068	BS-200E³	1.97	1.47	—	2.47	0.17
	BS-240⁴	0.867	0.651	—	1.082	0.072	BS-240⁴	2.22	1.67	—	2.77	0.18
	BS-240E⁵	0.758	0.568	—	0.949	0.063	BS-240E⁵	1.84	1.39	—	2.29	0.15
	BS-360E⁶	0.758	0.568	—	0.949	0.063	BS-360E⁶	1.84	1.39	—	2.29	0.15
	BS-380⁷	0.767	0.576	—	0.957	0.063	BS-380⁷	1.90	1.45	—	2.35	0.15
	BS-400⁸	0.767	0.576	—	0.957	0.063	BS-400⁸	2.02	1.52	—	2.52	0.17
	BS-430⁹	0.758	0.568	—	0.949	0.063	BS-430⁹	1.84	1.39	—	2.29	0.15
	BS-480¹⁰	0.758	0.568	—	0.949	0.063	BS-480¹⁰	1.84	1.39	—	2.29	0.15
	BS-600¹¹	0.758	0.568	—	0.949	0.063	BS-600¹¹	1.84	1.39	—	2.29	0.15
	BS-600M¹²	0.767	0.576	—	0.957	0.063	BS-600M¹²	2.02	1.52	—	2.52	0.17
	BS-800¹³	0.758	0.568	—	0.949	0.063	BS-800¹³	1.84	1.39	—	2.29	0.15
	BS-1000M¹⁴	0.767	0.576	—	0.957	0.063	BS-1000M¹⁴	2.02	1.52	—	2.52	0.17
	BS-2000¹⁵	0.783	0.588	—	0.979	0.065	BS-2000¹⁵	2.02	1.52	—	2.52	0.17
	BS-2800M¹⁶	0.767	0.576	—	0.957	0.063	BS-2800M¹⁶	2.02	1.52	—	2.52	0.17

mindray

ACE Control

For use on: BS-1000M, BS-2000,
BS-2800M, BS-600M,
BS-800

LOT 057725001

 2026-03-10

