

## Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

- 1.BS-240: BS-230, BS-240;  
 2.BS-240E: BS240E, BS240Pro;  
 3.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with "W9-" and software version starts with "35.00");  
 4.BS-380: BS-380, BS-390;  
 5.BS-400: BS-400, BS-420;  
 6.BS-430: BS-430, BS-450,BS-460;  
 7.BS-480: BS-480, BS-490;

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

- 8.BS-600: BS-600, BS-620;  
 9.BS-600M: BS-600M, BS-620M;  
 10.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;  
 11.BS-1000M: BS-1000M, BS-1100M;  
 12.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;  
 13.BS-2800M:BS-2600M, BS-2800M;  
 14.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction.

LOT : 057824001

: 2025-04-07

| English   | Abbreviated name         | Model   | Unit    | Assay Value        | Range (Assay Value $\pm$ 3SD)                    |
|-----------|--------------------------|---------|---------|--------------------|--------------------------------------------------|
| Русский   | сокращенное наименование | модель  | Прибор  | Результат анализа  | Диапазон(результат анализа $\pm$ 3CO)            |
| Português | Nome abreviado           | Modelo  | Unidade | Valores da análise | Faixa(Valores da análise $\pm$ 3SD)              |
| Español   | nombre abreviado         | modelo  | Unidad  | Valor de ensayo    | Rango(Valor de ensayo $\pm$ 3SD)                 |
| Italiano  | abbreviazione            | modelli | Unità   | Valori di dosaggio | Intervallo (valore di concentrazione $\pm$ 3 DS) |
| Türkçe    | kısaltılmış ad           | model   | Ünite   | Tayin Değeri       | Aralık (Tayin Değeri $\pm$ 3SD)                  |

## D-Dimer

|           |          |
|-----------|----------|
| English   | D-Dimer  |
| Русский   | D-димера |
| Português | D-Dimer  |
| Español   | Dímero D |
| Italiano  | D-Dimero |
| Türkçe    | D-Dimer  |

| Abbreviated name | Unit  | Level: L/N             |             |                         |   |      | Level: H/P |                        |                         |      |      |      |     |
|------------------|-------|------------------------|-------------|-------------------------|---|------|------------|------------------------|-------------------------|------|------|------|-----|
|                  |       | Model                  | Assay Value | Range (Assay Value±3SD) |   | 1 SD | Model      | Assay Value            | Range (Assay Value±3SD) |      | 1 SD |      |     |
| D-Dimer          | µg/mL | BS-240 <sup>1</sup>    | /           | /                       | / | /    | /          | BS-240 <sup>1</sup>    | /                       | /    | /    | /    | /   |
|                  |       | BS-240E <sup>2</sup>   | /           | /                       | / | /    | /          | BS-240E <sup>2</sup>   | /                       | /    | /    | /    | /   |
|                  |       | BS-360E <sup>3</sup>   | /           | /                       | / | /    | /          | BS-360E <sup>3</sup>   | /                       | /    | /    | /    | /   |
|                  |       | BS-380 <sup>4</sup>    | /           | /                       | / | /    | /          | BS-380 <sup>4</sup>    | /                       | /    | /    | /    | /   |
|                  |       | BS-400 <sup>5</sup>    | /           | /                       | / | /    | /          | BS-400 <sup>5</sup>    | /                       | /    | /    | /    | /   |
|                  |       | BS-430 <sup>6</sup>    | /           | /                       | / | /    | /          | BS-430 <sup>6</sup>    | /                       | /    | /    | /    | /   |
|                  |       | BS-480 <sup>7</sup>    | /           | /                       | / | /    | /          | BS-480 <sup>7</sup>    | /                       | /    | /    | /    | /   |
|                  |       | BS-600 <sup>8</sup>    | /           | /                       | / | /    | /          | BS-600 <sup>8</sup>    | /                       | /    | /    | /    | /   |
|                  |       | BS-600M <sup>9</sup>   | /           | /                       | / | /    | /          | BS-600M <sup>9</sup>   | /                       | /    | /    | /    | /   |
|                  |       | BS-800 <sup>10</sup>   | /           | /                       | / | /    | /          | BS-800 <sup>10</sup>   | /                       | /    | /    | /    | /   |
|                  |       | BS-1000M <sup>11</sup> | 1.40        | 0.98                    | — | 1.82 | 0.14       | BS-1000M <sup>11</sup> | 15.2                    | 12.8 | —    | 17.6 | 0.8 |
|                  |       | BS-2000 <sup>12</sup>  | /           | /                       | / | /    | /          | BS-2000 <sup>12</sup>  | /                       | /    | /    | /    | /   |
|                  |       | BS-2800M <sup>13</sup> | /           | /                       | / | /    | /          | BS-2800M <sup>13</sup> | /                       | /    | /    | /    | /   |

**mindray**  
**D-Dimer Control**

For use on: BS-1000M, BS-2000,  
BS-2800M, BS-600M,  
BS-800

LOT# 057824001

EXP 2025-04-07

