

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

- 1.BS-240: BS-230, BS-240;
 2.BS-240E: BS240E, BS240Pro;
 3.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with "W9-" and software version starts with "35.00");
 4.BS-380: BS-380, BS-390;
 5.BS-400: BS-400, BS-420;
 6.BS-430: BS-430, BS-450,BS-460;
 7.BS-480: BS-480, BS-490;

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

- 8.BS-600: BS-600, BS-620;
 9.BS-600M: BS-600M, BS-620M;
 10.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;
 11.BS-1000M: BS-1000M, BS-1100M;
 12.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;
 13.BS-2800M:BS-2600M, BS-2800M;
 14.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction.

LOT: 057824004

EXP: 2025-06-15

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value $\pm$ 3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа $\pm$ 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise $\pm$ 3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo $\pm$ 3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione $\pm$ 3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri $\pm$ 3SD)

D-Dimer

English	D-Dimer
Русский	D-димера
Português	D-Dimer
Español	Dímero D
Italiano	D-Dimero
Türkçe	D-Dimer

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
D-Dimer	µg/mL	BS-240 ¹	/	/	/	/	/	BS-240 ¹	/	/	/	/	/
		BS-240E ²	/	/	/	/	/	BS-240E ²	/	/	/	/	/
		BS-360E ³	/	/	/	/	/	BS-360E ³	/	/	/	/	/
		BS-380 ⁴	/	/	/	/	/	BS-380 ⁴	/	/	/	/	/
		BS-400 ⁵	/	/	/	/	/	BS-400 ⁵	/	/	/	/	/
		BS-430 ⁶	/	/	/	/	/	BS-430 ⁶	/	/	/	/	/
		BS-480 ⁷	/	/	/	/	/	BS-480 ⁷	/	/	/	/	/
		BS-600 ⁸	/	/	/	/	/	BS-600 ⁸	/	/	/	/	/
		BS-600M ⁹	/	/	/	/	/	BS-600M ⁹	/	/	/	/	/
		BS-800 ¹⁰	/	/	/	/	/	BS-800 ¹⁰	/	/	/	/	/
		BS-1000M ¹¹	1.40	0.98	—	1.82	0.14	BS-1000M ¹¹	15.2	12.8	—	17.6	0.8
		BS-2000 ¹²	/	/	/	/	/	BS-2000 ¹²	/	/	/	/	/
		BS-2800M ¹³	/	/	/	/	/	BS-2800M ¹³	/	/	/	/	/

mindray
D-Dimer Control

For use on: BS-1000M, BS-2000,
BS-2800M, BS-600M,
BS-800

 **LOT** 057824004

 **2025-06-15**

