

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

- 1.BS-200E: BS-200E, BS-220E;;
2.BS-240: BS-230, BS-240;
3.BS-240E: BS240E, BS240Pro;
4.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00");
5.BS-380: BS-380, BS-390;
6.BS-400: BS-400, BS-420;
7.BS-430: BS-430, BS-450,BS-460;
8.BS-480: BS-480, BS-490;

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

- 9.BS-600: BS-600, BS-620;
10.BS-600M: BS-600M, BS-620M;
11.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;
12.BS-1000M: BS-1000M, BS-1100M;
13.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;
14.BS-2800M:BS-2600M, BS-2800M;
15.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction.

LOT : 058023005 2024-10-20

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)
FER			IgE		
English		Ferritin			Immunoglobulin E
Русский		ферритин			Иммуноглобулин E
Português		Ferritina			Imunoglobulina E
Español		Ferritina			Inmunoglobulina E
Italiano		Ferritina			Immunoglobulina E
Türkçe		Ferritin			İmmünglobulin E
MYO					
English		Myoglobin			
Русский		Миоглобин			
Português		Mioglobina			
Español		Mioglobina			
Italiano		Mioglobina			
Türkçe		Myoglobin			

Abbreviated name	Unit	Level: L/N					Level: H/P				
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD
ng/mL		BS-200E ¹	/	/	/	/	BS-200E ¹	/	/	/	/
		BS-240 ²	/	/	/	/	BS-240 ²	/	/	/	/
		BS-240E ³	/	/	/	/	BS-240E ³	/	/	/	/
		BS-360E ⁴	/	/	/	/	BS-360E ⁴	/	/	/	/
		BS-380 ⁵	/	/	/	/	BS-380 ⁵	/	/	/	/
		BS-400 ⁶	/	/	/	/	BS-400 ⁶	/	/	/	/
		BS-430 ⁷	/	/	/	/	BS-430 ⁷	/	/	/	/
		BS-480 ⁸	/	/	/	/	BS-480 ⁸	/	/	/	/
		BS-600 ⁹	/	/	/	/	BS-600 ⁹	/	/	/	/
		BS-600M ¹⁰	/	/	/	/	BS-600M ¹⁰	/	/	/	/
		BS-800 ¹¹	/	/	/	/	BS-800 ¹¹	/	/	/	/

FER	pmol/L	BS-1000M ¹²	123	105	—	141	6	BS-1000M ¹²	459	390	—	528	23
		BS-2000 ¹³	/	/	/	/	/	BS-2000 ¹³	/	/	/	/	/
		BS-2800M ¹⁴	/	/	/	/	/	BS-2800M ¹⁴	/	/	/	/	/
		BS-200E ¹	/	/	/	/	/	BS-200E ¹	/	/	/	/	/
		BS-240 ²	/	/	/	/	/	BS-240 ²	/	/	/	/	/
		BS-240E ³	/	/	/	/	/	BS-240E ³	/	/	/	/	/
		BS-360E ⁴	/	/	/	/	/	BS-360E ⁴	/	/	/	/	/
		BS-380 ⁵	/	/	/	/	/	BS-380 ⁵	/	/	/	/	/
		BS-400 ⁶	/	/	/	/	/	BS-400 ⁶	/	/	/	/	/
		BS-430 ⁷	/	/	/	/	/	BS-430 ⁷	/	/	/	/	/
		BS-480 ⁸	/	/	/	/	/	BS-480 ⁸	/	/	/	/	/
		BS-600 ⁹	/	/	/	/	/	BS-600 ⁹	/	/	/	/	/
		BS-600M ¹⁰	/	/	/	/	/	BS-600M ¹⁰	/	/	/	/	/
		BS-800 ¹¹	/	/	/	/	/	BS-800 ¹¹	/	/	/	/	/
		BS-1000M ¹²	276	236	—	317	13	BS-1000M ¹²	1031	876	—	1186	52
		BS-2000 ¹³	/	/	/	/	/	BS-2000 ¹³	/	/	/	/	/
		BS-2800M ¹⁴	/	/	/	/	/	BS-2800M ¹⁴	/	/	/	/	/
Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
IgE	IU/mL	BS-240 ²	/	/	/	/	BS-240 ²	/	/	/	/		
		BS-240E ³	/	/	/	/	BS-240E ³	/	/	/	/		
		BS-360E ⁴	/	/	/	/	BS-360E ⁴	/	/	/	/		
		BS-380 ⁵	/	/	/	/	BS-380 ⁵	/	/	/	/		
		BS-400 ⁶	/	/	/	/	BS-400 ⁶	/	/	/	/		
		BS-480 ⁸	/	/	/	/	BS-480 ⁸	/	/	/	/		
		BS-600 ⁹	/	/	/	/	BS-600 ⁹	/	/	/	/		
		BS-600M ¹⁰	/	/	/	/	BS-600M ¹⁰	/	/	/	/		
		BS-800 ¹¹	/	/	/	/	BS-800 ¹¹	/	/	/	/		
		BS-1000M ¹²	80.9	56.6	—	105.2	8.1	BS-1000M ¹²	540	378	—	702	54
		BS-2000 ¹³	/	/	/	/	/	BS-2000 ¹³	/	/	/	/	/
		BS-2800M ¹⁴	/	/	/	/	/	BS-2800M ¹⁴	/	/	/	/	/
Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
MYO	ng/mL	BS-240 ²	/	/	/	/	BS-240 ²	/	/	/	/		
		BS-240E ³	/	/	/	/	BS-240E ³	/	/	/	/		
		BS-360E ⁴	/	/	/	/	BS-360E ⁴	/	/	/	/		
		BS-380 ⁵	/	/	/	/	BS-380 ⁵	/	/	/	/		
		BS-400 ⁶	/	/	/	/	BS-400 ⁶	/	/	/	/		
		BS-430 ⁷	/	/	/	/	BS-430 ⁷	/	/	/	/		
		BS-480 ⁸	/	/	/	/	BS-480 ⁸	/	/	/	/		
		BS-600 ⁹	/	/	/	/	BS-600 ⁹	/	/	/	/		
		BS-600M ¹⁰	/	/	/	/	BS-600M ¹⁰	/	/	/	/		
		BS-800 ¹¹	/	/	/	/	BS-800 ¹¹	/	/	/	/		
		BS-1000M ¹²	79.8	67.8	—	91.8	4.0	BS-1000M ¹²	289	247	—	331	14
		BS-2000 ¹³	/	/	/	/	/	BS-2000 ¹³	/	/	/	/	/
		BS-2800M ¹⁴	/	/	/	/	/	BS-2800M ¹⁴	/	/	/	/	/

mindray

Multimmun Control

For use on: BS-1000M, BS-2000,
BS-2800M, BS-600M,
BS-800

LOT 058023005

 2024-10-20

