

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

1.BS-200E: BS-200E, BS-220E,;

2.BS-240: BS-230, BS-240;

3.BS-240E: BS240E, BS240Pro;

4.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-

330E(V35.00)(Serial Number starts with "W8-" and software version

starts with "35.00"), BS-350E(V35.00) (Serial Number starts with "W9-"

and software version starts with "35.00");

5.BS-380: BS-380, BS-390;

6.BS-400: BS-400, BS-420;

7.BS-430: BS-430, BS-450,BS-460;

8.BS-480: BS-480, BS-490;

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

9.BS-600: BS-600, BS-620;

10.BS-600M: BS-600M, BS-620M;

11.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-

1800plus;

12.BS-1000M: BS-1000M, BS-1100M;

13.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

14.BS-2800M:BS-2600M, BS-2800M;

15.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction.

LOT : 058023008

EXP: 2025-01-06

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value $\pm$ 3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа $\pm$ 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise $\pm$ 3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo $\pm$ 3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione $\pm$ 3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri $\pm$ 3SD)

FER				IgE	
English		Ferritin		Immunoglobulin E	
Русский		ферритин		Иммуноглобулин E	
Português		Ferritina		Imunoglobulina E	
Español		Ferritina		Inmunoglobulina E	
Italiano		Ferritina		Immunoglobulina E	
Türkçe		Ferritin		İmmünglobulin E	
MYO					
English		Myoglobin			
Русский		Миоглобин			
Português		Mioglobina			
Español		Mioglobina			
Italiano		Mioglobina			
Türkçe		Myoglobin			

Abbreviated name	Unit	Level: L/N					Level: H/P				
		Model	Assay Value	Range (Assay Value $\pm$ 3SD)		1 SD	Model	Assay Value	Range (Assay Value $\pm$ 3SD)		1 SD
ng/mL		BS-200E ¹	/	/	/	/	BS-200E ¹	/	/	/	/
		BS-240 ²	/	/	/	/	BS-240 ²	/	/	/	/
		BS-240E ³	/	/	/	/	BS-240E ³	/	/	/	/
		BS-360E ⁴	/	/	/	/	BS-360E ⁴	/	/	/	/
		BS-380 ⁵	/	/	/	/	BS-380 ⁵	/	/	/	/
		BS-400 ⁶	/	/	/	/	BS-400 ⁶	/	/	/	/
		BS-430 ⁷	/	/	/	/	BS-430 ⁷	/	/	/	/
		BS-480 ⁸	/	/	/	/	BS-480 ⁸	/	/	/	/
		BS-600 ⁹	/	/	/	/	BS-600 ⁹	/	/	/	/
		BS-600M ¹⁰	/	/	/	/	BS-600M ¹⁰	/	/	/	/
		BS-800 ¹¹	/	/	/	/	BS-800 ¹¹	/	/	/	/

FER	pmol/L	BS-1000M ¹²	123	105	—	141	6	BS-1000M ¹²	459	390	—	528	23
		BS-2000 ¹³	/	/	/	/	/	BS-2000 ¹³	/	/	/	/	/
		BS-2800M ¹⁴	/	/	/	/	/	BS-2800M ¹⁴	/	/	/	/	/
		BS-200E ¹	/	/	/	/	/	BS-200E ¹	/	/	/	/	/
		BS-240 ²	/	/	/	/	/	BS-240 ²	/	/	/	/	/
		BS-240E ³	/	/	/	/	/	BS-240E ³	/	/	/	/	/
		BS-360E ⁴	/	/	/	/	/	BS-360E ⁴	/	/	/	/	/
		BS-380 ⁵	/	/	/	/	/	BS-380 ⁵	/	/	/	/	/
		BS-400 ⁶	/	/	/	/	/	BS-400 ⁶	/	/	/	/	/
		BS-430 ⁷	/	/	/	/	/	BS-430 ⁷	/	/	/	/	/
		BS-480 ⁸	/	/	/	/	/	BS-480 ⁸	/	/	/	/	/
		BS-600 ⁹	/	/	/	/	/	BS-600 ⁹	/	/	/	/	/
		BS-600M ¹⁰	/	/	/	/	/	BS-600M ¹⁰	/	/	/	/	/
		BS-800 ¹¹	/	/	/	/	/	BS-800 ¹¹	/	/	/	/	/
		BS-1000M ¹²	276	236	—	317	13	BS-1000M ¹²	1031	876	—	1186	52
		BS-2000 ¹³	/	/	/	/	/	BS-2000 ¹³	/	/	/	/	/
		BS-2800M ¹⁴	/	/	/	/	/	BS-2800M ¹⁴	/	/	/	/	/
Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
IgE	IU/mL	BS-240 ²	/	/	/	/	BS-240 ²	/	/	/	/	/	/
		BS-240E ³	/	/	/	/	BS-240E ³	/	/	/	/	/	/
		BS-360E ⁴	/	/	/	/	BS-360E ⁴	/	/	/	/	/	/
		BS-380 ⁵	/	/	/	/	BS-380 ⁵	/	/	/	/	/	/
		BS-400 ⁶	/	/	/	/	BS-400 ⁶	/	/	/	/	/	/
		BS-480 ⁸	/	/	/	/	BS-480 ⁸	/	/	/	/	/	/
		BS-600 ⁹	/	/	/	/	BS-600 ⁹	/	/	/	/	/	/
		BS-600M ¹⁰	/	/	/	/	BS-600M ¹⁰	/	/	/	/	/	/
		BS-800 ¹¹	/	/	/	/	BS-800 ¹¹	/	/	/	/	/	/
		BS-1000M ¹²	80.9	56.6	—	105.2	8.1	BS-1000M ¹²	540	378	—	702	54
		BS-2000 ¹³	/	/	/	/	BS-2000 ¹³	/	/	/	/	/	/
		BS-2800M ¹⁴	/	/	/	/	BS-2800M ¹⁴	/	/	/	/	/	/
Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
MYO	ng/mL	BS-240 ²	/	/	/	/	BS-240 ²	/	/	/	/	/	/
		BS-240E ³	/	/	/	/	BS-240E ³	/	/	/	/	/	/
		BS-360E ⁴	/	/	/	/	BS-360E ⁴	/	/	/	/	/	/
		BS-380 ⁵	/	/	/	/	BS-380 ⁵	/	/	/	/	/	/
		BS-400 ⁶	/	/	/	/	BS-400 ⁶	/	/	/	/	/	/
		BS-430 ⁷	/	/	/	/	BS-430 ⁷	/	/	/	/	/	/
		BS-480 ⁸	/	/	/	/	BS-480 ⁸	/	/	/	/	/	/
		BS-600 ⁹	/	/	/	/	BS-600 ⁹	/	/	/	/	/	/
		BS-600M ¹⁰	/	/	/	/	BS-600M ¹⁰	/	/	/	/	/	/
		BS-800 ¹¹	/	/	/	/	BS-800 ¹¹	/	/	/	/	/	/
		BS-1000M ¹²	79.8	67.8	—	91.8	4.0	BS-1000M ¹²	289	247	—	331	14
		BS-2000 ¹³	/	/	/	/	BS-2000 ¹³	/	/	/	/	/	/
		BS-2800M ¹⁴	/	/	/	/	BS-2800M ¹⁴	/	/	/	/	/	/

mindray

Multimmun Control

For use on: BS-1000M, BS-2000,
BS-2800M, BS-600M,
BS-800

LOT 058023008

 2025-01-06

