

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-200E: BS-200E, BS-220E,;

9.BS-600: BS-600, BS-620;

2.BS-240: BS-230, BS-240;

10.BS-600M: BS-600M, BS-620M;

3.BS-240E: BS240E, BS240Pro;

11.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

4.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-

12.BS-1000M: BS-1000M, BS-1100M;

330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with "W9-" and software version starts with "35.00");

13.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

5.BS-380: BS-380, BS-390;

14.BS-2800M:BS-2600M, BS-2800M;

6.BS-400: BS-400, BS-420;

15.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction.

7.BS-430: BS-430, BS-450,BS-460;

8.BS-480: BS-480, BS-490;

LOT : 059025004

: 2026-07-13

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value $\pm$ 3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа $\pm$ 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise $\pm$ 3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo $\pm$ 3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione $\pm$ 3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri $\pm$ 3SD)

CysC

English	Cystatin C
Русский	цистатин C
Português	Cistatina C
Español	cistatina C
Italiano	cistatina C
Türkçe	Cystatin C

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
CysC II	mg/L	BS-200E ¹	0.748	0.637	—	0.859	0.037	BS-200E ¹	1.30	1.09	—	1.51	0.07
		BS-240 ²	0.733	0.622	—	0.844	0.037	BS-240 ²	1.27	1.09	—	1.45	0.06
		BS-240E ³	0.774	0.657	—	0.891	0.039	BS-240E ³	1.27	1.09	—	1.45	0.06
		BS-360E ⁴	0.748	0.637	—	0.859	0.037	BS-360E ⁴	1.30	1.09	—	1.51	0.07
		BS-380 ⁵	0.769	0.655	—	0.883	0.038	BS-380 ⁵	1.30	1.09	—	1.51	0.07
		BS-400 ⁶	0.769	0.655	—	0.883	0.038	BS-400 ⁶	1.30	1.09	—	1.51	0.07
		BS-430 ⁷	0.748	0.637	—	0.859	0.037	BS-430 ⁷	1.30	1.09	—	1.51	0.07
		BS-480 ⁸	0.769	0.655	—	0.883	0.038	BS-480 ⁸	1.30	1.09	—	1.51	0.07
		BS-600 ⁹	0.769	0.655	—	0.883	0.038	BS-600 ⁹	1.30	1.09	—	1.51	0.07
		BS-600M ¹⁰	0.777	0.660	—	0.894	0.039	BS-600M ¹⁰	1.32	1.11	—	1.53	0.07
		BS-800 ¹¹	0.748	0.637	—	0.859	0.037	BS-800 ¹¹	1.30	1.09	—	1.51	0.07
		BS-1000M ¹²	0.769	0.655	—	0.883	0.038	BS-1000M ¹²	1.30	1.09	—	1.51	0.07
		BS-2000 ¹³	0.762	0.648	—	0.876	0.038	BS-2000 ¹³	1.31	1.10	—	1.52	0.07
		BS-2800M ¹⁴	0.769	0.655	—	0.883	0.038	BS-2800M ¹⁴	1.30	1.09	—	1.51	0.07

mindray
Cystatin C Control

For use on: BS-1000M, BS-2000,
BS-2800M, BS-600M,
BS-800

LOT 059025004

2026-07-13

