

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.	Русский: Данные совпадают во всех группах.
Português: A dados de cada grupo é a mesma.	Español: la datos de cada grupo es la misma.
Italiano: la dati di ogni gruppo è la stessa.	Türkçe: her grubun veri aynıdır.
1.BS-120: BS-120, BS-130, BS-180, BS-190;	10.BS-480: BS-480, BS-490;
2.BS-200: BS-200, BS-220, BS-330, BS-350;	11.BS-600: BS-600, BS-620;
3.BS-240: BS-230, BS-240;	12.BS-600M: BS-600M, BS-620M;
4.BS-240E: BS240E, BS240Pro;	13.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;
5.BS-330E:BS-330E(Serial Number starts with "XQ-"), BS-350E(Serial Number starts with "XS-")	14.BS-1000M: BS-1000M, BS-1100M;
6.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00");	15.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;
7.BS-380: BS-380, BS-390;	16.BS-2800M:BS-2600M, BS-2800M;
8.BS-400: BS-400, BS-420;	17.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction.
9.BS-430: BS-430, BS-450,BS-460;	

LOT : 060823007

: 2024-10-19

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

HCY	
English	Homocysteine (HCY)
Русский	Гомоцистеин (HCY)
Português	Homocisteína (HCY)
Español	Homocisteína (HCY)
Italiano	Omocisteina (HCY)
Türkçe	Homosistain (HCY)

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)			1 SD	Model	Assay Value	Range (Assay Value±3SD)			1 SD
HCY II	μmol/L	BS-120 ¹	/	/	—	/	/	BS-120 ¹	/	/	—	/	/
		BS-200 ²	/	/	—	/	/	BS-200 ²	/	/	—	/	/
		BS-240 ³	/	/	—	/	/	BS-240 ³	/	/	—	/	/
		BS-240E ⁴	/	/	—	/	/	BS-240E ⁴	/	/	—	/	/
		BS-330E ⁵	/	/	—	/	/	BS-330E ⁵	/	/	—	/	/
		BS-360E ⁶	/	/	—	/	/	BS-360E ⁶	/	/	—	/	/
		BS-380 ⁷	/	/	—	/	/	BS-380 ⁷	/	/	—	/	/
		BS-400 ⁸	/	/	—	/	/	BS-400 ⁸	/	/	—	/	/
		BS-430 ⁹	/	/	—	/	/	BS-430 ⁹	/	/	—	/	/
		BS-480 ¹⁰	/	/	—	/	/	BS-480 ¹⁰	/	/	—	/	/
		BS-600 ¹¹	/	/	—	/	/	BS-600 ¹¹	/	/	—	/	/
		BS-600M ¹²	/	/	—	/	/	BS-600M ¹²	/	/	—	/	/
		BS-800 ¹³	/	/	—	/	/	BS-800 ¹³	/	/	—	/	/
		BS-1000M ¹⁴	16.1	12.8	—	19.4	1.1	BS-1000M ¹⁴	32.1	25.8	—	38.4	2.1
		BS-2000 ¹⁵	/	/	—	/	/	BS-2000 ¹⁵	/	/	—	/	/

mindray

HCY Control

For use on: BS-1000M, BS-2000,
BS-2800M, BS-600M,
BS-800

LOT 060823007

 2024-10-19

