

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-120: BS-120, BS-130, BS-180, BS-190;

10.BS-480: BS-480, BS-490;

2.BS-200: BS-200, BS-220, BS-330, BS-350;

11.BS-600: BS-600, BS-620;

3.BS-240: BS-230, BS-240;

12.BS-600M: BS-600M, BS-620M;

4.BS-240E: BS240E, BS240Pro;

5.BS-330E: BS-330E (Serial Number starts with "XQ-"), BS-350E (Serial Number starts with "XS-")

13.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

6.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-

14.BS-1000M: BS-1000M, BS-1100M;

330E(V35.00) (Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with "W9-" and software version starts with "35.00");

15.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

7.BS-380: BS-380, BS-390;

16.BS-2800M: BS-2600M, BS-2800M;

8.BS-400: BS-400, BS-420;

17.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction.

9.BS-430: BS-430, BS-450, BS-460;

LOT: 060825006

EXP: 2026-09-26

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value $\pm$ 3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа $\pm$ 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise $\pm$ 3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo $\pm$ 3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione $\pm$ 3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri $\pm$ 3SD)

HCY

English	Homocysteine (HCY)
Русский	Гомоцистеин (HCY)
Português	Homocisteína (HCY)
Español	Homocisteína (HCY)
Italiano	Omocisteina (HCY)
Türkçe	Homosistain (HCY)

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
HCY II	μmol/L	BS-120 ¹	15.4	12.4	—	18.4	1.0	BS-120 ¹	31.9	25.6	—	38.2	2.1
		BS-200 ²	15.4	12.4	—	18.4	1.0	BS-200 ²	31.6	25.3	—	37.9	2.1
		BS-240 ³	15.4	12.4	—	18.4	1.0	BS-240 ³	31.9	25.6	—	38.2	2.1
		BS-240E ⁴	15.5	12.5	—	18.5	1.0	BS-240E ⁴	31.9	25.6	—	38.2	2.1
		BS-330E ⁵	15.5	12.5	—	18.5	1.0	BS-330E ⁵	30.7	24.7	—	36.7	2.0
		BS-360E ⁶	15.5	12.5	—	18.5	1.0	BS-360E ⁶	31.9	25.6	—	38.2	2.1
		BS-380 ⁷	15.4	12.4	—	18.4	1.0	BS-380 ⁷	32.0	25.7	—	38.3	2.1
		BS-400 ⁸	15.4	12.4	—	18.4	1.0	BS-400 ⁸	32.0	25.7	—	38.3	2.1
		BS-430 ⁹	15.5	12.5	—	18.5	1.0	BS-430 ⁹	31.9	25.6	—	38.2	2.1
		BS-480 ¹⁰	15.4	12.4	—	18.4	1.0	BS-480 ¹⁰	32.0	25.7	—	38.3	2.1
		BS-600 ¹¹	15.4	12.4	—	18.4	1.0	BS-600 ¹¹	32.0	25.7	—	38.3	2.1
		BS-600M ¹²	15.4	12.4	—	18.4	1.0	BS-600M ¹²	32.0	25.7	—	38.3	2.1
		BS-800 ¹³	15.5	12.5	—	18.5	1.0	BS-800 ¹³	31.9	25.6	—	38.2	2.1
		BS-1000M ¹⁴	15.4	12.4	—	18.4	1.0	BS-1000M ¹⁴	32.0	25.7	—	38.3	2.1
		BS-2000 ¹⁵	15.7	12.7	—	18.7	1.0	BS-2000 ¹⁵	32.2	25.9	—	38.5	2.1

BS-2800M ¹⁶	15.4	12.4	—	18.4	1.0	BS-2800M ¹⁶	32.0	25.7	—	38.3	2.1
------------------------	------	------	---	------	-----	------------------------	------	------	---	------	-----

mindray

HCY Control

For use on: BS-1000M, BS-2000,
BS-2800M, BS-600M,
BS-800

LOT 060825006

 **2026-09-26**

