

## Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-120: BS-120, BS-130, BS-180, BS-190;

11.BS-430: BS-430, BS-450,BS-460;

2.BS-200: BS-200, BS-220, BS-330, BS-350;

12.BS-480: BS-480, BS-490;

3.BS-200E: BS-200E, BS-220E,;

13.BS-600: BS-600, BS-620;

4.BS-240: BS-230, BS-240;

14.BS-600M: BS-600M, BS-620M;

5.BS-240E: BS240E, BS240Pro;

15.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

6.BS-300: BS-300, BS-320;

16.BS-1000M: BS-1000M, BS-1100M;

7.BS-330E:BS-330E(Serial Number starts with "XQ-"), BS-350E(Serial Number starts with "XS-")

17.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

8.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-

18.BS-2800M:BS-2600M, BS-2800M;

330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00");

9.BS-380: BS-380, BS-390;

19.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction.

10.BS-400: BS-400, BS-420;

LOT : 061123008

2025-06-06

| English   | Abbreviated name         | Model   | Unit    | Assay Value        | Range (Assay Value $\pm$ 3SD)                       |
|-----------|--------------------------|---------|---------|--------------------|-----------------------------------------------------|
| Русский   | сокращенное наименование | модель  | Прибор  | Результат анализа  | Диапазон(результат анализа $\pm$ 3CO)               |
| Português | Nome abreviado           | Modelo  | Unidade | Valores da análise | Faixa(Valores da análise $\pm$ 3SD)                 |
| Español   | nombre abreviado         | modelo  | Unidad  | Valor de ensayo    | Rango(Valor de ensayo $\pm$ 3SD)                    |
| Italiano  | abbreviazione            | modelli | Unità   | Valori di dosaggio | Intervallo<br>(valore di concentrazione $\pm$ 3 DS) |
| Türkçe    | kisaltılmış ad           | model   | Ünite   | Tayin Değeri       | Aralık (Tayin Değeri $\pm$ 3SD)                     |

## Lp(a) II

|           |                      |
|-----------|----------------------|
| English   | Lp(a) II Calibrator  |
| Русский   | Калибратор Lp(a) II  |
| Português | Calibrador Lp(a) II  |
| Español   | Calibrador Lp(a) II  |
| Italiano  | Calibratore Lp(a) II |
| Türkçe    | Lp(a) II Kalibratörü |

| Abbreviated name | Unit   | Level: L/N            |             |                         |   |   | Level: H/P |                       |             |                         |   |   |      |
|------------------|--------|-----------------------|-------------|-------------------------|---|---|------------|-----------------------|-------------|-------------------------|---|---|------|
|                  |        | Model                 | Assay Value | Range (Assay Value±3SD) |   |   | 1 SD       | Model                 | Assay Value | Range (Assay Value±3SD) |   |   | 1 SD |
| Lp(a) II         | nmol/L | BS-120 <sup>1</sup>   | /           | /                       | / | / | /          | BS-120 <sup>1</sup>   | /           | /                       | / | / | /    |
|                  |        | BS-200 <sup>2</sup>   | /           | /                       | / | / | /          | BS-200 <sup>2</sup>   | /           | /                       | / | / | /    |
|                  |        | BS-200E <sup>3</sup>  | /           | /                       | / | / | /          | BS-200E <sup>3</sup>  | /           | /                       | / | / | /    |
|                  |        | BS-240 <sup>4</sup>   | /           | /                       | / | / | /          | BS-240 <sup>4</sup>   | /           | /                       | / | / | /    |
|                  |        | BS-240E <sup>5</sup>  | /           | /                       | / | / | /          | BS-240E <sup>5</sup>  | /           | /                       | / | / | /    |
|                  |        | BS-300 <sup>6</sup>   | /           | /                       | / | / | /          | BS-300 <sup>6</sup>   | /           | /                       | / | / | /    |
|                  |        | BS-330E <sup>7</sup>  | /           | /                       | / | / | /          | BS-330E <sup>7</sup>  | /           | /                       | / | / | /    |
|                  |        | BS-360E <sup>8</sup>  | /           | /                       | / | / | /          | BS-360E <sup>8</sup>  | /           | /                       | / | / | /    |
|                  |        | BS-380 <sup>9</sup>   | /           | /                       | / | / | /          | BS-380 <sup>9</sup>   | /           | /                       | / | / | /    |
|                  |        | BS-400 <sup>10</sup>  | /           | /                       | / | / | /          | BS-400 <sup>10</sup>  | /           | /                       | / | / | /    |
|                  |        | BS-430 <sup>11</sup>  | /           | /                       | / | / | /          | BS-430 <sup>11</sup>  | /           | /                       | / | / | /    |
|                  |        | BS-480 <sup>12</sup>  | /           | /                       | / | / | /          | BS-480 <sup>12</sup>  | /           | /                       | / | / | /    |
|                  |        | BS-600 <sup>13</sup>  | /           | /                       | / | / | /          | BS-600 <sup>13</sup>  | /           | /                       | / | / | /    |
|                  |        | BS-600M <sup>14</sup> | /           | /                       | / | / | /          | BS-600M <sup>14</sup> | /           | /                       | / | / | /    |

|                        |      |      |   |      |     |                        |     |    |   |     |    |
|------------------------|------|------|---|------|-----|------------------------|-----|----|---|-----|----|
| BS-800 <sup>15</sup>   | /    | /    | / | /    | /   | BS-800 <sup>15</sup>   | /   | /  | / | /   | /  |
| BS-1000M <sup>16</sup> | 44.6 | 29.6 | — | 59.6 | 5.0 | BS-1000M <sup>16</sup> | 124 | 82 | — | 166 | 14 |
| BS-2000 <sup>17</sup>  | /    | /    | / | /    | /   | BS-2000 <sup>17</sup>  | /   | /  | / | /   | /  |
| BS-2800M <sup>18</sup> | /    | /    | / | /    | /   | BS-2800M <sup>18</sup> | /   | /  | / | /   | /  |

# **mindray**

## **Lipoprotein (a) Control**

For use on: BS-1000M, BS-2000,  
BS-2800M, BS-600M,  
BS-800

**LOT** 061123008

 **2025-06-06**

